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Alcohols / alcools / Alkohole

Redox, eVALC

These specially designed filters shall be used to provide protection against chemical perturbation. Mineral acids, halogens, powerful oxidizing agents and corrosive substances may damage the cells and reduce their lifetime. Chemical filters shall be replaced, when a recalibration of the cell becomes necessary more often, when replacing a cell, but latest after one year. The eVALC needs not to be replaced, if not mechanically damaged.

The higher the concentration of an aggressive substance, the less effective will be the protection capability of a chemical filter, so the cell lifetime still can be reduced. If the concentration exceeds certain limits, no protection is guaranteed.

-> Only practical testing can answer the question about exchange cycles of a filter.

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Novasina FILTER SYSTEMS

eVC-21 Chemical Filter

acetic acid (eg. mustard)
 formic acid
 butyric acid
 lactic acid
 sulphuric acid
 sulfur dioxide
 hydrogen sulfide
 chlorine
 carboxylic acids, aliphatic
 hydrogen halides
 hydrogen peroxide
 formaldehyde
 jodine
 chlorine, hypochloric acid

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eVC-26 Universal Filter

oil fog
 fine dusts
 amines
 ozone
 nicotine
 aromatic hydrocarbons
 perchlor-, trichlor-ethylene
 chloroform
 nitrogen oxides

Art.nr. 1111003

Oelnebel	
Feinstaub	
Amine	
Ozon	(O ₃)
Nikotin	
arom. Kohlenwasserstoffe	(Benzol, Kerosen, Toluol)
Perchlor-, Trichlorethylen	(C ₂ HCl ₃)
Chloroform	(CHCl ₃)
Stickoxyde	(NO _x)

Redox Filter

primary, secondary alcohols
 glycerol
 ethanol (ethyl alcohol)
 isopropanol (isopropyl alcohol)
 propanol (propyl alcohol)
 propylen-glycol
 polyethylen-glycol
 methanol
 phenol
 yeast and hop aromas
 coffee aromas
 artificial fruit aromas
 onions, garlic, horseraddish
 furanes (bread aromas)
 ketones (eg. acetone or MEK)
 partly *amines*
 perfumes, flavours, fragrances

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