

Scan[®]

Automatic colony counters & inhibition zone readers



interscience



interscience

Our quality for your lab

- Designer and manufacturer for microbiological analyses
- Made in France
- R&D leadership with innovative & reliable products
- Worldwide distribution network, immediate delivery

Scan[®] colony counters

High quality analyses,
full traceability

Scan[®] 300, Scan[®] 500 & Scan[®] 1200

High technology automatic colony counters


With a digital camera and high technology software, they can be linked to a PC via a FireWire connection. They count all colonies on a Petri dish in less than 1 second and provide a complete, fast, accurate and traceable reading of the result.

■ Bacterial enumeration

- **Food analyses**
- **Total flora analyses:** bacterial enumeration, aerobic & anaerobic, yeasts, lactobacillus...
- **Pathogenic bacteria research**
- **Environmental research**
- **Pharmaceutical analyses**
- **Medical analyses**
- **Cosmetics analyses**

■ Inhibition zones

- **Pharmaceutical industry, medical research & hospitals**
(antibiograms, resistance tests to pathogenic microbes, medical diagnoses...)
- **Food industry**
(Tests on lactic ferments & for dairy ingredients industry...)




Scan[®] colony counters

- Automatic colony counters: No settings
- Inhibition zone readers*
- Data traceability and full report



High performance

- > Count colonies of numerous media
- > Read chromogenic Petri dishes*: Colored differentiation of colonies (up to 7 different colors on the same dish)
- > Inhibition zone measurement



Live image

- > Fits any type of dish: automatic adjustment of contrast and lighting
- > High-definition color image
- > Each colony is marked with a cross
- > Powerful zoom: up to x28 for **Scan[®] 1200**

Sample N°	Count	Dilution
E COLI	174	1 2
E COLI	353	1 5
SPIRAL	47	1/1000 9
PETRIFILM	89	1 8
RIDA COUNT	179	1 1
Filteing Membrane	111	1 1

Instant results

- > 1000 colonies detected in 1 second
- > Counts 30 dishes in 5 minutes (in real condition with presetting)
- > Reproducible and standardized results
- > **Scan[®]** results: instant and automatic



* on Scan[®] 500 & Scan[®] 1200



Easy-to-use

- Counting in 1 click <
- All functions in 1 single window <
- Custom parameters: day, users, project... <




Dark Field technology

- Display of every colony <
- Optimized lighting & contrast <
- Long lasting LED lighting <
- 6 lighting combinations <



Traceability & reporting

- Automatic archiving and printing of data: <
 - pictures, comments & results
- Export to EXCEL™, PDF, JPEG, BMP <
- Barcode reader <
- Connection to LIMS network <



**DOWNLOAD
SCAN® SOFTWARE**

www.interscience.fr




** Free update of the software during guarantee period

Efficiency & time saving

Instant results

Thanks to the live image display of the Petri dish on your computer, count more than **1000 CFU/sec.** on all media. Each counted colony is marked with a cross and the result is automatically saved.



No settings

Choose your pre-set parameters for Petri dishes:




Also available on **Scan® 1200**:


- RIDA™ COUNT - Sanita-kun™: AC, CC, EC/CC, SA
- Petrifilm™: AC, ETB, CC, EC/CC, EC
- Compact Dry™: TC, CF, EC, ETB

High-performance colony counters


Scan® works for every kind of colony. The minimum size is **0.05 mm for Scan® 1200** and **0.1 mm for Scan® 300** and **Scan® 500**. **Scan®** colony counter automatically separates confluent colonies, allows you to create polygonal exclusion areas and ignores agar flaws and air bubbles. You can also add or remove colonies manually. Every change is automatically saved in your report.



Scan® read all the colonies, even the smallest




Automatic separation of confluent colonies



Automatic elimination of counting grids



Cross on each counted colony



Polygonal exclusion areas




Automated software with manual control

Color detection & chromogenic media

Scan® 500 and **Scan® 1200** can read chromogenic agar and differentiate colonies by color: **up to 7 different colors on the same Petri dish**. Color selection can be made directly from the color of the bacteria and a cursor allows you to set the sensitivity.

Chromogenic media reading allows the detection of Salmonella on XLD media and E.Coli on TBX media, for example.



INCLUDED IN THE SOFTWARE*

* on Scan® 500 & Scan® 1200

Scan[®] : 3 models adapted to your needs

Scan[®] 300 Essential

- CCD color camera, zoom x7
- 6 combinations of lighting and backgrounds
- Brightness, contrast and sensitivity are automatically optimized by the software
- Minimum size of detected colony: 0.1 mm
- Automatic separation of clustered colonies
- Creation of polygonal exclusion zones
- Long lasting LED lighting
- Export of data to Excel[™], PDF
- Software in French, English, Japanese, Chinese, Russian, Spanish

Count these supports



Petri Dish



Spiral[®] plating



Ref 436 300

Scan[®] 500 Efficient

- Color CMOS camera, zoom x7
- 6 combinations of lighting and backgrounds
- Brightness, contrast and sensitivity are automatically optimized by the software
- Minimum size of detected colony: 0.1 mm
- Automatic separation of clustered colonies
- Creation of polygonal exclusion zones
- Long lasting LED lighting
- Export of data to Excel[™], PDF
- Software in French, English, Japanese, Chinese, Russian, Spanish
- **Detects and counts up to 7 colors on the same dish**
- **Adapts itself to most kinds of agar**

Count these supports



Petri Dish



Spiral[®] plating



Chromogenic
Petri Dish



Inhibition
zone



Ref 436 000

Scan[®] 1200 High-Resolution

- 6 combinations of lighting and backgrounds
- Detects and counts up to 7 colors on the same dish
- Brightness, contrast and sensitivity are automatically optimized by the software
- Minimum size of detected colony: 0.05 mm
- Automatic separation of clustered colonies
- Creation of polygonal exclusion areas
- Long lasting LED lighting
- Export of data to Excel[™], PDF
- Software in French, English, Japanese, Chinese, Russian, Spanish
- **HD color CCD camera, zoom x28**
- **Adapts itself to most kinds of agar**
- **Motorized background color**
- **Lights controlled by Scan[®] Software**

Count these supports



Petri Dish



Spiral[®] plating



Chromogenic
Petri Dish



Inhibition zone



Petrifilm[™]



Compact Dry[™]



RIDA[™] Count
Sanita-kun[™]



Filtration
membrane



Ref 437 000

Inhibition zone from paper discs, agar wells & peni cylinders

Performance and flexibility

Scan® 500 and **Scan® 1200** allow efficient work flow because you can create and edit a list of antibiotics, useful for routine analysis.

Measured by **Scan®**, inhibition zones guarantee repeatability and reproducibility of analysis and diagnosis reliability.

- Rapid detection: up to 8 antibiotic sensitivities in 1 click.
- Paper discs, agar wells and peni cylinders may be manually added or deleted. Inhibition zones may be manually resized.




The result of sensitivity in contact with the antibiotic is fast and visualization of results is clear:

- > Red (resistant)
- > Yellow (intermediate)
- > Green (susceptible)

Medical analysis


Inhibition zone measurement allows you to test the efficiency of antibiotics on micro-organisms to accelerate the diagnosis in order to choose precisely an appropriate antibiotic treatment for a patient. **Scan®** has a built-in antibiotic database from the **French Society of Microbiology (CA-SFM)** and of the **European Committee on Antimicrobial Susceptibility Testing (EUCAST)** which determines the sensitivity of the bacteria to the antibiotic. This database is fully editable.



Precision of inhibition zone radius measurement from paper discs : 0.1 mm

Pharmaceutical analysis

In the pharmaceutical industry, **Scan®** allows you to test the quality of an antibiotic during its manufacturing process by measuring the inhibition zones. To evaluate the action of an antibiotic, antibiotic diffusion from paper disc, agar well or peni cylinder is supported.



Precision of inhibition zone radius measurement from agar wells : 0.1 mm



* on Scan® 500 & Scan® 1200

Comfort of use

> High definition live image


This feature enables total control of colony counting.

Optimum visualization

Enjoy comfortable viewing of the colonies with the unequalled **Dark Field technology** (double oblique tangent and crossed light), high definition live image and with the automatic optimization of the image (lighting, contrast and sensitivity). You can also check key areas thanks to the digital zoom.



Dark Field: LED are disposed in a circle for optimal contrast



Scan® automatically optimizes contrast, luminosity and sensitivity



Digital zoom with the mouse wheel (up to x28)




Image with Scan® 1200

Easy-to-use

All **Scan®** functions are in **one single window** and colonies are counted in one click.

The **Scan®** easy commands (visualization, settings and results) allow quick access to both ongoing and archived work sessions.

Scan® software is available in **6 languages** (English, French, Chinese, Russian, Japanese and Spanish) and is updated regularly. The intuitive use of **Scan®** **does not require any special training.**



Fast communication, total traceability

Results harmonization

Using the **Scan**[®] allows more reliable analyses and harmonizes the results within a team.

You can save as many settings as you wish and customize the settings according to the type of dishes and agar you use.

The automatic archiving of data, photos, comments and results ensures total traceability.

Print your results

You can export your results to your PC, archive it in Excel[™], PDF or SCA format. You can also export pictures from the camera in JPEG, PNG and BMP format.



Printed report example

Analysis result

Comments

Sample information

Add your own logo in the reports

Petri dish before counting

Petri dish after counting

The screenshot shows a printed report with the following sections:


- Logo:** A vertical green box with the word "interscience" in white.
- Petri dish before counting:** A circular image of a petri dish with a blue agar surface and numerous small white colonies.
- Petri dish after counting:** A circular image of a petri dish with a blue agar surface and colonies that are color-coded in green and orange. A small purple box in the top right corner of the image indicates a diameter of "Φ 80.8 mm".
- Sample information table:**

Sample N°:	1	Count :	915	Dilution :	1/10	CFU/mL :	1.14E+004
Parameters :	GL Lyon	Area (%) :	81 %	Sensitivity :	37		
Date Time :	15/05/2009 15:47:28						
- Comments:** A text box containing "Comment : 1 CFU added manually --".

Reproducibility of results

Automatic counting is a guarantee of **regularity and standardization** of analyses, which is the key to ensure accurate and reliable results. **Reproducibility** of results is guaranteed whatever the day, conditions and user.

A scientific study has proved 98% precision for **Scan®** colony counters. This study is available on request.



Manual counting:

Random results over time and different users by manual counting of colonies

Automatic counting:


Standardized and reproducible results by automatic counting

Internal traceability


Thanks to the **LIMS** connection and the barcode reader, photos of counted plates are **saved and traceable**. The images are accessible and recountable at any time.



Barcode




Archiving



LIMS connection


Secure your sessions

Sessions are secured with a **security code** (one per operator) and the impossibility to alter each saved counting. **Scan®** use allows the compliance with **CFR 21 part 11**: system securization, operational controls and documentation management.



> External traceability

Scan® software provides numerous possibilities to easily and quickly export your results.



Work sessions saving




PDF export



JPEG, PNG & BMP formats export



Export results to Excel™ for ensure traceability



Print report from Scan®

Plate & Count System[®] with dataLink[™]

INCREASE
your lab capacity with
easySpiral[®] & Scan[®]!



More info on
www.interscience.fr

Plate & Count System[®] enables automatic plating and colony counting with full traceability!

- **INCREDIBLE SAVINGS:** Save up to 75% in time, consumables and bench space
- **FAST:** Full plating cycle in 25 seconds and counting in 1 click
- **RELIABLE:** 98% repeatable and reproducible results
- **TOTAL TRACEABILITY WITH dataLink[™]:** Automatic data saving and reporting

Plate & Count System® includes:

- **easySpiral®**: Automatic Spiral platers
- **Scan®**: Automatic colony counters



Plate and count your Petri dishes: UP TO 75% SAVINGS

Once the Spiral® plated dish is out of the incubator, it is ready to be counted by **Scan®** automatic colony counters. Results are immediately displayed and saved.

easySpiral® and **Scan®** guarantee the regularity and standardization of the analyses, save time, consumables and bench space of up to 75%.

Plating with easySpiral®: FROM 30 TO 10⁷ CFU/ML ON ONE SINGLE PETRI DISH

In compliance with **AFNOR V08-100**, **ISO 4833-2** and **ISO 7218** standards, **easySpiral®** automatically plates a sample in **8 seconds**: from 30 to 10⁷ CFU/mL on a **single Petri dish** without prior sample dilution. Once the sample is plated and incubated, you can analyze it manually or automatically by counting all or part of its colonies.

ISO 4833-2	ISO 7218	AFNOR V08-100
----------------------	--------------------	-------------------------

Manual plating method

This method requires repetitive actions: at least **four dilutions** and **four successive platings** are necessary to obtain one good and readable Petri dish.

Automatic Spiral® method

With this method, make your analysis on **1 single Petri dish!**

All dilutions on 1 Petri dish!

- From 30 to 10 million CFU/mL on 1 single Petri dish
- Up to 75% less consumables
- Full plating cycle in 25 seconds!

Technical specifications

		Scan® 300	Scan® 500	Scan® 1200
	Reference	436 300	436 000	437 000
Image	Camera	CCD color camera	CMOS color camera	HD CCD color camera
	Digital zoom	x 7		x 28
	Resolution	640 x 480 pixels		1280 x 960 pixels
	Counting time	1000 colonies per second		
	Minimum size of colonies	0.1 mm		0.05 mm
	Lighting technology	Long-life white LEDs / Dark Field		
	Lighting system	6 combinations, top and/or bottom light white or black background		Automatic 6 combinations, top and/or bottom light white or black background
Counting	Counting	Automatic, with manual control		
	Results / data export	Scan® recountable file, PDF report, JPEG, PNG, BMP, Excel™		
	Color detection	-	7 colors on the same Petri dish	
	Chromogenic medium	-	✓	✓
	Inhibition zones	-	✓	✓
	Petri dishes	✓ (55-90 mm)	✓ (55-90 mm)	✓ (55-90 mm)
	Spiral® plating	✓	✓	✓
	Petrifilm™	-	-	✓
	RIDA™ Count	-	-	✓
	Compact Dry™	-	-	✓
	Filtration membrane	-	-	✓
	LIMS connection	✓	✓	✓
	Languages	English, French, Japanese, Chinese, Russian, Spanish		
Hardware	Dimensions (L x W x H)	27 x 27 x 29 cm		27 x 27 x 38 cm
	Weight	6.6 kg		8.7 kg
	Body	Stainless steel		
	Computer connection	Firewire A		
	Power	100-240 V~ 50/60 Hz		
PC requirements	Operating systems	Windows™: Vista, 7, 8		
	Processor	Intel (recommended) or others, 2 GHz or better		Intel Core (recommended) or AMD Phenom and superior, 2 GHz
	RAM	2 Go		3 Go
	Equipment	FireWire connector or free PCI slot / CD ROM drive		
	Screen	1280 X 1024 pixels and more		
	Computer	Desktop computer recommended / Laptop not recommended		

Delivered with : Scan® software CD-ROM/ user guide/ FireWire adapters: PCI express, Expresscard, adapter 6/4 pins



Guarantee

3 years Guarantee / 3 years free Software update (after registration of the guarantee form)

PC requirements are subject to change. Please check our website www.interscience.fr for current updates and additional informations.

Scan® accessories



Barcode reader
Ref. : 522 000



Petrifilm™ adapter
Ref. : 436 002



RIDA™ Count /
Sanita-kun™ adapter
Ref. : 436 001




Compact Dry™ adapter
Ref. : 436 004



55 mm dishes adapter
Ref. : 436 005

Discover our complete range for microbiology



www.interscience.fr



DOWNLOAD
SCAN® SOFTWARE

www.interscience.fr

On our website, find and download:

- Latest information on our products
- Demonstration photos and videos



interscience INTERNATIONAL

30, ch. Bois Arpents - 78860 St Nom - FRANCE
T: +33 (0)1 34 62 62 61 - F: +33 (0)1 34 62 43 03
info@interscience.fr - www.interscience.fr

interscience USA & CANADA

199 Weymouth ST. - ROCKLAND - MA 02370 - USA
T: +1 781 792 2133 - F: +1 781 792 2134
info@intersciencelab.com - www.intersciencelab.com

interscience CHINA

上海市徐汇区吴兴路277号锦都大厦718室-200030
电话: +86 (0)21-64739390 - 手机: +86 13040616422
sales.china@interscience.fr - www.interscience.fr/china

Distributor