① Conexión BNC.
Connessione BNC.
Connexion BNC.
BNC connection.

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Ag electrode 50 46

Introduction

The **50 46** is a silver electrode for argentometric titrations.

Electrode preparation

Remove the storage protector containing electrolyte by unscrewing it ⑤.

Remove the cap from the filling port ③.

Measurements. Recommendations

- The sample must cover the electrode's diaphragms ⑦.
- After measurement, clean the electrode with distilled water.

Maintenance

- Periodically refill the electrode with electrolyte. Maintain the level as high as possible ④.
- Diaphragm cleaning: use the diaphragm cleaning solution, code 96 10.
- Silver ring cleaning: can be ribbed with soft abrasive, for example tooth paste.

Troubleshooting

... the reading is always the same.

Short circuit in the electrode or connector.

Substitute the electrode.

... unstable reading.

Broken cable.

Substitute the cable.

... very slow response.

Dirt on the diaphragm or the silver ring.

Reference system worn out.

Clean, regenerate or substitute the electrode.

Storage

Place the electrode in the storage protector ⑤ with appropriate electrolyte KNO_3 1 M.

Close the refilling port with the stopper ③.

Life of an electrode

The life expectancy of the electrode could be indefinite when the suitable maintenance is applied.

Specifications

Measuring range	± 2000 mV
Indicator metal	silver
Operating temp.	0 ... 80 °C
Reference	encapsulated Ag/AgCl
Diaphragm	ceramic
Electrolyte	KNO_3 1M
Body material	glass
Min. puncture depth	12 mm

Applications:

Argentometric titrations.

Limits

Products with colloids or solids in suspension.

⚠ Warning.

It uses a special electrolyte, KNO_3 1 M, which prevents clogging of the diaphragm caused by AgCl precipitation.

Warranty

CRISON guarantees this electrode against manufacturing defects.

Validity: 6 months from shipment date.

Limitations:

- Accidental damage.
- Inadequate applications.
- Non-fulfilment of the recommendations.
- Normal wear and tear of the electrode.

QUALITY CERTIFICATE

The attached electrode **50 46** has passed the quality test and follows the next specifications:

- Deviation $< \pm 10$ mV in AgNO_3 standard.
- Response time in AgNO_3 standard < 20 s.

CRISON INSTRUMENTS, S.A.

Riera Principal 34 - 36
08328 Alella BARCELONA
SPAIN

Service

Tel.: +34 935 550 318
Fax: +34 935 400 857
E-mail: service@crison.es

Tel.: +39 059 651 922
Fax: +39 059 652 011
E-mail: crison@crison.it

Tel.: +33 0555 639 726
Fax: +33 0555 639 727
E-mail: info@crison.fr

www.crisoninstruments.com