50 51 T

 **CRISON**

50 51 T



Fig. A

Puntos clave. *Punti chiavi.*

Points clé. *Key points.*

1

Conector: MP-5, BNC, DIN, etc.
Connettore: MP-5, BNC, DIN, ecc.
Connecteur: MP-5, BNC, DIN, etc.
Connector: MP-5, BNC, DIN, etc.

2

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Introduction

The **50 51 T** is a pH low maintenance electrode, as it does not require monitoring of electrolyte level, with built-in A.T.C. Pt 1000. The open diaphragm allows direct contact between the sample and the reference electrolyte. This avoids the classical problem of diaphragm clogging.

Electrode preparation

Remove the storage protector containing electrolyte by unscrewing it and put the working protector on its place (see Fig. A).

⚠ *MP-5 connector (IP 67). Connect and disconnect by pushing. Do not unscrew it.*

Calibration

⚠ *See pH-meter instruction manual.*

Measurements. Recommendations

- During calibration and measurement the buffers and the samples must be slightly stirred. The stirring speed should be identical in both cases. In field measurements, these types of electrodes allow stirring the sample directly holding the electrode by the handle.
- The sample must cover the electrode's diaphragm ⑦.
- After measurement, immediately rinse the electrode with distilled water. Dry it with soft tissue, only by contact, without rubbing.
- Between measurements the electrode's membrane should be immersed in aqueous solution. The most appropriate place is the storage protector with some electrolyte.

⚠ WARNING

- *Do not leave the electrode immersed in distilled water.*
- *Do not leave the electrode immersed in solutions with pH < 2.*
- *Do not bend the cable on the handle level.*

Maintenance

- Maintain the membrane's surface ⑧ clean.
- When an electrode gets contaminated, the simply rinsing with distilled water is not sufficient. The most suitable cleaning solution is which acts in a most selective way over the contamination.

Proteins: Immerse the electrode in 2% NaOH solution during approx. 20 min.

Fats, oils: Wash the electrode with ethanol or soapy tap water at approx. 40°C.

Afterwards, rinse the electrode with plenty of distilled water and leave for rest in the reference electrolyte.

Life of an electrode

The average life expectancy of a pH electrode is approximately one year depending on the working conditions.

An electrode used for more than six months that does not response appropriately after regeneration, must be replaced by a new one.

⚠ *If an electrode has too short life span, less than three months, ask us about it. Probably there are more suitable electrode for the application.*

Storage

Place the electrode in the storage protector ④ with appropriate electrolyte (Crisolyt G, (KCl with glycerine).

Troubleshooting

- ... the pH-meter always measures the same pH value.
Short circuit in the electrode, cable or connector.
Substitute the electrode.
- ... fortuitous reading.
Broken cable.
Substitute the electrode.
- ... impossible to perform calibration.
- ... very slow response.
Dirt on the diaphragm and/or the membrane.
De-hydrated membrane.
Electrode ageing or wear.
Clean, regenerate or replace the electrode.

Specifications

pH measuring range	2 ... 14
Operating temp.	0 ... 80 °C
Reference element	sleeved Ag wire coated with AgCl
Diaphragms	open, circular
Electrolyte	solid
Body material	glass
Min. immersion depth	15 mm 

Applications

Food industry, viscous media, dirty samples, etc.

Limitations

- Samples with pH < 2.
- Temperatures exceeding 80 °C.
- Samples with low conductivity, distilled water.

Warranty

CRISON guarantees this electrode against manufacturing defects.

Validity: 6 months from shipment date.

Limitations:

- Accidental damage.
- Inadequate applications.
- Non-fulfilment of the recommendations.
- Normal wear and tear of the electrode.

QUALITY CERTIFICATE

The attached electrode **50 51 T** has passed the quality test and follows the next specifications:

- Asymmetry potential < ± 15 mV.
- Sensitivity, pH 4...7 (at 25 °C) > 98%.
- Response time, pH 4...7 < 20 s
- Temperature error $\leq 0.5^\circ\text{C}$

CRISON INSTRUMENTS, S.A.

Riera Principal 34 - 36
08328 Alella BARCELONA
SPAIN

Service

Tel.: +34 935 550 318
Fax: +34 935 400 857
E-mail: service@crison.es

Tel.: +39 059 651 922
Fax: +39 059 652 011
E-mail: crison@crison.it

www.crisoninstruments.com