

52 60

mV

CRISON

Puntos clave. *Punti chiave.* Points clé. *Key points.*

① Tapón protector.
Tappo protettore.
Bouchon protecteur.
Protective cap.

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Ag electrode 52 60

Introduction

The **52 60** is a silver electrode for argentometric titrations.

Electrode preparation

Remove the storage protector containing electrolyte by unscrewing it (5).

Remove the cap from the filling port (3).

Measurements. Recommendations

- The sample must cover the electrode's diaphragms (7).
- After measurement, clean the electrode with distilled water.

Maintenance

- Periodically refill the electrode with electrolyte. Maintain the level as high as possible (4).
- Diaphragm cleaning: use the diaphragm cleaning solution, code 96 10.
- Indicator metal cleaning: can be ribbed with soft abrasive, for example tooth paste.

Troubleshooting

... the reading is always the same.

Short circuit in the electrode or connector.

Substitute the electrode.

... unstable reading.

Broken cable.

Substitute the cable.

... very slow response.

Dirt on the diaphragm or the silver ring.

Reference system worn out.

Clean, regenerate or substitute the electrode.

Storage

The connector should be protected from humidity by the protective cap (1).

Place the electrode in the storage protector (5) with appropriate electrolyte KNO_3 1M.

Close the refilling port with the stopper (3).

Life of an electrode

The life expectancy of the electrode could be indefinite when the suitable maintenance is applied.

Specifications

Measuring range	± 2000 mV
Indicator metal	silver
Operating temp.	0 ... 80 °C
Reference	encapsulated Ag/AgCl
Diaphragm	ceramic
Electrolyte	KNO_3 1M
Body material	glass
Min. puncture depth	12 mm 

Applications:

Argentometric titrations.

Limits

Products with colloids or solids in suspension.

⚠ Warning.

It uses a special electrolyte, KNO_3 1M, which prevents clogging of the diaphragm caused by $AgCl$ precipitation.

Warranty

CRISON guarantees this electrode against manufacturing defects.

Validity: 6 months from shipment date.

Limitations:

- Accidental damage.
- Inadequate applications.
- Non-fulfilment of the recommendations.
- Normal wear and tear of the electrode.

QUALITY CERTIFICATE

The attached electrode 52 60 has passed the quality test and follows the next specifications:

- Deviation $< \pm 10$ mV in $AgNO_3$ standard.
- Response time in $AgNO_3$ standard < 20 s.

CRISON INSTRUMENTS, S.A.

Riera Principal 34 - 36
08328 Alella BARCELONA
SPAIN

Service

Tel.: +34 935 550 318
Fax: +34 935 400 857
E-mail: service@crison.es

Tel.: +39 059 651 922
Fax: +39 059 652 011
E-mail: crison@crison.it

www.crisoninstruments.com