

*insta***BAG**[®]

Filter blender bag with dehydrated media

Your diluent ready in minutes!



interscience

Ready-to-use media for easy preparation

instaBAG® is a filter bag which includes a dehydrated pre-dosed media, allowing a 2-in-1 process. Skip the time-consuming media preparation.



1 Place the sample in the *instaBAG*®



2 Add sterile and deionized water



3 Blend the *instaBAG*® with the BagMixer®

➔ The sample is ready to be pipetted

All the advantages of a filter bag

Debris are blocked by the filter during blending

Long shelf life

Easy pouch opening

Zip pouch with desiccant: optimal conservation after opening

Light and humidity proof aluminum pouch



Just add water

Ready-to-use bag with pre-dosed dehydrated media

Perfect dissolution during blending

Non-woven filter identical to the bag filter

Available in BPW 90 mL and BPW 225 mL



PATENTS
PENDING

SHELF LIFE
12 MONTHS

SRT
6°C to 30°C

* Store at Room Temperature

Space saving




100
instaBAG



100
glass bottles

No waste

100 %
used
3 *instaBAG*
for 3 analyses.



13 %
used
1 (5L) pouch
for 3 analyses.
The leftover liquid
is thrown away.

Less garbage



instaBAG
+ diluent + sample




Media bottle
+ bag + diluent + sample

Easy preparation



with *instaBAG*.



with media preparation

Green

instaBAG : no water transported



Diluent pouches




Long shelf life

instaBAG



Media bottle (BPW)



Technical specifications

	instaBAG® BPW 90	instaBAG® BPW 225 mL
Reference	114 090	114 225
Type of bag	Bag with lateral filter	
Composition of the bag	MultiLayer®: reinforced multicoated complex	
Type of filter	Non-woven	
Filter porosity	< 250 microns	
Type of dehydrated media	Buffered peptone water	
Weight of dehydrated media	1.8 g	4.5 g
Sample weight	10 g (± 5 %)	25 g (± 5 %)
Volume of water to be added (sterile and deionized)	90 mL (± 2 %)	225 mL (± 2 %)
Blending time	1 mn	
Bag dimensions	190 x 300 mm	
Aluminum pouch of	10 bags	
Shelf life	12 months	
Storage conditions	Temperature from + 6°C to + 30°C	
Traitment	Gamma ray (10 to 25 KGy)	
Standards	ISO 11133 - ISO 11290-2 - ISO 215281 - ISO 7218 - ISO 6579 - ISO 6887- FDA-BAM	
Box of	100 bags	

MSDS available on request

Ask for free samples

Certified production

Products made for INTERSCIENCE by Interlab, an ISO 9001 certified company.



The winning combination



DiluFlow® Gravimetric dilutors



instaBAG®



BagMixer® Lab blenders

Your local distributor

interscience

30, chemin du Bois des Arpents - 78860 St Nom - FRANCE
T: +33 (0)1 34 62 62 61 - Email: info@interscience.com

interscience USA & CANADA

32 Cummings Park - Woburn, MA 01801 - USA
P: +1 781 937 0007 - F: +1 781 937 0017 - Email: sales.usa@interscience.com

interscience CHINA

上海市徐汇区吴兴路277号锦都大厦588室 - 200030
电话: +86 (0)21-64739390 - +86 189 3097 0733 - 邮址: sales.china@interscience.com

interscience SOUTH-EAST ASIA

The Centropod, 80 Changi Rd - 05-07 Singapore 419715
T: +65 6909 0825 - M: +65 8118 5178 - E-mail: sales.asia@interscience.com

www.interscience.com