

52 11

pH

 **CRISON**

# Puntos clave. *Punti chiave.* Points clé. *Key points.*

① Tapón protector.  
*Tappo protettore.*  
Bouchon protecteur.  
*Protective cup.*

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# pH electrode 52 11

## Introduction

The **52 11** is a pH electrode with screw cap S7, suitable for samples containing proteins. Large, semi-spheric membrane, very robust and fast responding.

## Electrode preparation

Remove the storage protector containing electrolyte by unscrewing it (5).

Verify that there are not air bubbles inside the membrane. They can be removed by shaking the electrode as a clinical thermometer. Remove the cup from the electrolyte filling port (3).

## Calibration

**⚠** See pH-meter instruction manual.

## Measurements. Recommendations

- During calibration and measurement the buffers and the samples must be slightly stirred. The stirring speed should be identical in both cases.
- The sample must cover the electrode's diaphragms (7).
- After measurement, immediately rinse the electrode with distilled water. Dry it with soft tissue, only by contact, without rubbing.
- Between measurements the electrode's membrane should be immersed in aqueous solution. The most appropriate place is the storage protector with some electrolyte.

**⚠** Do not leave the electrode immersed in distilled water.

## Maintenance

- Periodically refill the electrode with electrolyte. Maintain the level as high as possible (4).
- Maintain the membrane's surface clean (8). Protect it from friction or knocking.
- When an electrode gets contaminated, the simple rinsing with distilled water is not sufficient. The most suitable cleaning solution is one which acts in a more selective way over the contamination (look for CRISON regeneration solutions).
- **Re-hydration.** A "dry" electrode recovers by immersing it in CRISON cleaning electrode solution, code 96 00, or in diluted HCl solution for several hours. It is recommended to read the information supplied with the electrode's regenerating solutions.

## Troubleshooting

... the pH-meter always measures the same pH value.

Short circuit in the electrode or connector.

**Substitute the electrode.**

... unstable reading.

Air bubble in the membrane.

**Shake the electrode as a thermometer.**

... impossible to perform calibration.

... very slow response.

Dirt on the diaphragm or the membrane.

De-hydrated membrane.

Electrode ageing or wear.

**Clean, regenerate or substitute the electrode.**

## Storage

The connector should be protected from humidity by the protective cap ①.

Place the electrode in the storage protector ⑤ with Crisolyt G (KCl + glycerine).

Close the refilling port with the stopper ③.

## Life of an electrode

The average life expectancy of a pH electrode is approximately one year depending on the working conditions.

An electrode used for more than six months that does not respond appropriately after regeneration, must be replaced by a new one.

To avoid accidental breakage of the glass body electrodes, the use of the working protector, code 91 62, is recommended.

**⚠** *If an electrode has too short lifespan, less than three months, ask us about it. Probably there is a more suitable electrode for the application.*

## Specifications

Measuring range	0... 12
Operating temp.	0... 100 °C
Reference element	encapsulated Ag/AgCl crystals
Diaphragms	3 ceramic
Electrolyte	solid
Body material	glass
Min. immersion depth	20 mm

## Applications

It is appropriate for samples with proteins and also for high temperatures up to 100 °C.

## Limitations

Very viscous solutions or products with colloids able to clog the three diaphragms.

## Warranty

CRISON guarantees this electrode against manufacturing defects.

**Validity:** 6 months from shipment date.

**Limitations:**

- Accidental damage.
- Inadequate applications.
- Non-fulfilment of the recommendations.
- Normal wear and tear of the electrode.

## QUALITY CERTIFICATE

The attached electrode 52 11 has passed the quality test and follows the next specifications:

- Asymmetry potential  $< \pm 15$  mV.
- Sensitivity, pH 4...7 (at 25 °C)  $> 98\%$ .
- Response time, pH 4...7  $< 20$  s.

**CRISON INSTRUMENTS, S.A.**  
Riera Principal 34 - 36  
08328 Alella - BARCELONA  
SPAIN

**Service**

Tel.: +34 935 550 318  
Fax: +34 935 400 857  
E-mail: [service@crison.es](mailto:service@crison.es)

Tel.: +39 059 651 922  
Fax: +39 059 652 011  
E-mail: [crison@crison.it](mailto:crison@crison.it)

[www.crisoninstruments.com](http://www.crisoninstruments.com)