
52 33

pH

**CRISON**

# 52 33

# Puntos clave. *Punti chiave.* Points clé. *Key points.*

1 Tapón protector.  
*Tappo protettore.*  
Bouchon protecteur.  
*Protective cup.*

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# pH electrode 52 33

## Introduction

The **52 33** is a puncture electrode with plastic body (POM). Suitable for measurements in semisolid samples.

## Electrode preparation

Remove, with care, the storage protector containing electrolyte by unscrewing it (6).

## Calibration

pH electrodes must be calibrated with buffer solutions before measuring.

**⚠** See pH-meter instruction manual.

## Measurements. Recommendations

- Do not introduce the electrode directly in thick samples or samples containing solid materials, bones, etc. It is necessary to open a path using a stainless steel perforator, introduce the electrode and take the measurement.
- During measurements, hold the measured product still.
- The sample must cover the electrode's diaphragms (4).
- After measurement, clean the electrode with cotton soaked in alcohol.
- Between measurements the electrode's membrane should be immersed in aqueous solution. The most appropriate place is the storage protector with some electrolyte.
- This electrode allows measure with the membrane in any position.

## **⚠** WARNING

- Do not leave the electrode immersed in distilled water.
- Do not bend the electrode. The internal body is made from glass.

## Maintenance

- Maintain the membrane's surface (5) clean. Do not rub it. Avoid knocks and scratches.
- After measurement, clean the electrode with cotton soaked in alcohol.
- At the end, clean with lukewarm soapy water.

**⚠** Never use cleaning solutions with  $\text{pH} < 2$ . It provokes alteration in the reference electrolyte and damages irreversibly the electrode.

## Troubleshooting

... the pH-meter always measures the same pH value.

Short circuit in the electrode or connector.

*Substitute the electrode.*

... unstable reading.

Broken cable.

*Substitute the cable.*

... impossible to perform calibration.

... very slow response.

Dirt on the diaphragm or the membrane.

De-hydrated membrane.

Electrode ageing or wear.

*Clean, regenerate or substitute the electrode.*

## Storage

The connector should be protected from humidity by the protective cap ①.

Place the electrode in the storage protector ⑥ with appropriate electrolyte (Crisolyt G, KCl + glycerine).

## Life of an electrode

The average life expectancy of a pH electrode is approximately one year depending on the working conditions.

An electrode used for more than six months that does not respond appropriately after regeneration, must be replaced by a new one.

 *If an electrode has too short life span, less than three months, ask us about it.*

## Specifications

Measuring range	2 ... 14 pH
Operating temp.	0 ... 60 °C
Reference	encapsulated Ag/AgCl
Diaphragms	1 ceramic + 1 open
Electrolyte	solid
Body material	glass / plastic
Min. puncture depth	12 mm

## Applications:

- Cheeses, meats, fish, fruits, bread doughs, etc. It can be used for measurements in aqueous solutions too.

## Limits

pH values < 2.

Temperatures > 60 °C.

## Warranty

CRISON guarantees this electrode against manufacturing defects.

**Validity:** 6 months from shipment date.

### Limitations:

- Accidental damage.
- Inadequate applications.
- Non-fulfilment of the recommendations.
- Normal wear and tear of the electrode.

## QUALITY CERTIFICATE

The attached electrode 52 33 has passed the quality test and follows the next specifications:

- Asymmetry potential <  $\pm 15$  mV.
- Sensitivity, pH 4...7 (at 25 °C) > 98%.
- Response time, pH 4...7 < 20 s.

**CRISON INSTRUMENTS, S.A.**

Riera Principal 34 - 36  
08328 Alella - BARCELONA  
SPAIN

**Service**

Tel.: +34 935 550 318  
Fax: +34 935 400 857  
E-mail: [service@crison.es](mailto:service@crison.es)

Tel.: +39 059 651 922  
Fax: +39 059 652 011  
E-mail: [crison@crison.it](mailto:crison@crison.it)

[www.crisoninstruments.com](http://www.crisoninstruments.com)