50 52 T

 **CRISON**

# 50 52 T



Fig. A

## Puntos clave. *Punti chiavi.* Points clé. *Key points.*

1 Conector MP-5.  
*Connettore MP-5.*  
Connecteur MP-5.  
*Connector MP-5.*

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# pH electrode 50 52 T

## Introduction

The **50 52 T** is a pH electrode, with large ring diaphragm of PTFE and built-in A.T.C. Pt 1000, for use with CRISON portable PH 25 and PH 25 DL. These are low maintenance electrodes, as they do not require monitoring of electrolyte levels. This model is supplied with a "working protector", protecting the electrode against knocks.

## Electrode preparation

Remove the storage protector containing electrolyte by unscrewing it and put the working protector on its place (see Fig. A).

**⚠** *MP-5 connector (IP 67). Connect and disconnect by pushing. Do not unscrew it.*

## Calibration

All pH electrodes must be calibrated with buffer solutions before measuring.

**⚠** *See pH-meter instruction manual.*

## Measurements. Recommendations

- During calibration and measurement the buffers and the samples must be slightly stirred. The stirring speed should be identical in both cases. In field measurements, these types of electrodes allow stirring the sample directly holding the electrode by the handle.
- The sample must cover the electrode's diaphragm ⑦.
- After measurement, immediately rinse the electrode with distilled water. Dry it with soft tissue, only by contact, without rubbing.
- Between measurements the electrode's membrane should be immersed in aqueous solution. The most appropriate place is the storage protector with some electrolyte.

## **⚠** WARNING

- *Do not leave the electrode immersed in distilled water.*
- *Do not bend the cable on the handle level.*
- *Never unscrew the MP-5 connector.*

## Maintenance

- Maintain the membrane's surface ⑨ clean.
- When an electrode gets contaminated, the simply rinsing with distilled water is not sufficient. The most suitable cleaning solution is which acts in a most selective way over the contamination (see CRISON regeneration solutions).
- **Re-hydration.** A "dry" electrode recovers by immersing it in CRISON cleaning electrode solution, code 96 00, or in diluted HCl solution for several hours. It is recommended to read the information supplied with the electrodes regenerated solutions.

## Life of an electrode

The average life expectancy of a pH electrode is approximately one year depending on the working conditions.

An electrode used for more than six months that does not response appropriately after regeneration, must be replaced by a new one.

**⚠** *If an electrode has too short life span, less than three months, ask us about it. Probably there are more suitable electrode for the application.*

## Storage

Place the electrode in the storage protector ④ with appropriate electrolyte (Crisolyt, KCl 3M).

## Troubleshooting

- ... the pH-meter always measures the same pH value.  
Short circuit in the electrode, cable or connector.  
**Substitute the electrode.**
- ... fortuitous reading.  
Broken cable.  
**Substitute the electrode.**
- ... impossible to perform calibration.
- ... very slow response.  
Dirt on the diaphragm and/or the membrane.  
De-hydrated membrane.  
Electrode ageing or wear.  
**Clean, regenerate or replace the electrode.**

## Specifications

pH measuring range	0 ... 14
Operating temp.	0 ... 100 °C
Reference element	sleeved Ag wire coated with AgCl
Diaphragms	ring of porous PTFE
Electrolyte	gel
Body material	glass
Temperature sensor	Pt 1000
Min. immersion depth	12 mm

## Applications:

- Measurements in "difficult" samples, creams, dirty or viscous samples, etc.
- Measurements at high temperature.
- Distilled water.

## Warranty

CRISON guarantees this electrode against manufacturing defects.

**Validity:** 6 months from shipment date.

## Limitations:

- Accidental damage.
- Inadequate applications.
- Non-fulfilment of the recommendations.
- Normal wear and tear of the electrode.

## QUALITY CERTIFICATE

The attached electrode **50 52 T** has passed the quality test and follows the next specifications:

- Asymmetry potential  $< \pm 15$  mV.
- Sensitivity, pH 4...7 (at 25 °C)  $> 98\%$ .
- Response time, pH 4...7  $< 20$  s
- Temperature error  $\leq 0.5^\circ\text{C}$

**CRISON INSTRUMENTS, S.A.**

Riera Principal 34 - 36  
08328 Alella BARCELONA  
SPAIN

**Service**

Tel.: +34 935 550 318  
Fax: +34 935 400 857  
E-mail: [service@crison.es](mailto:service@crison.es)

Tel.: +39 059 651 922  
Fax: +39 059 652 011  
E-mail: [crison@crison.it](mailto:crison@crison.it)

[www.crisoninstruments.com](http://www.crisoninstruments.com)