COW Milk Express

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**Latex test for the detection of raw and heated
cow's milk in milk from other species**

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Development | Manufacture | R&D and OEM Services

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**This document must be fully read before running the kit
for the first time.**

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TABLE OF CONTENTS

1.	FIELD OF APPLICATION	9
2.	ASSAY PRINCIPLE	9
3.	KIT CONTENT	9
4.	STORAGE CONDITIONS / SHELF LIFE	10
5.	STORAGE AND PREPARATION OF THE SAMPLES.....	10
6.	INSTRUCTIONS	10
7.	READING AND INTERPRETATION.....	11
8.	TECHNICAL SPECIFICATIONS	12
9.	SAFETY RECOMMENDATIONS.....	12
	REFERENCE	13

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1. FIELD OF APPLICATION

The Cow Milk Express is a latex test for the rapid detection of cow's milk in milk of other species (except buffalo milk).

The test is based upon the specific detection of bovine beta-lactoglobulin (Blg) in the sample. Since the Blg is resistant to thermal denaturation, the test can be applied for the detection of raw and heated (pasteurized and sterilized) cow's milk. The test can also be used to detect the presence of powder of bovine milk or whey in the sample.

The limits of detection of the test are 0,25% of raw milk and 1% of sterilized milk in only 2 minutes.

The Cow Milk Express can be used with whole milk, semi-skimmed, skimmed milk and whey.

2. ASSAY PRINCIPLE

The Cow Milk Express used a monoclonal antibody specific for the bovine beta-lactoglobulin (1). This purified antibody has been used to sensitise blue latex particles.

A drop of the milk is mixed with a drop of the latex test reagent. If agglutination of the latex reagent occurs, then at least 1% of sterilized bovine milk is present in the sample tested.

3. KIT CONTENT

- a flask of latex particles: latex particules sensitised with the anti-bovine beta-lactoglobulin monoclonal antibody.
- disposable micro-pipettes: used for the collection of milk sample.
- disposable plastic spatula: used for mixing latex particules with the milk sample.
- test cards. Each card has 6 test circles. If only two circles are to be used, then the card may be cut and the unused portion saved for later use.
- an instruction manual.

	Reference L_CME_30	Reference L_CME_120
latex particles	1 flask of 1,3mL	1 flask of 4,9mL
micro-pipettes	30	120
plastic spatula	30	120
test card	5	20

4. STORAGE CONDITIONS / SHELF LIFE

Store the flask of latex at +2°C to +8°C. **Do not freeze.**

All reagents will remain stable until the expiry date shown on the bottle labels. Do not use beyond the expiry date.

5. STORAGE AND PREPARATION OF THE SAMPLES

Milk and whey samples can be stored at +2°C to +8°C for 1 to 3 days.

For long term storage, freeze the samples (- 20°C). After thawing, leave the samples equilibrate to room temperature before performing the assay. If the samples have a granular appearance after thawing, filter them (0.45 µm filter) or centrifuge in Eppendorf tube for few minutes.

For the whole milks, take care to remove the maximum of fat from the sample before running the assay.


6. INSTRUCTION

1. Bring the latex reagent to room temperature. **Mix the latex suspensions by vigorous shaking.** Expel any excess latex from the dropper flask to ensure complete mixing. Mix the samples. Identify the assay circles on the test card.
2. Use a disposable micro-pipette to collect milk sample. Load one drop of the milk sample onto one of the test circles reaction card. Discard the micro-pipette.
3. Dispense one free-falling drop of the test latex beside the milk sample. Take care not to touch the milk sample with the dropper of the latex flask.
4. Mix the two drops together using a plastic spatula and spread out to cover most of the test circle area. Use a new plastic spatula for each sample.
5. Incubate 2 minutes by slowly rocking (right, up, left and down) the test card.

7. READING AND INTERPRETATION

Positif Result : A result is positive (presence of bovine milk – circle nb#1 on the picture below) if agglutination of the latex occurs. Agglutination corresponds to a discoloration of the latex solution with the apparition of dark blue aggregates.

Negative result : A result is negative when agglutination do not occurs after 2 minutes of rocking. The mix remains of an homogenous blue color (circle nb#2 on the picture below).

Further indications :

Samples that contains more than 1% of raw bovine milk leads to an immediate agglutination.

The reference incubation time is 2 minutes. By extending the incubation time, it is possible to increase the sensitivity of the test. Do not exceed 4 minutes for the incubation.

8. TECHNICAL SPECIFICATIONS

The Cow Milk Express has been tested with a wide range of raw and heat treated milks (pasteurized, sterilized). The technical data presented below were obtained by dilution of a mixture of three commercial UHT cow's milk in a large mixing UHT goat milk.

The reference time of incubation to proceed to the interpretation of results is 2 minutes.

	Type of milk	Result
100%	sterilised Cow's milk	++
10 %	sterilised Cow's milk	++
1 %	sterilised Cow's milk	+
1%	raw Cow's milk	++
0,5%	raw Cow's milk	+
0,25%	raw Cow's milk	+
100%	raw Goat's milk	-
100%	sterilised Goat's milk	-
100%	raw Sheep's milk	-
100%	sterilised Sheep's milk	-

Legend :

- ++ => instant agglutination, after less than 2 minutes
- + => agglutination in 2 minutes
- => no agglutination

9. SAFETY RECOMMENDATION

This product is for laboratory use only. Do not freeze. Always use reagents with the same lot number. Reagents contains 0,09% sodium azide as preservative.

REFERENCES

- 1 - **Production and epitopic characterization of monoclonal antibodies against bovine beta-lactoglobulin.** Venien A et al. ; J Dairy Sci 1997 Sep ; 80 (9) : 1977-87.

IDBiotech

ImmunoDiffusion Biotechnologies SARL

Rue Marie Curie

63500 Issoire – FRANCE

Tel: +33 (0)4 73 54 95 01

Fax: +33 (0)4 73 54 43 94

Web : www.idbiotech.com

e-mail : contact@idbiotech.com