
50 11 T

 **CRISON**

50 11 T

1 Conector BNC + banana.
Connettore BNC + banana.
Connecteur BNC + banana.
Connector BNC + banana.

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Introduction

The **50 11 T** is a pH low maintenance electrode, as it does not require monitoring of electrolyte level, with built-in A.T.C. Pt 1000. The open diaphragm allows direct contact between the sample and the reference electrolyte. This avoids the classical problem of diaphragm clogging.

Electrode preparation

Remove the storage protector containing electrolyte by unscrewing it ③.

Verify that there are not air bubbles inside the membrane ⑦. They can be removed by shaking the electrode as a clinical thermometer.

Calibration

⚠ See *pH-meter instruction manual*.

Measurements. Recommendations

- During calibration and measurement the buffers and the samples must be slightly stirred. The stirring speed should be identical in both cases.
- The sample must cover the electrode's diaphragm ⑥.
- After measurement, rinse the electrode with distilled water. Dry it with soft tissue, only by contact, without rubbing.
- Between measurements the electrode's membrane should be immersed in aqueous solution. The most appropriate place is the storage protector with some electrolyte.

⚠ Do not leave the electrode immersed in distilled water.

Do not leave the electrode immersed in solutions with $\text{pH} < 2$.

Maintenance

- Maintain the membrane's surface clean ⑦. Protect it from friction or knocking.
- When an electrode gets contaminated, the simple rinsing with distilled water is not sufficient. The most suitable cleaning solution is one which acts in a more selective way over the contamination.

Proteins: Immerse the electrode in 2% NaOH solution during approx. 20 min.

Fats, oils: Wash the electrode with ethanol or soapy tap water at approx. 40°C.

Afterwards, rinse the electrode with plenty of distilled water and leave for rest in the reference electrolyte.

Troubleshooting

... the pH-meter always measures the same pH value.

Short circuit in the electrode or connector.

Substitute the electrode.

... unstable reading.

Air bubble in the membrane.

Shake the electrode as a thermometer.

... impossible to perform calibration.

... very slow response.

Dirt on the membrane.

De-hydrated membrane.

Electrode ageing or wear.

Clean, regenerate or substitute the electrode.

Storage

Place the electrode in the storage protector ③ with Crisolyt G (KCl with glycerine).

Life of an electrode

The average life expectancy of a pH electrode is approximately one year depending on the working conditions.

An electrode used for more than six months that does not respond appropriately after regeneration, must be replaced by a new one.

⚠ *If an electrode has too short life span, less than three months, ask us about it. Probably there is a more suitable electrode for the application.*

Specifications

pH measuring range	2 ... 14
Operating temp.	0 ... 80 °C
Reference element	sleeved Ag wire coated with AgCl
Diaphragm	open, circular
Electrolyte	solid
Body material	glass
Min. immersion depth	15 mm 
Temperature probe	Pt 1000

Applications

Food industry, viscous media, dirty samples, etc.

Limitations

Samples with pH < 2.

Temperatures exceeding 80 °C.

Samples with low conductivity, distilled water.

Warranty

CRISON guarantees this electrode against manufacturing defects.

Validity: 6 months from shipment date.

Limitations:

- Accidental damage.
- Inadequate applications.
- Non-fulfilment of the recommendations.
- Normal wear and tear of the electrode.

QUALITY CERTIFICATE

The attached electrode 50 11 T has passed the quality test and follows the next specifications:

- Asymmetry potential < ± 15 mV.
- Sensitivity, pH 4...7 (at 25 °C) > 98%.
- Response time, pH 4...7 < 20 s.
- Temperature error $\leq 0.5^\circ\text{C}$

CRISON INSTRUMENTS, S.A.

Riera Principal 34 - 36
08328 ALELLA (Barcelona)

Service

Tel.: +34 935 550 318
Fax: +34 935 400 857
E-mail: service@crison.es

Tel.: +39 059 651 922
Fax: +39 059 652 011
E-mail: crison@crison.it

www.crisoninstruments.com