PAL TEST

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**Specific immunocapture assay kit for the detection
of bovine alkaline phosphatase in cow milk and
cow milk cheese**

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Specific immunocapture assay kit for the detection of bovine alkaline phosphatase in cow milk and cow milk cheese

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This document must be fully read before running the kit for the first time.

Quality Management System **ISO9001:2008** certified by

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TABLE OF CONTENTS

1.	INTRODUCTION	20
2.	FIELD OF APPLICATION	20
3.	ASSAY PRINCIPLE	20
4.	KIT CONTENT	21
5.	STORAGE CONDITIONS	21
6.	REAGENTS AND EQUIPMENT REQUIRED	22
7.	REAGENT PREPARATION	22
8.	SAMPLE PREPARATION	23
9.	STANDARD PREPARATION	24
10.	INSTRUCTIONS	25
11.	RESULTS	27
12.	TECHNICAL ADVICES	28
13.	SAFETY RECOMMENDATIONS	28
14.	SPECIFICATIONS	28
	ANNEXE / ANNEX	34

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1. INTRODUCTION

The ALP Test is an immunocapture assay used for measuring the levels of alkaline phosphatase (ALP) in cow milk and cow milk cheese. It is based on a monoclonal antibody specific to the ALP of cow milk which allows the measurement of the ALP enzymatic activity as an indicator of the effects of heat treatment applied to milk & dairy products. Since the antibody is only specific to the ALP of cow milk, this assay does not detect exogenous alkaline phosphatases of bacterial, microbial or fungal origin. Furthermore the assay is not sensitive to the presence of ALP inhibitors.

This assay can run 42 samples in duplicate, simultaneously, and can also be automated.

This test cannot replace the current European standard methods used for the official testing of heat treated dairy products.

2. FIELD OF APPLICATION

The ALP TEST can be used for the following products:

- Whole or skimmed cow milk
- Cow milk cheeses of the following types:
 - Blue-veined cheeses
 - Uncooked hard cheeses
 - Hard cooked cheeses
 - Soft-ripened cheeses
 - Fresh cheeses

3. ASSAY PRINCIPLE

The ALP TEST is presented in the format of an ELISA plate coated with a monoclonal antibody specific to the alkaline phosphatase found in cow milk. The assay technique is based on a 2-step immunocapture:

- After extraction with butanol-1, the calibration standards and the diluted samples are distributed in the wells of the microplate.
- After rinsing, the ALP captured by the antibodies is detected by the addition of a chromogenic substrate specific of the ALP (pNPP).

After addition of the stop solution, the optical densities are measured at 405 nm with a microplate reader.

The results are expressed in equivalent-milliunit per litre (Eq.mU/l) for milk and in equivalent-milliunit per gram (Eq.mU/g) for cheese (*).

() The standard unit for the alkaline phosphatase assay is (the) milliunit per litre (mU/l) for milk and (the) milliunit per gram (mU/g) for cheese. These units are routinely applied when the ALP analyses are carried out following the two standard methods described in ISO 11816-1:2006 and ISO 11816-2:2003, respectively.*

The ALP TEST evaluates the amount of active enzyme captured by the monoclonal antibody. For these types of tests, the results are usually expressed in terms of concentration, based on a calibration range containing known quantities of the purified analyte. In the absence of reference material expressing known concentrations of the alkaline phosphatase and since it is not possible to obtain the purified enzyme, the ALP TEST has been calibrated in relation to the results obtained using the standardised method (ISO11816-1 :2006), and then the results are expressed as “equivalent” to the results obtained with this method. However, as it has been shown during a comparative study with both methods, based on statistical analysis of the accuracy profile, the ALP TEST cannot be substituted to the standard method for determining ALP in cheeses.

4. KIT CONTENT

The kit is delivered in a box containing:

- 6 16–wells strips (totally 96 wells) coated with a monoclonal antibody specific to cow milk ALP, wrapped in re-sealable aluminium pouches.
- 1 strip holder.
- 6 adhesive covers for 16-well strips.
- 1 bag containing 40 disposable transfer pipettes for pipeting sample extracts.
- Reagent **R1**: Dilution buffer for the calibration range and for the samples - 2 x 30 ml vials.
- Reagent **R2**: Standard solution at 600,000 Eq.mU/l – 1 x 0.5ml, microtube
- Reagent **R3**: Wash buffer, 20-fold concentrated – 1 x100 ml vial.
- Reagent **R4**: Substrate dilution buffer - 1 x 32 ml vial.
- pNPP substrate: 6 x 5mg tablets in individual hermetically sealed pouches – **to be stored at -20°C upon receipt.**
- Reagent **R5**: Stop solution, 1.5M NaOH - 1 x 20 ml vial.

5. STORAGE CONDITIONS

Store all reagents but the pNPP substrate at a temperature ranging from +2°C to +8°C.

Upon receipt, store the pouches containing the pNPP substrate at -20°C.

Do not use the reagents after their expiry date.

Prevent all the reagents, especially the substrate, from light exposure during storage and during assay incubation steps.

6. REAGENTS AND EQUIPMENT REQUIRED

- Distilled or deionised water
- 1-butanol (n-butanol): $\text{CH}_3(\text{CH}_2)_3\text{OH}$, FW 74.12 (e.g. product code 360465 Sigma-Aldrich)
- Precision micropipettes (50 – 1000 μl) with corresponding disposable tips
- 8-channel micropipette with variable volume (50-200 μl) and corresponding disposable tips
- 300 μl mono or multichannel dispenser or manual or automatic plate washer
- 12 ml disposable glass haemolysis tubes (for sample extraction)
- Disposable 1.2ml \varnothing 8.8mm x 45mm microtitration tubes for sample dilution
- Disposable plastic 50 ml conical bottom centrifuge tubes with stoppers
- Reservoir (for multichannel pipette)
- Electric mincer
- Polytron®-type homogeniser
- Orbital shaker (type IKA-Vibrax®) with tray for microplate
- Water bath (42 - 45°C)
- Microplate reader with 405 nm filter
- Data treatment software including a 4 parameter logistics regression model $Y = (A-D)/(1+(X/C)^B) + D$
- Centrifuge with 12 ml haemolysis tube holder.
- Incubator set to 35-38°C
- Scalpels
- Positive and negative control samples

7. REAGENT PREPARATION

- **Before use bring all reagents and the aluminium pouch containing the plates to room temperature**
- Homogenise each vial manually by turning upside down several times before use.
- Do not mix reagents from different lots.

Preparation of 1X wash buffer

For an assay using the 6 strips (96 wells), proceed as follows:

Into a clean new vial, mix the 100 ml of 20X wash buffer (R3) with 1900 ml of distilled or deionised water.

1. Homogenise without generating foam
2. The solution can be stored at +2°C/+8°C during 15 days

If the kit is used strip per strip, prepare the necessary volume of wash buffer (50 ml 1X wash buffer per strip is required, plus the expected dead volume required for the washing system).

8. SAMPLE PREPARATION

Milk samples can be kept frozen at -20°C for 1 month before analysis or at $+2^{\circ}\text{C}/+8^{\circ}\text{C}$ for a maximum of 48h.

Cheese samples can be stored at $+2^{\circ}\text{C}/+8^{\circ}\text{C}$. Do not freeze.

Bring milk and cheese samples to room temperature before extraction.

8.1. Milk samples

Shake the sample vigorously to ensure homogeneity, then:

1. Pipet 3ml milk using an automatic pipette and pour into a disposable glass tube. Add 3ml butanol-1. Put a stopper on the tube then shake vigorously by hand for 5 to 10 seconds. Vortex for 30 seconds. Centrifuge for 30 minutes between 2500g and 3500g.
2. Collect the aqueous phase under the fat layer with a transfer pipette (expelling the air slowly while crossing the alcohol phase and the fat layer so as not to clog the end of the transfer pipette).
3. The recommended dilution for a milk extract is 1/200. This dilution factor must be adapted in relation to the thermal treatment applied to the sample. The minimum possible dilution of the extract is 1/5.
For a dilution of 1/200, proceed in two steps as follows:
Dilution 1/10: 50 μl extract + 450 μl dilution buffer R1. Vortex
Dilution 1/200: 25 μl dilution 1/10 + 475 μl dilution buffer R1. Vortex.

8.2. Cheese samples

1. Remove the crust or the surface of the cheese sample with a clean scalpel. Take a representative amount of the cheese and homogenise it in a grinder. Wash the bowl carefully with hot water containing detergent and dry it between each sample.
2. Weigh about 5g precisely of the homogenate and place it in a 50ml centrifuge tube with a conical bottom. Add twice the same weight of deionised water. Place the tube in a water bath at $42-45^{\circ}\text{C}$ for 15 minutes.
3. Homogenise the sample in a Polytron[®] at 12000 rpm for 30 seconds. Rinse the blade thoroughly with hot water containing some detergent and dry it between each sample in order to avoid any cross-contamination. Return the tubes back to the water bath at $42-45^{\circ}\text{C}$ for 15 min. Mix thoroughly by inversion.

NB: in the case of fresh curd cheeses, check the pH of the extract. This must be higher than 6.0 to ensure the validity of the results.

4. Transfer 3ml of the mixture into a disposable glass tube. Add 3ml butanol-1. Put a stopper on the tube then shake vigorously by hand. Vortex for 30 seconds. Centrifuge for 30 minutes at 2900 G.
5. Collect the aqueous phase under the emulsion layer with a transfer pipette (expelling the air slowly while crossing the butanol phase and the emulsion layer so as not to clog the end of the transfer pipette).
6. The recommended dilution for an extract of cheese made from crude milk is 1/200. This dilution factor must be adapted in relation to the thermal treatment applied to the sample. The minimum possible dilution of the extract is 1/5.
For a dilution of 1/200, proceed in two steps as follows:
Dilution 1/10: 50 µl extract + 450 µl dilution buffer R1. Vortex
Dilution 1/200: 25µl dilution 1/10 + 475 µl buffer R1. Vortex.

NB: *The aqueous phases after extraction with butanol can be kept undiluted at +2°C/+8°C for a maximum of 10 days before analysis. This can notably allow a sufficient number of extracts to be prepared and analysed simultaneously. The dilutions must be carried out on the day of analysis.*

9. STANDARD PREPARATION

A working solution at 15,000 Eq.mU/l is prepared extemporaneously by diluting the stock solution R2 to 1/40 in the dilution buffer R1 as follows:

Working sol. 15,000 Eq.mU/l 50µl R2 + 1950µl R1

The standard solutions are then prepared extemporaneously from the working solution at 15,000 Eq.mU/l as follows:

- STD1	5000 Eq.mU/l:	125µl (15,000 Eq.mU/l) + 250 µl R1
- STD2	3000 Eq.mU/l:	75 µl (15,000 Eq.mU/l) + 300 µl R1
- STD3	1000 Eq.mU/l:	25 µl (15,000 Eq.mU/l) + 350 µl R1
- STD4	500 Eq.mU/l:	50 µl (STD1) + 450 µl R1
- STD5	300 Eq.mU/l:	50 µl (STD2) + 450 µl R1
- STD6	100 Eq.mU/l:	50 µl (STD3) + 450 µl R1

Sample and standard dilutions can be done directly in microtitration tubes. The volumes indicated here allow for duplicates of each sample or standard solution.

10. INSTRUCTIONS

Bring all reagents to room temperature before use. It is recommended that all samples, controls and standards are analysed in duplicate.

1. Prepare the reagents, samples, controls and standards as described previously. Remove the required number of strips from the aluminium pouch and place them on the strip holder. Carefully reseal the aluminium pouch and return it at +2°C/+8°C. Be careful not to touch the underneath of the plates.
2. Distribute 300 µl of wash buffer (**R3 1X**) into each well. Then eliminate the washing liquid by turning the plate upside down or by aspiration using a pipet. Repeat the operation 4 times. After the final wash, eliminate any residual droplets by tapping the overturned plate onto clean and dry absorbent paper. *NB: Do not let the wells drying out between consecutive steps.*
3. Distribute, using a multichannel pipette, **100 µl/well** of standards, controls and samples in the corresponding wells. Cover the plate with an adhesive cover and shake gently for a few seconds on a plate shaker. Incubate for **1 hour** at room temperature (18-25°C) or **30 minutes under gentle shaking**.
4. During the incubation, prepare extemporaneously the volume of substrate necessary according to the table below:

Number of strips used	Number of pNPP tablets	Volume of R4 buffer (ml)
1	1	5
2 to 5	2	10
6	3	15

Preheat the substrate by placing it in the incubator at 37°C for about 20 minutes.

5. Empty the liquid from the wells by turning the plate upside down or by aspiration. Repeat the washing operation as described in step 2.
6. Using a multichannel pipette, distribute 100 µl of substrate into each well. Cover the plate with the adhesive cover and shake gently for a few seconds on the plate shaker. Incubate for 2 hours at 35-38°C or 1 hour at 35-38°C if step 3 was carried out with agitation. A yellow colour develops.
7. Using a multichannel pipette, distribute 50 µl of stop solution (R5) into each well. Cover with the adhesive cover and shake gently for a few seconds on the plate shaker.
8. Remove adhesive cover and read the absorbance at 405nm using a microplate reader.

NB : in order to optimize the incubation time of substrate (to obtain OD max corresponding to OD max of the microplate reader), it's possible to read the absorbance at 405nm during step 6, before the addition of stop solution.

Instructions summary

Sample preparation

MILK

↓

CHEESE

Take a representative sample and

grind

↓

Weigh 5 g + twice weight of deionised water

↓

Water bath 45°C - \approx 15min

↓

Homogenise in Polytron®

↓

Water bath 45°C- \approx 15min

↓

3ml in glass tube + 3ml butanol-1

↓

Shake vigorously by hand + **vortex** - \approx 30 seconds

↓

Centrifuge at 2900g - \approx 30min

↓

Take aqueous phase

↓

Dilute aqueous phase

↓

Prepare working solution 15,000 Eq.mU/l based on stock solution R2

Prepare standards based on working solution 15,000 Eq.mU/l

Sample analysis

Wash 4 times (4 x 300 μ l)

↓

Distribute 100 μ l/well of each calibration and sample in duplicate

↓

Incubate 1 hour at room temperature or 30 min under agitation

↓

Prepare and preheat substrate 20 min at 37°C

↓

Wash 4 times (4 x 300 μ l)

↓

Distribute 100 μ l/well of substrate

↓

Incubate 2 hours at +37°C or 30 min to 1 hour if previous incubation carried out with agitation

↓

Distribute 50 μ l/well of stop solution

↓

Read absorbance at 405nm

↓

Interpret results


11. RESULTS

The calibration curve is obtained by plotting on the Y-axis the optical density (OD) read for each standard and on the X-axis their values in Eq.mU/l (logarithmic scale).

The calibration curve can be generated by using ideally a logistic 4 parameter model ($Y = (A-D)/(1+(X/C)^B) + D$) using appropriate software to analyse the data.

In the absence of this type of software or this regression model, it is also possible to use the graphic functions of the Excel software or similar software to calculate the results using polynomial regression of the order of 4. *A typical model is available on request from IDBiotech.*

An example of the calibration curve is given below:

Values in Eq.mU/l or Eq.mU/g of samples are then calculated as follows:

- **Milk samples:**

Value in Eq.mU/l = value in Eq.mU/l sample x dilution factor

- **Cheese samples:**

Value in Eq.mU/g =
$$\frac{\text{Value in Eq.mU/l sample} \times (\text{weight of water in g})}{(\text{weight of cheese in g} \times 1000)}$$

12. TECHNICAL SPECIFICATIONS

Any change at the level of the operator, pipette technique or washing, time of incubation and/or temperature variations may lead to variations in the final results.

To obtain precise and reproducible results, the calibration and sample distribution must be done carefully. Only use precision pipettes that have been metrologically checked and for which the verification is still valid.

Avoid creating bubbles in the solutions during dilutions.

To avoid cross-contaminations, change the pipette tips between each sample, calibration and reagent.

Use individual flasks and reservoirs labelled for each reagent.

Respect the number of washing cycles as the washing stages are very important.

It is imperative to have equivalent timing between the steps of reagent distribution and washing in each well.

The stop solution must be distributed in the wells in the same order and with the same timing as the substrate.

13. SAFETY RECOMMENDATIONS

Be careful to respect all hygiene and safety recommendations concerning the handling of reagents and notably the organic solvent used for extraction (work in fume cupboard). The work premises must be well ventilated. Avoid the presence of ignition or heat sources (flames, sparks, solar radiation, etc.). Avoid inhaling aerosol vapours and wear appropriate protective clothing to avoid all contact with skin and eyes.

14. TECHNICAL SPECIFICATIONS


The technical specifications of the ALP Test were determined following the main directions of the document "Validation of analytical procedures: text and methodology Q2(R1) November 2005" published by the ICH (International Conference on Harmonization).

14.1. Linearity

Two regression models were studied.

- **4 parameter logistic model $Y = (A-D)/(1+(X/C)^B) + D$**
-


	Absorbance 405 nm (OD)	OD/ODmax (%)	Mean of calculated concentration (Eq. mU/l)	CV (%)
STD1	3.029 2.818 3.030 2.891 2.982 3.035 2.900 2.971 2.863 2.897	-	5003.9	2.7
STD2	1.893 1.741 1.854 1.747 1.788 1.835 1.737 1.815 1.746 1.717	61	2990.1	3.5
STD3	0.705 0.645 0.680 0.650 0.661 0.673 0.654 0.667 0.649 0.611	22	1014.6	4.3
STD4	0.402 0.369 0.394 0.369 0.377 0.386 0.369 0.384 0.368 0.340	13	506.2	6.2
STD5	0.285 0.260 0.275 0.260 0.267 0.275 0.262 0.270 0.262 0.235	9	303.6	8.1
STD6	0.154 0.147 0.152 0.148 0.149 0.150 0.147 0.158 0.147 0.147	5	86.6	8.3

$Y = (A-D)/(1+(X/C)^B) + D$				
A	B	C	D	R ²
0.107	1.04	1.37E+05	91.5	0.998

- Polynomial regression of the order of 4: $Y = A + BX + CX^2 + DX^3 + EX^4$

	Absorbance 405 nm (OD)	OD/ODmax %	Mean of calculated concentration (Eq. mU/l)	CV (%)
STD1	3.029 2.818 3.030 2.891 2.982 3.035 2.900 2.971 2.863 2.897	-	5000.3	2.7
STD2	1.893 1.741 1.854 1.747 1.788 1.835 1.737 1.815 1.746 1.717	61	2999.9	3.4
STD3	0.705 0.645 0.680 0.650 0.661 0.673 0.654 0.667 0.649 0.611	22	1000.5	4.4
STD4	0.402 0.369 0.394 0.369 0.377 0.386 0.369 0.384 0.368 0.340	13	497.5	6.1
STD5	0.285 0.260 0.275 0.260 0.267 0.275 0.262 0.270 0.262 0.235	9	303.2	7.7
STD6	0.154 0.147 0.152 0.148 0.149 0.150 0.147 0.158 0.147 0.147	5	102.9	6.3

$Y = A + BX + CX^2 + DX^3 + EX^4$					
A	B	C	D	E	R ²
0.0903	0.000581	-1.76E-08	5.70E-12	-5.24E-16	0.998

4 parameter logistic model $Y = (A-D)/(1+(X/C)^B) + D$ as well as the order 4 polynomial regression $Y = A + BX + CX^2 + DX^3 + EX^4$ describe the optimal “concentration / OD” relation with a correlation coefficient R^2 superior or equal to 0.99.

The 4 parameter logistic model was used to determine the technical specifications presented in the rest of this document.

14 .2. Detection limit – Quantitation limit

The **LOD** (Limit of Detection) was determined by interpolating the mean based on 36 Optical Densities plus 3 standard deviations. It is **50 Eq.mU/l**.

The **LOQ** (Limit of Quantification) was determined by interpolating the mean based on 36 Optical Densities plus 10 standard deviations.

It is **250 Eq.mU/l**.

- **Application to milk testing**

The **minimum** dilution factor (**1/5**) of samples must be taken into account. Thus the LOD for **milk samples** is **250 Eq.mU/l** and the LOQ is **1250 Eq.mU/l** (i.e. **about 0.1% of raw milk in pasteurised milk**).

- **Application to cheese testing**

The extraction of the sample (weight of cheese and volume of buffer) must be taken into account, as well as the minimum dilution factor (**1/5**) of samples. Thus the LOD for **cheese samples** is **2.5 Eq.mU/g** and the LOQ is **12.5 Eq. mU/g**.

14.3. Repeatability – Intermediate reliability

Whole milk - Study data:

ALP (Eq.mU/l)	Intra-assay repeatability (CV%)	Inter-assay intermediate reliability (CV%)
2454	2.0	13.0
837	3.8	12.3
456	5.8	15.6

Soft cheese - Study data:

ALP (Eq.mU/g)	Intra-assay repeatability (CV%)	Inter-assay intermediate reliability (CV%)
1866	5.3	18.0
1080	5.6	20.0
633	5.8	21.0

Fresh curd cheese - Study data:

ALP (Eq.mU/g)	Intra-assay repeatability (CV%)	Inter-assay intermediate reliability (CV%)
1326	5.5	11.0
278	3.2	23.0

Blue-veined cheese - Study data:

ALP (Eq.mU/g)	Intra-assay repeatability (CV%)	Inter-assay intermediate reliability (CV%)
1655	2.8	7.5
159	2.2	11.9

Uncooked pressed cheese - Study data:

ALP (Eq.mU/g)	Intra-assay repeatability (CV%)	Inter-assay intermediate reliability (CV%)
5508	3.1	9.9
44	10.1	12.4

14.4. Precision (Whole milk)

Theoretical concentration (Eq.mU/l)	Measured concentration (Eq.mU/l)	Recovery (%)
-	659 805	-
494 854	464 916	94
329 903	310 698	94
164 951	150 442	91
65 981	57 054	86

14.5. Specificity

The monoclonal antibody used in this ALP Test has not shown any cross reactivity with the following alkaline phosphatases:

- Bovine intestinal alkaline phosphatase (SIGMA P6774 – FLUKA 79390)
- Fungal alkaline phosphatase (penicillium roqueforti)
- *E.coli* alkaline phosphatase (SIGMA P4377)
- Alkaline phosphatase from goat and sheep milk

Comparison with fluorometric method

This study has been organised and funded by the “Centre National Interprofessionnel de l'Economie Laitière” (CNIEL).

The objective consisted of comparing the results obtained for milk and cheese by the ALP TEST with those obtained using standard fluorometric methods ISO 11816-1 and ISO 11816-2 (in a version under revision). The statistical method of interpretation used was the procedure with precision profile in accordance with standard NF V03-110. The results of this study are available from the CNIEL.

ANNEXE / ANNEX

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	1	2	3	4	5	6	7	8	9	10	11	12
A												
B												
C												
D												
E												
F												
G												
H												

IDBiotech

*ImmunoDiffusion Biotechnologies SARL
Rue Marie Curie
63500 Issoire – FRANCE*

Tel: +33 (0)4 73 54 95 01

Fax: +33 (0)4 73 54 43 94

Web : www.idbiotech.com

e-mail : contact@idbiotech.com