50 71

 **CRISON**

50 71

1

Conector: telefónico (8c).
Connettore: telefonico (8c).
Connecteur: téléphonique (8c).
Connector: telephonic (8c).

2

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Introduction

The **50 71** is a platinum cell for very low conductivities measurements, designed for use with the **CRISON BASIC** and **GLP “+” series** conductivity meters.

It has 2 platinum electrodes and an integrated Pt 1000 temperature sensor.

Cell preparation

Before use for the first time immerse the electrodes in ethanol for 15 s, and rinse with distilled water. Apply the same procedure for a cell not being in use for long time.

Calibration

Before measurement, calibrate with standard solutions.

 See *User's Manual of conductivity meter*.

Measurements. Recommendations

- Both, calibration and measurement must be performed under slight stirring of standards and samples. The stirring speed in the two cases should be identical.
- The sample must cover the measuring electrodes ⑥.
- Ensure that during measurement there are not air bubbles over the measuring electrodes.
- After measurement, rinse the cell with distilled water or with sample to be measured.

Warning.

The cell is made of glass, any abrupt manipulation could break it.

Maintenance

Keep the electrodes surface clean ⑥.

If the thorough rinse with distilled water is not sufficient, suitable solvent for the elimination of sample rests can be used. This solvent must not produce chemical attack over the cell materials.

Inorganic contamination: Immerse the cell in diluted hydrochloric acid or in chromic acid. Afterwards rinse thoroughly with distilled water.

Organic contamination: Use suitable solvent (see cell's materials), next utilize alcohol and distilled water as in Cell preparation is described.

 Do not apply any mechanical cleaning over the platinum cell.

Storage

Store the cell dry.

Life of the cell

The cell lifespan is unlimited if suitable maintenance is applied and the cell does not break.

Troubleshooting

... Reading 0 independently of the measured sample.

Broken cable.

Replace the cell.

... Unstable reading.

... Impossible calibration.

Dirty measuring electrodes or air bubbles on their surface.

Clean or shake slightly the cell.

... The obtained readings differ a lot from the real values.

Verify that the conductivity is covered by the measuring interval of the cell.

Re-calibrate using standards in good conditions.

.... slow cell response

Regenerate the cell.

Specifications

Approx. constant	0,1 cm ⁻¹
Measuring range	0,05 µS/cm... 30 mS/cm
Operating temp.	-30... 85 °C
Body materials	glass
Electrode material	platinum
Min. immersion depth	45 mm

Applications

Low conductivity samples, highly purified and de-mineralised waters, etc.

Ultrapure water, using as accessory the glass chamber code 91 18.

Limits

Conductivities of over 30 mS/cm.

Very dirty samples.

Warning.

The measuring range is only guaranteed when working with CRISON conductivity meters

Warranty

CRISON guarantees this cell against manufacturing defects.

Validity: 6 months from shipment date.

Limitations:

- Accidental damage.
- Inadequate applications.
- Non-fulfilment of the recommendations.
- Normal wear and tear of the cell.

QUALITY CERTIFICATE

The attached cell **50 71** has passed the quality test and fulfils the following specifications:

- % error ≤ 5 % between 0,05 µS/cm...30 mS/cm.
- Temperature error ≤ 0.5 °C.

CRISON INSTRUMENTS, S.A.
Riera Principal 34 - 36
08328 ALELLA (Barcelona)

Service

Tel.: +34 935 550 318
Fax: +34 935 400 857
E-mail: service@crison.es

Tel.: +39 059 651 922
Fax: +39 059 652 011
E-mail: crison@crison.it

Tel.: +33 0555 639 726
Fax: +33 0555 639 727
E-mail: info@crison.fr

www.crisoninstruments.com