Citrate  $\xrightarrow{\text{CL}}$  Oxaloacetate + Acetate

Oxaloacetate + NADH + H<sup>+</sup>  $\xrightarrow{\text{MDH}}$  L-Malate + NAD<sup>+</sup>

Pyruvate + NADH + H<sup>+</sup>  $\xrightarrow{\text{LDH}}$  L-Lactate + NAD<sup>+</sup>

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#### Sample preparation:

The concentration of citrate in the sample used in the assay procedure has to be between 0.04 to 0.4 g/l.

#### Assay procedure:

Wavelength: 340nm / Optical path: 1cm / Temperature: 20-37°  
Measurement: against water or air.

	Blank	Sample
R 1	1 ml	1 ml
R 2	0.2 ml	0.2 ml
R 3	0.02 ml	0.02 ml
Water	1 ml	0.9 ml
Sample	0	0.1ml
Mix and read DO	<b>DO 1</b>	<b>DO 1</b>
R 4	0.02 ml	0.02 ml
Mix and read DO after 10 mn	<b>DO 2</b>	<b>DO 2</b>

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Dilution factor of the sample has to be considered in the calculation.

#### Storage instructions and reagent stability

The reagents are stable up to the indicated month of expiry, if stored at 2-8 °C, contamination is avoided.

The dissolved R4 is stable about 5 days, if stored at 2-8 °C.

#### Warnings and precautions

Do not swallow. Avoid contact with the skin and mucous membranes.  
Take necessary precautions for the use of laboratory reagents.

#### General information on sample preparation

- Use liquid samples, clear, colorless and with pH practically 7.8, directly, or after dilution.
- Crush or homogenize solid or semi-solid samples
- Deproteinize samples containing proteins with perchloric acid; alternatively clarify with Carrez reagents.
- Extract samples containing fat with hot water

#### Assay control

citrate reagents must be validated with the use of the standard included in the kit. The standard is ready-to-use.

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