

50 21 T

50 21 T

1

Conexiones: BNC para pH + Banana para sonda Pt 1000.
Connessioni: BNC per pH + Banana per sonda Pt 1000.
Connexions: BNC pour pH + Banana pour capteur Pt 1000.
Connections: BNC for pH + Banana for Pt 1000 probe.

2

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pH electrode 50 21 T

Introduction

The **50 21 T** is a pH electrode with built-in A.T.C. Pt 1000, suitable for difficult samples.

Electrode preparation

Remove the storage protector containing electrolyte ⑤.

Remove the air bubbles that are inside the membrane ⑧. They can be removed by shaking the electrode as a clinical thermometer.

Remove the cap from the electrolyte filling port ③.

Before measuring, slide the plastic ring and return it to its original position, by pressing softly, to ensure the electrical contact, see fig A. If it is necessary, refill with reference electrolyte.

Calibration

⚠ See pH-meter instruction manual.

Measurements. Recommendations

- During calibration and measurement the buffers and the samples must be slightly stirred. The stirring speed should be identical in both cases.
- The sample must cover the electrode's diaphragm ⑦.
- After measurement, immediately rinse the electrode with distilled water. Dry it with soft tissue, only by contact, without rubbing.
- Between measurements the electrode's membrane should be immersed in aqueous solution. The most appropriate place is the storage protector with some electrolyte.

⚠ Do not leave the electrode immersed in distilled water.

Maintenance

- Periodically refill the electrode with electrolyte. Maintain the level as high as possible ④.
- Maintain the membrane's surface clean ⑧. Protect it from friction or knocking.
- When an electrode gets contaminated, the simple rinsing with distilled water is not sufficient. The most suitable cleaning solution is one which acts in a more selective way over the contamination (look for CRISON regeneration solutions).
- **Re-hydration.** A "dry" electrode recovers by immersing it in CRISON cleaning electrode solution, code 96 00, or in diluted HCl solution for several hours. It is recommended to read the information supplied with the electrode's regenerating solutions.

Troubleshooting

... always measures the same pH value.

Short circuit in the electrode or cable.

Substitute the electrode.

... unstable reading.

Air bubble in the membrane.

Shake the electrode as a thermometer.

... impossible to perform calibration.

... very slow response.

Dirt on the diaphragm or the membrane.

De-hydrated membrane.

Electrode ageing or wear.

Clean, regenerate or substitute the electrode.

Storage

Place the electrode in the storage protector ⑤ with Crisolyt (KCl 3M).

Close the refilling port with the stopper ③.

Life of an electrode

The average life expectancy of a pH electrode is approximately one year depending on the working conditions.

An electrode used for more than six months that does not respond appropriately after regeneration, must be replaced by a new one.

⚠ *If an electrode has too short lifespan, less than three months, ask us about it. Probably there is a more suitable electrode for the application.*

Specifications

pH measuring range	0 ... 14
Operating temp.	0 ... 60 °C
Reference element	cartridge with Ag ⁺ ion barrier
Diaphragm	sleeve
Electrolyte	CRISOLYT
Body material	glass
Min. immersion depth	15 mm 

Applications

Samples with very low conductivity.

Distilled water.

Samples with colloids, wines, paints, etc.

Very viscous sample, emulsions, creams...

Limitations

Temperatures > 60 °C.

Warranty

CRISON guarantees this electrode against manufacturing defects.

Validity: 6 months from shipment date.

Limitations:

- Accidental damage.
- Inadequate applications.
- Non-fulfilment of the recommendations.
- Normal wear and tear of the electrode.

QUALITY CERTIFICATE

The attached electrode 50 21 T has passed the quality test and follows the next specifications:

- Asymmetry potential < ± 15 mV.
- Sensitivity, pH 4...7 (at 25 °C) > 98%.
- Response time, pH 4...7 < 20 s.
- Temperature error ≤ 0.5 °C.

CRISON INSTRUMENTS, S.A.

Riera Principal 34 - 36
08328 Alella BARCELONA
SPAIN

Service

Tel.: +34 935 550 318
Fax: +34 935 400 857
E-mail: service@crison.es

Tel.: +39 059 651 922
Fax: +39 059 652 011
E-mail: crison@crison.it

Tel.: +33 0555 639 726
Fax: +33 0555 639 727
E-mail: info@crison.fr

www.crisoninstruments.com