ACHTUNG / ATTENTION / ATTENZIONE / LET OP / NB

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(EN) Important information for the evaluation!

Without hydrolysis, only the (dissolved) **orthophosphate** is measured.
The result of the orthophosphate measurement can be expressed as: mg/L PO₄-P (e.g. for process analysis), mg/L PO₄ (e.g. for analyses of drinking water or boiler water), mg/L P₂O₅ (e.g. for soils analyses)

With hydrolysis, all of the **phosphorus** (Total-P, P_{total}) is measured.
The result of the total phosphorus measurement can be expressed as: **mg/L P_{tot} = Display mg/L PO₄-P** (e.g. for monitoring threshold values in waste water), mg/L PO₄ (e.g. for analyses of drinking water or boiler water), mg/L P₂O₅ (e.g. for soils analyses).

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LCK 350

Phosphorus total / Phosphate ortho

! Please check the "Edition Date" (see data table). Safety advice and expiry date on package.

Principle

Phosphate ions react with molybdate and antimony ions in an acidic solution to form an antimonyl phosphomolybdate complex, which is reduced by ascorbic acid to phosphomolybdenum blue.

Range of Application

Waste water, drinking water, boiler water, surface water, process analysis

Interferences

The ions listed in **T1** have been individually checked up to the given concentrations and do not cause interference. We have not determined cumulative effects and the influence of other ions.

The measurement results must be subjected to plausibility checks (dilute and/or spike the sample).

Removal of Interferences

If phosphonic acids are present the time for hydrolysis in the thermostat must be increased to **2 h** at **100°C** in order to prevent low-bias results (see procedure for the determination of total phosphorus).

pH sample 2–10
Temperature sample/reagents 15–25 °C
In case of not working at the right recommended temperature an incorrect result may be obtained.

Note

Inverting the cuvette after hydrolysis improves the reliability of the result.

T1
5000 mg/L: SO ₄ ²⁻
2000 mg/L: Cl ⁻
1000 mg/L: K ⁺ , Na ⁺ , Ca ²⁺
500 mg/L: Mg ²⁺ , NO ₃ ⁻
50 mg/L: Co ²⁺ , Fe ²⁺ , Fe ³⁺ , Zn ²⁺ , Cu ²⁺ , Ni ²⁺ , I ⁻ , NO ₂ ⁻ , Cd ²⁺ , Sn ⁴⁺ , NH ₄ ⁺ , Mn ²⁺ , Al ³⁺ , Hg ²⁺ , Pb ²⁺ , SiO ₂
25 mg/L: Ag ⁺
10 mg/L: Cr ³⁺
5 mg/L: Cr ⁶⁺

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EN

1. **Carefully** remove the foil from the screwed-on **DosiCap Zip**.
2. Unscrew the **DosiCap Zip**.
3. Pipette **0.4 mL** sample.
4. Screw the **DosiCap Zip** back **tightly**; fluting at the top.
5. Shake firmly.
6. Heat in the thermostat.
HT 200 S: in standard program HT for **15 min**
Thermostat: 60 min at 100°C
7. Allow to cool to room temperature. Shake firmly.
8. Pipette into the cooled cuvette: **0.5 mL** Reagent **B** (LCK 350 B). Close Reagent B **immediately** after use.
9. Screw a **grey DosiCap C** (LCK 350 C) onto the cuvette.
10. Invert a few times. After **10 min** invert a few times more, thoroughly clean the outside of the cuvette and evaluate.

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