McFarland Standard

IVD

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IVD

SUMMARY AND EXPLANATION

McFarland Standard is a series of standards of different opacity, allowing the estimation of the density of bacterial suspensions. Its use is necessary for the standardization of microbiological methods. These standards are designated by the tube numbers of the scale initially described by McFarland.

PRINCIPLE

The density of the bacterial suspension is checked by comparison with a suspension of known opacity contained in an ampule of the same diameter.

CONTENT OF THE KIT

- 6 McFarland standards (0.5, 1, 2, 3, 4, 5), 17.75 mm in diameter (standard bioMérieux ampule)
- 1 package insert

COMPOSITION OF THE STANDARDS

	BaSO ₄		BaSO ₄
Standard 0.5	2.40 10 ⁻⁵ mol/l	Standard 3	1.44 10 ⁻⁴ mol/l
Standard 1	4.80 10 ⁻⁵ mol/l	Standard 4	1.92 10 ⁻⁴ mol/l
Standard 2	9.60 10 ⁻⁵ mol/l	Standard 5	2.40 10 ⁻⁴ mol/l

WARNINGS AND PRECAUTIONS

- For *in vitro* diagnostic use and microbiological control.
- For professional use only.
- Do not use standards past the expiration date.
- Before use, check that the ampules are intact.
- All specimens, microbial cultures and inoculated products should be considered infectious and handled appropriately. Aseptic technique and usual precautions for handling the bacterial group studied should be observed throughout this procedure. Refer to "NCCLS M29-A, *Protection of Laboratory Workers from Instrument Biohazards and Infectious Disease Transmitted by Blood, Body Fluids, and Tissue; Approved Guideline - December 1997*". For additional handling precautions, refer to "Biosafety in Microbiological and Biomedical Laboratories, HHS Publication No. (CDC) 93-8395, 3rd Edition (May 1993)", or to the regulations currently in use in each country.

STORAGE CONDITIONS

- The ampules should be stored in the dark at 2-30°C until the expiration date indicated on the packaging.
- Do not open the McFarland Standard ampules: the standard suspensions are used in the original, sealed ampules.

INSTRUCTIONS FOR USE

- Take the standard chosen as the reference.
- Prepare a bacterial suspension in an ampule of same diameter; mix well.
- Strongly shake the standard ampule.
- Immediately compare both opacities by placing the ampules against a black background.
- If necessary, adjust the bacterial suspension and shake the standard a second time before performing a new comparison.

INTERPRETATION

Equivalences

Standard	Bacterial Concentration (1) x 10 ⁶ /ml	Theoretical Optical Density (2) at 550 nm
0.5	150	0.125
1	300	0.25
2	600	0.50
3	900	0.75
4	1200	1.00
5	1500	1.25

1. Bacterial concentration depends on microorganism size. The numbers represent an average value valid for bacteria.
For yeasts, which are larger in size, these numbers should be divided by about 30.
2. Values correspond to optical densities of bacterial suspensions. The BaSO₄ solutions do not have the same optical density, because size and particle forms differ from those of bacteria and light is diffracted differently.

LIMITATIONS OF THE METHOD

- Do not use these ampules to test the DENSIMAT.

WASTE DISPOSAL

It is the responsibility of each laboratory to handle waste and effluents produced according to their type and degree of hazardousness and to treat and dispose of them (or have them treated and disposed of) in accordance with any applicable regulations.

WARRANTY

bioMérieux disclaims all warranties, express or implied, including any implied warranties of MERCHANTABILITY AND FITNESS FOR A PARTICULAR USE. bioMérieux shall not be liable for any incidental or consequential damages. IN NO EVENT SHALL BIOMERIEUX'S LIABILITY TO CUSTOMER UNDER ANY CLAIM EXCEED A REFUND OF THE AMOUNT PAID TO BIOMERIEUX FOR THE PRODUCT OR SERVICE WHICH IS THE SUBJECT OF THE CLAIM.

LITERATURE REFERENCES p. I
INDEX OF SYMBOLS p. II

bioMérieux® sa
au capital de 11 879 045 €
673 620 399 RCS LYON
69280 Marcy-l'Etoile / France
Tel. 33 (0)4 78 87 20 00
Fax 33 (0)4 78 87 20 90
<http://www.biomerieux.com>

bioMérieux, Inc
Box 15969,
Durham, NC 27704-0969 / USA
Tel. (1) 919 620 20 00
Fax (1) 919 620 22 11

Printed in France

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The logo is a registered and protected trademark of bioMérieux sa or one of its subsidiaries