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NCFB : National Collection for Food Bacteria (=NCFB), Institute of Food Research, Reading Laboratory, Earley Gate, Reading RG6 6BZ, ENGLAND

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Carbohydrates

SUMMARY AND EXPLANATION

API 50 CH is a standardized system, associating 50 biochemical tests for the study of the carbohydrate metabolism of microorganisms. API 50 CH is used in conjunction with API 50 CHL Medium for the identification of *Lactobacillus* and related genera and with API 50 CHB/E Medium for the identification of *Bacillus* and related genera, *Enterobacteriaceae* and *Vibrionaceae*. The complete list of those organisms that it is possible to identify with this system can be found in the Identification Table at the end of the package inserts of the associated media.

PRINCIPLE

The API 50 CH strip consists of 50 microtubes used to study fermentation of substrates belonging to the carbohydrate family and its derivatives (heterosides, polyalcohols, uronic acids).

The fermentation tests are inoculated with API 50 CHL Medium or API 50 CHB/E Medium, which rehydrates the substrates.

During incubation, **fermentation** is revealed by a **color change** in the **tube**, caused by the anaerobic production of acid and detected by the pH indicator present in the chosen medium. The first tube, which does not contain any active ingredient, is used as a negative control.

NOTE : The API 50 CH strip may be used to test two other pathways :

- **oxidation** which is revealed by a **color change** in the **cupule**, caused by the aerobic production of acid and detected by the pH indicator present in the chosen medium.
- **assimilation** which is revealed by **growth** of the organism in the **cupule** when the substrate is used as the only available source of carbon.

In this case, the choice of medium to be used for inoculation of the strips will depend on the metabolism and nutritional requirements of the microbial group to be tested (see literature references).

CONTENT OF THE KIT (Kit for 10 tests)

- 10 API 50 CH strips
- 10 incubation boxes
- 10 result sheets
- 1 package insert

COMPOSITION OF THE STRIP

The composition of the API 50 CH strip is given below in the list of tests :

Strip 0 - 9

Tube	Test	Active ingredients	QTY (mg/cup.)
0		CONTROL	-
1	GLY	GLYcerol	1.64
2	ERY	ERYthritol	1.44
3	DARA	D-ARAbinose	1.4
4	LARA	L-ARAbinose	1.4
5	RIB	D-RIBose	1.4
6	DXYL	D-XYLose	1.4
7	LXYL	L-XYLose	1.4
8	ADO	D-ADOnitol	1.36
9	MDX	Methyl-βD-Xylopyranoside	1.28

Strip 10 - 19

Tube	Test	Active ingredients	QTY (mg/cup.)
10	GAL	D-GALactose	1.4
11	GLU	D-GLUcose	1.56
12	FRU	D-FRUctose	1.4
13	MNE	D-MaNnosE	1.4
14	SBE	L-SorBosE	1.4
15	RHA	L-RHAmnose	1.36
16	DUL	DULcitol	1.36
17	INO	INOsitol	1.4
18	MAN	D-MANnitol	1.36
19	SOR	D-SORbitol	1.36

Strip 20 - 29

Tube	Test	Active ingredients	QTY (mg/cup.)
20	MDM	Methyl-αD-Mannopyranoside	1.28
21	MDG	Methyl-αD-Glucopyranoside	1.28
22	NAG	N-AcetylGlucosamine	1.28
23	AMY	AMYgdalin	1.08
24	ARB	ARButin	1.08
25	ESC	ESCUlin ferric citrate	1.16 0.152
26	SAL	SALicin	1.04
27	CEL	D-CELlobiose	1.32
28	MAL	D-MALtose	1.4
29	LAC	D-LACTose (bovine origin)	1.4

Strip 30 - 39

Tube	Test	Active ingredients	QTY (mg/cup.)
30	MEL	D-MELibiose	1.32
31	SAC	D-SACcharose (sucrose)	1.32
32	TRE	D-TREhalose	1.32
33	INU	INUlin	1.28
34	MLZ	D-MeLeZitose	1.32
35	RAF	D-RAFfinose	1.56
36	AMD	AmiDon (starch)	1.28
37	GLYG	GLYcoGen	1.28
38	XLT	XyLiToI	1.4
39	GEN	GENtiobiose	0.5

Strip 40 - 49

Tube	Test	Active ingredients	QTY (mg/cup.)
40	TUR	D-TURanose	1.32
41	LYX	D-LYXose	1.4
42	TAG	D-TAGatose	1.4
43	DFUC	D-FUCose	1.28
44	LFUC	L-FUCose	1.28
45	DARL	D-ARAbitoL	1.4
46	LARL	L-ARAbitoL	1.4
47	GNT	potassium GlucoNaTe	1.84
48	2KG	potassium 2-KetoGluconate	2.12
49	5KG	potassium 5-KetoGluconate	1.8

The quantities indicated may be adjusted depending on the titer of the raw materials used.

REAGENTS AND MATERIAL REQUIRED BUT NOT PROVIDED

Reagents :

- Inoculation medium :
 - API 50 CHL Medium (Ref. 50 410)
 - API 50 CHB/E Medium (Ref. 50 430)
 - (+ products mentioned in the package inserts of these media)
 - or any other suitable medium
- Mineral oil (Ref. 70 100)
- McFarland Standard (Ref. 70 900) or DENSIMAT (Ref. 99 234) or ATB® Densitometer
- Identification software (consult bioMérieux)

Material :

- Pipettes or PSlpettes
- Ampule rack
- Ampule protector
- Swabs
- General microbiology laboratory equipment

WARNINGS AND PRECAUTIONS

- **For *in vitro* diagnostic use and microbiological control.**
- **For professional use only.**
- This kit contains products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not totally guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious, and handled observing the usual safety precautions (do not ingest or inhale).
- All specimens, microbial cultures and inoculated products should be considered infectious and handled appropriately. Aseptic technique and usual precautions for handling the bacterial group studied should be observed throughout this procedure. Refer to "NCCLS M29-A, *Protection of Laboratory Workers from Instrument Biohazards and Infectious Disease Transmitted by Blood, Body Fluids, and Tissue; Approved Guideline* - December 1997". For additional handling precautions, refer to "Biosafety in Microbiological and Biomedical Laboratories, HHS Publication No. (CDC) 93-8395, 3rd Edition (May 1993)", or to the regulations currently in use in each country.
- Do not use reagents past the expiration date.
- Before use, check that the packaging of the various components is intact.
- Do not use strips which have been damaged : cupules deformed, desiccant sachet open, ...
- The performance data were obtained using the procedure indicated in this package insert. Any change or modification in the procedure may affect the results.
- Interpretation of the test results should be made taking into consideration the patient history, the source of the specimen, colonial and microscopic morphology of the strain and, if necessary, the results of any other tests performed, particularly the antimicrobial susceptibility patterns.

STORAGE CONDITIONS

The strips should be stored at 2-8°C until the expiration date indicated on the packaging.

SPECIMENS (COLLECTION AND PREPARATION)

API 50 CH is not for use directly with clinical or other specimens.

The microorganisms to be identified must first be isolated on a suitable culture medium according to standard microbiological techniques.

INSTRUCTIONS FOR USE

Depending on which medium is used, API 50 CHL Medium or API 50 CHB/E Medium, carefully read the corresponding package insert.

Preparation of the strips

Each strip is made up of 5 smaller strips each containing 10 numbered tubes.

- Prepare an incubation box (tray and lid).
- Record the reference of the strain on the elongated flap of the tray. (Do not record the reference on the lid as it may be misplaced during the procedure.)
- Distribute about 10 ml of distilled water or demineralized water [or any water without additives or chemicals which may release gases (e.g. Cl₂, CO₂, etc.)] into the honeycombed wells of the tray to create a humid atmosphere.
- Remove the 2 strips (0-19 and 20-39) from their packaging, separate into 4 smaller strips (0-9, 10-19, 20-29 and 30-39) and place all 4 in the incubation tray.
- Take the remaining smaller strip (40-49) out of the packaging and place it next to the others in the incubation tray to complete the strip.

Preparation of the inoculum

- Culture the microorganism on a medium adapted to its growth.
- Check the purity of the strain.
- Harvest all the bacteria from a solid medium using a swab or from a liquid medium using centrifugation.
- Prepare the inoculum in the appropriate medium (see the API 50 CHL Medium and API 50 CHB/E Medium package inserts).
This suspension must be used immediately after preparation.

Inoculation of the strips

Distribute the bacterial suspension using a sterile pipette into the 50 tubes as follows :

- Tilt the incubation box slightly forwards.
- Avoid the formation of bubbles by placing the tip of the pipette against the side of the cupule.
- When only the tube is to be inoculated, do not exceed the top of the tube so as to maintain anaerobic conditions.
- When the tube and the cupule are to be completely filled, avoid the formation of a concave or convex meniscus.
- Incubate the strips at the optimum temperature for growth of the group of microorganisms being tested : 30°C, 37°C or 55°C.

READING AND INTERPRETATION

Reading the strips

The strips are read after the stipulated incubation times (e.g., 24 hrs., 48 hrs.), depending on the microorganism and the type of reaction studied.

Interpretation

Interpret each test (positive (+), negative (-), doubtful (?)) and record the results on the result sheet.

The results form a biochemical profile which, when entered in the identification software, provides the identification of *Lactobacillus* and related genera or *Bacillus* and related genera, *Enterobacteriaceae* and *Vibrionaceae*.

NOTE : The results may be used for other purposes :

- Epidemiological grouping into types of the microorganism.
- Taxonomical analysis of a group of microorganisms.
- Classification of an unknown bacterial population into homogeneous groups.

QUALITY CONTROL

The strips are systematically controlled at various stages of their manufacture. For those users who wish to perform their own quality control tests with the strip, the following strains may be used :

For *Lactobacillus* : ***Lactobacillus paracasei ssp paracasei* NCFB 206** or **ATCC BAA-52** (with API 50 CHL Medium)

	0	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49
24	-	-	-	-	+	-	-	-	-	+	+	+	+	-	-	V	-	+	V	-	+	+	-	V	+	V	+	+	-	-	+	+	+	+	-	-	-	-	-	-	+	-	+	-	-	-	-	V	-	-
48	-	-	-	-	+	-	-	-	-	+	+	+	+	-	-	+	-	+	+	+	+	+	V	+	+	+	+	+	-	-	+	+	+	+	-	-	-	-	V	+	-	+	-	-	-	-	+	-	-	

For *Bacillus* : ***Bacillus polymyxa* (*) ATCC 43865** (with API 50 CHB/E Medium)

	0	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49			
24	-	+	-	-	+	+	+	-	-	V	+	+	+	+	-	V	-	-	+	-	-	+	-	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+		
48	-	+	-	-	+	+	+	-	-	+	+	+	+	+	-	V	-	-	+	-	-	+	-	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+

(*) *Bacillus polymyxa* identified as *Paenibacillus polymyxa* with API 50 CH and API 50 CHB/E Medium.

Results obtained after incubation at 30°C.

For *Enterobacteriaceae* : ***Klebsiella pneumoniae ssp pneumoniae* ATCC 35657** (with API 50 CHB/E Medium)

	0	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49		
24	-	+	-	-	+	+	+	-	+	-	+	+	+	+	-	+	-	+	+	+	+	-	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	
48	-	+	-	V	+	+	+	-	+	-	+	+	+	+	-	+	-	+	+	+	+	-	+	+	V	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+

ATCC : American Type Culture Collection, 10801 University Boulevard, Manassas, VA 20110-2209, USA.

NCFB : National Collection for Food Bacteria (=NCDO), Institute of Food Research, Reading Laboratory, Earley Gate, Reading RG6 6BZ, ENGLAND

It is the responsibility of the user to perform Quality Control in accordance with any local applicable regulations.

LIMITATIONS OF THE METHOD

- The user assumes full responsibility for the identification of any species not included in the API 50 CHL and API 50 CHB/E databases.
- Only pure cultures of a single organism should be used.

RANGE OF EXPECTED RESULTS

Consult the Identification Tables in the package inserts of the media associated with this strip for the range of expected results for the various biochemical reactions.

PERFORMANCE

Refer to the Performance section in the package inserts of the media associated with this strip.


WASTE DISPOSAL

It is the responsibility of each laboratory to handle waste and effluents produced according to their type and degree of hazardousness and to treat and dispose of them (or have them treated and disposed of) in accordance with any applicable regulations.

WARRANTY

bioMérieux disclaims all warranties, express or implied, including any implied warranties of MERCHANTABILITY AND FITNESS FOR A PARTICULAR USE. bioMérieux shall not be liable for any incidental or consequential damages. IN NO EVENT SHALL BIOMERIEUX'S LIABILITY TO CUSTOMER UNDER ANY CLAIM EXCEED A REFUND OF THE AMOUNT PAID TO BIOMERIEUX FOR THE PRODUCT OR SERVICE WHICH IS THE SUBJECT OF THE CLAIM.

LITERATURE REFERENCES p. I
 INDEX OF SYMBOLS p. II

bioMérieux® sa
 au capital de 11 879 045 €
 673 620 399 RCS LYON
 69280 Marcy-l'Etoile / France
 Tel. 33 (0)4 78 87 20 00
 Fax 33 (0)4 78 87 20 90
<http://www.biomerieux.com>

bioMérieux, Inc
 Box 15969,
 Durham, NC 27704-0969 / USA
 Tel. (1) 919 620 20 00
 Fax (1) 919 620 22 11

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The logo is a registered and protected trademark of bioMérieux sa or one of its subsidiaries