NBB[®]-A, NBB[®]-B, NBB[®]-C Manual

WE BRING
IDEAS TO LIFE.

NATURAL INGREDIENTS
INGREDIENT SYSTEMS
INTEGRATED SOLUTIONS

Table of Contents

	Page
1. Introduction: Microbial quality control in breweries with NBB®-A, NBB®-B and NBB®-C	5
2. NBB®-A, NBB®-B and NBB®-C: Packaging formats	7
3. Storage	8
4. Intended purpose	8
5. Safety requirements	8
6. Quality control	9
7. Your benefits with the ready-to-use media: NBB®-A, NBB®-B and NBB®-C	9
8. Materials needed for implementation (not included)	9
9. Important information	13
10. Protocols NBB®-A Agar	13
10.1. NBB®-A Agar: Preparation of agar plates	14
10.1.1. Liquefaction of NBB®-A Agar	14
10.1.2. Plate pouring with NBB®-A Agar	15
10.2. NBB®-A Agar: Detection of beer spoiling bacteria with membrane filtration	16
10.2.1. Membrane filtration of clear beer samples	16
10.2.2. Incubation of NBB®-A Agar plates	17
10.2.3. Evaluation of NBB®-A Agar plates with membrane filter	17
10.3. NBB®-A Agar: Detection of beer spoiling bacteria using plating or spreading method	18
10.3.1. Plating or spreading on NBB®-A Agar	18
10.3.2. Incubation of NBB®-A Agar plates	19
10.3.3. Evaluation of NBB®-A Agar plates	19
10.4. NBB®-A Agar: Detection of beer spoiling bacteria from air samples	20
10.4.1. Air sampling with NBB®-A Agar	20
10.4.2. Incubation of NBB®-A Agar plates	20
10.4.3. Evaluation of NBB®-A Agar plates	20
10.5. NBB®-A Agar: Detection of beer spoiling bacteria from CO ₂ and compressed air	21
10.5.1. CO ₂ and compressed air sampling with NBB®-A Agar	21
10.5.2. Membrane filtration	21
10.5.3. Incubation of NBB®-A Agar plates	22
10.5.4. Evaluation of NBB®-A Agar plates with membrane filter	22
10.6. NBB®-A Agar: Detection of beer spoiling bacteria using the pour-plate method	23
10.6.1. Liquefaction of NBB®-A Agar	23
10.6.2. Pouring plates with NBB®-A Agar	24
10.6.3. Incubation of NBB®-A Agar pour plates	25
10.6.4. Evaluation of NBB®-A Agar pour plates	25

	Page
11. Protocols NBB[®]-B Broth	25
11.1. NBB[®]-B Broth: Detection of beer spoiling bacteria from yeast samples in tubes	25
11.1.1. Preparation of NBB [®] -B Broth tubes	25
11.1.2. Yeast sampling in NBB [®] -B Broth tubes	26
11.1.3. Incubation of NBB [®] -B Broth tubes	27
11.1.4. Evaluation of NBB [®] -B Broth tubes	27
11.2. NBB[®]-B Broth: Detection of beer spoiling bacteria from yeast samples in swing top bottles	27
11.2.1. Addition of yeast into swing top bottles	27
11.2.2. Addition of NBB [®] -B Broth in swing top bottles	28
11.2.3. Incubation of NBB [®] -B Broth swing top bottles	28
11.2.4. Evaluation of NBB [®] -B Broth swing top bottles	28
11.3. NBB[®]-B Broth: Detection of beer spoiling bacteria from CO₂ and compressed air	29
11.3.1. Preparation of NBB [®] -B Broth tubes	29
11.3.2. CO ₂ and compressed air sampling with NBB [®] -B Broth	29
11.3.3. Incubation of NBB [®] -B Broth tubes	30
11.3.4. Evaluation of NBB [®] -B Broth tubes	30
12. Protocols NBB[®]-C Concentrate	31
12.1. NBB[®]-C Concentrate: Enrichment of beer spoiling bacteria from turbid beer samples	31
12.1.1. NBB [®] -C Concentrate: Preparation of swing top bottles	31
12.1.2. Initiating enrichment with NBB [®] -C Concentrate in swing top bottles or beer bottles	32
12.1.3. Incubation of NBB [®] -C Concentrate	33
12.1.4. Evaluation of NBB [®] -C Concentrate	34
12.2. NBB[®]-C Concentrate: Enrichment of beer spoiling bacteria from wort samples	35
12.2.1. NBB [®] -C Concentrate: Preparation of swing top bottles	35
12.2.2. NBB [®] -C Concentrate: Enrichment of wort samples in swing top bottles	35
12.2.3. Incubation of NBB [®] -C Concentrate	36
12.2.4. Evaluation of NBB [®] -C Concentrate	36
12.3. NBB[®]-C Concentrate: Enrichment of beer spoiling bacteria from yeast samples (alternative protocol)	37
12.3.1. NBB [®] -C Concentrate: Preparation of sample enrichment	37
12.3.2. Incubation of NBB [®] -C Concentrate with yeast samples	37
12.3.3. Evaluation of NBB [®] -C Concentrate with yeast samples	38

	Page
13. FAQs on application and evaluating the results	39
13.1. NBB [®] -A Agar	39
13.2. NBB [®] -B Broth	40
13.3. NBB [®] -C Concentrate	40
14. Appendix:	42
14.1. Brief overview: Detection of beer spoiling bacteria using NBB [®] -A Agar	42
14.2. Brief overview: Detection of beer spoiling bacteria using NBB [®] -B Broth	43
14.3. Brief overview: Detection of beer spoiling bacteria using NBB [®] -C Concentrate	43
14.4. Glossary (in alphabetical order)	44
15. Supply Sources	45
16. Information	46
16.1. Trademarks	46
16.2. Distribution	46
16.3. Recommended further reading	46

1. Introduction: Microbial quality control in breweries with NBB[®]-A, NBB[®]-B and NBB[®]-C

Producing microbiologically immaculate beer is the key role of a brewery. However, beer spoiling bacteria can impair fermentation or have a negative impact on the beer's flavour. Microbiological control of the brewing process is therefore extremely important and critical. The **NBB[®]** products make it possible to carry out microbiological control on all samples of the brewing process in a simple, fast and comprehensive way (see Figure 1).

Gram-positive lactic acid bacteria like *Lactobacillus brevis*, *Lactobacillus lindneri* and *Pediococcus damnosus* are responsible for 70–80% of contamination in breweries. In addition, gram-negative bacteria like *Pectinatus* spp. and *Megasphaera* spp. cause around 3–7% of microbial contamination and, together with lactic acid bacteria, form the most significant group of obligate beer spoiling bacteria.

With **NBB[®]-A** Agar, beer spoiling bacteria can be quantitatively detected in rinsing water, in beer, in CO₂ or compressed air systems and in the ambient air. With **NBB[®]-B** Broth, the brewing yeast can be qualitatively analysed for microbiological contamination, allowing beer spoiling bacteria to be spotted before the start of fermentation. The microbiological detection with **NBB[®]-A** Agar and **NBB[®]-B** Broth can be easily and reliably seen through a change in the colour from red to yellow, even without any previous knowledge.

NBB[®]-C Concentrate is used to enrich beer spoiling bacteria from turbid beer samples. Additionally, these bacteria can be enriched by wort samples by adding sterile beer and **NBB[®]-C** Concentrate. The samples can be analysed with a microscope after incubation. The combination of **NBB[®]-A** Agar, **NBB[®]-B** Broth and **NBB[®]-C** Concentrate enables complete microbiological control of the brewing process and can even detect microbial trace contaminations. This reduces the risk of beers with microbial impairment and enhances efficiency in the brewery.

The **NBB[®]** products have been used successfully in quality assurance in breweries for 30 years and have established themselves in various collections of standard methods, such as MEBAK and EBC.


Figure 1: Overview of the microbial analysis options for the brewing process with the **NBB®** product portfolio.

Overview of the detection options for beer spoiling bacteria with the **NBB®** product portfolio:

Labelling	Sample types	Methods	NBB® product
1 2	Clear beer samples/rinsing water	Membrane filtration & pour-plate method	NBB®-A Agar
3	Yeast cultures	Spreading or plating on agar	NBB®-A Agar
3	Air	Airborne bacteria collection on agar	NBB®-A Agar
3 2	CO ₂ and compressed air	Blowing into sterile water and membrane filtration on agar or direct gassing in broth	NBB®-A Agar NBB®-B Broth
1 2	Yeast samples	Enrichment with broth or enrichment with concentrate	NBB®-B Broth NBB®-C Concentrate
3 4	Turbid beer samples	Enrichment with concentrate	NBB®-C Concentrate
1	Wort samples	Enrichment with concentrate (sterile beer addition)	NBB®-C Concentrate

2. NBB[®]-A, NBB[®]-B and NBB[®]-C: Packaging formats

2.04709.782

2.04723.646 / 2.04710.782

2.04711.782

Art. no.	Döhler Microsafety Design [®] product	Description	Target bacteria	Packaging unit
2.04709.782	NBB[®]-A Agar (bottle)	For the analysis of clear beer samples (membrane filtration), pour-plate method, plating & swab processes, CO ₂ and compressed air (first blowing into sterile water). Ready-to-use NBB[®]-A Agar for the detection of beer spoiling bacteria through colour change.	Beer spoiling bacteria: e.g. <i>Lactobacillus</i> spp. and <i>Pediococcus</i> spp. as well as <i>Pectinatus</i> spp. and <i>Megasphaera</i> spp.	1 cardboard box (9x 250 ml NBB[®]-A Agar bottle)
2.04723.646	NBB[®]-B Broth (tubes)	For the analysis of pitching and harvesting yeast as well as CO ₂ and compressed air. Ready-to-use NBB[®]-B Broth		1 polystyrene box (20x 10 ml NBB[®]-B Broth tubes)
2.04710.782	NBB[®]-B Broth (bottle)	for the detection of beer spoiling bacteria through colour change.		1 cardboard box (9x 250 ml NBB[®]-B Broth bottle)
2.04711.782	NBB[®]-C Concentrate (bottle)	For the enrichment of (yeast-) turbid beer samples. Ready-to-use NBB[®]-C Concentrate for the detection of beer spoiling bacteria through microscopy.		1 cardboard box (9x 250 ml NBB[®]-C Concentrate bottle)

3. Storage

NBB®-A Agar, **NBB®-B** Broth and **NBB®-C** Concentrate must be stored in a dark, dry place at a temperature of +4 to +8° C. The products must not be frozen. The expiration date can be found on the label on the packaging (see product label).

4. Intended purpose

NBB®-A Agar, **NBB®-B** Broth and **NBB®-C** Concentrate are intended for the detection of beer spoiling bacteria in breweries.

When using the products, we recommend that you pay attention to the safety requirements outlined below and work particularly cautiously and hygienically to prevent secondary contamination during sampling.

5. Safety requirements

The basic rules for microbiological operations shall be observed when using this product. This includes using clean equipment such as lab coats, safety goggles and gloves. Such equipment not only contributes to your own personal safety, but also prevents any secondary contaminations being passed on by the user. The safety data sheets can be found at www.doehler-dmd.com.

NBB®-A Agar, **NBB®-B** Broth and **NBB®-C** Concentrate are not suitable for consumption.

Special care is required when using Bunsen burners while wearing latex gloves. Always maintain a sufficiently safe distance from the flame. If the gloves catch fire, it can result in extremely severe injuries.

Fire hazard: Never use alcohol for sterilisation or decontamination purposes while working with an open flame (Bunsen burner).

6. Quality control

The quality of **NBB®-A** Agar, **NBB®-B** Broth and **NBB®-C** Concentrate is carefully tested. In particular, the functionality of the culture media is specifically tested and secured using beer spoiling bacteria and regularly checked at the Research Center Weihenstephan at Technical University Munich. You can obtain a certificate of analysis from Doehler GmbH.

7. Your benefits with the ready-to-use media: **NBB®-A**, **NBB®-B** and **NBB®-C**

NBB®-A Agar, **NBB®-B** Broth and **NBB®-C** Concentrate are ready-to-use culture media, eliminating the need for extensive and costly preparations such as weighing and dissolving chemicals, in addition to adjusting the pH values and conducting subsequent autoclaving processes. Quality control has also already been conducted on the culture media's function. This saves you both time and money so you can focus on your microbiological quality control straight away.

8. Materials needed for implementation (not included)

When using **NBB®-A** Agar, **NBB®-B** Broth and **NBB®-C** Concentrate, the following items are recommended for carrying out the experiments:

General laboratory items	Use
Mobile Bunsen burner	Bunsen burner flames create a sterile work area (approx. 50 cm around the flame) in which microbiological experiments can be carried out. The heat kills off any bacteria in the immediate atmosphere. The Bunsen burner is also a useful tool for collecting sterile samples.

General laboratory items	Use
A sterile (serological) pipette (1 – 10 ml) with pipetting aid (e.g. Peleus ball or pipetting aid made of plastic)	<p>(Serological) pipettes are suitable for transferring sterile liquids.</p> <p> Warning: Never pipette fluids with your mouth! Always use a pipetting aid!</p>
Incubator with thermostat	<p>For incubating samples at defined temperatures.</p>
Disposable gloves	<p>Disposable gloves to prevent secondary contamination by the user.</p>
Sterile workbench	<p>Provides a sterile working environment. Essential lab equipment for sterile work.</p>
Laboratory vessels (e. g. Duran® glass bottles)	<p>Inert and autoclavable vessels for sampling.</p>
Anaerobiosis jar (incl. CO₂ gassing)	<p>For incubating samples under anaerobic conditions.</p>
Autoclave	<p>Standard laboratory device for sterilising laboratory materials.</p>
70 % alcohol solution	<p>For decontaminating surfaces and devices. Composed of 30 % distilled, deionised water and 70 % ethanol.</p>

Laboratory items for NBB [®] -A Agar	Use
Water bath	Water bath with temperature control for gentle liquefaction of NBB[®]-A Agar.
Air sampler (as needed)	Airborne microorganism collector that can be equipped with Petri dishes. It should be possible to adjust the air volume variably.
Membrane filtration system	For the filtration of liquids with sterile membrane filters with different pore sizes (e. g. 0.45 µm).
Sterile tweezers	Stainless steel tweezers that can be autoclaved in glass vessels. For transferring membrane filters.
Drigalski spatula or inoculation loop	Disposable or reusable laboratory items for plating or spreading microorganisms. Reusable spatulas (stainless steel) or loops (platinum wire) must be sterilised by flaming before use.
Sterile Petri dishes with vents	Sterile Petri dishes can be obtained from specialist laboratory suppliers. The vents enable gas exchange.

Laboratory items for NBB [®] -B Broth	Use
Sterile swab	Sterile (wooden) swab with cotton tip, usually available sterile and individually-packaged (see Chapter 15, p. 45)
Transparent swing top bottle	Swing top bottles are available in various volumes. 50 ml swing top bottles are needed for NBB[®]-B Broth. Transparent bottles are needed in order to recognise the colour change. Stainless steel wire should be used for autoclaving. (Procurement source: specialist laboratory suppliers e.g. VWR)

Laboratory items for NBB [®] -B Broth	Use
Screw-top test tubes	Inert and autoclavable vessels available in various volumes. The screw top allows loose to tight sealing. Not required for use with NBB[®]-B-Am Broth tubes (2.04723.646).
Test tube racks	NBB[®]-B Broth tube holder

Laboratory items for NBB [®] -C Concentrate	Use
Transparent swing top bottles	Swing top bottles are available in various volumes. 200 ml swing top bottles are needed for NBB[®]-C Concentrate. Stainless swing tops should be used for autoclaving. (Procurement source: specialist laboratory suppliers e.g. VWR)
Shaker	Laboratory device whose rotational frequency can be controlled to shake microbial culture vessels.
Erlenmeyer flasks with cotton stoppers	Inert and autoclavable laboratory vessels with a capacity of at least 200 ml.

9. Important information

Clean microbiological work is the basis of all analysis. In order to obtain reliable results, we strongly recommend following the safety instructions and the work steps below (see Protocols, from p. 13).

NBB®-A Agar, **NBB®-B Broth** and **NBB®-C Concentrate** are not suitable for consumption.

Following analysis, the incubated culture media should be deactivated before being disposed of. Autoclaving the culture media or disposing of them as hazardous waste is recommended for deactivating the bacteria.

10. Protocols NBB®-A Agar

DMD® product	Preparation	Application	Frequency	Protocol
NBB®-A Agar	Preparation of agar plates (see Protocol 10.1.)	Membrane filtration of clear beer samples or water	Regular (every production cycle)	10.2.
		Spreading and plating method	As needed	10.3.
		Airborne microorganism collection	Monthly/as needed	10.4.
		CO ₂ and compressed air sampling	Occasional (for analysis of contamination source)	10.5.
	No preparation	Pour-plate method	As needed	10.6.

10.1. NBB[®]-A Agar: Preparation of agar plates

10.1.1. Liquefaction of NBB[®]-A Agar

Liquefy an appropriate number of **NBB[®]-A Agar** bottles for your tests in a water bath at approx. 95 °C.

Liquefaction of a bottle of **NBB[®]-A Agar** takes approx. 40–45 min.

Allow the **NBB[®]-A Agar** to cool to approx. 45–50 °C, so that you can touch the bottles without danger of being burned. Make sure that the agar does not solidify again as it cools.

Warning:

Never keep the liquefied NBB[®]-A Agar in a water bath at 45 °C for more than 4 h. This damages the agar structure and changes its solidification properties irrevocably. Therefore, avoid long liquefaction or multiple liquefactions at 95 °C, as the growth components could be damaged otherwise.

Do not use a microwave to liquefy NBB[®]-A Agar. Danger of fire due to the bottle's metal lid and local overheating also destroys some ingredients.

10.1.2. Plate pouring with NBB[®]-A Agar

Take the bottles, now cooled to 45 °C, into a sterile workbench and arrange the required number of empty Petri dishes with vents. Consider the diameter of the Petri dishes needed for the tests.

Homogenise the **NBB[®]-A** Agar by gently swirling the bottle. Pour into the **NBB[®]-A** Agar plates as required. If possible, pour a whole bottle of **NBB[®]-A** Agar into Petri dishes, as heating the **NBB[®]-A** Agar bottle repeatedly damages the ingredients and leads to poorer results.

Petri dish diameter	Fill volume with NBB [®] -A Agar	Number of plates from one bottle of NBB [®] -A Agar
5.5 cm	10 ml	25 pc.
9 cm	31 ml	8 pc.
14 cm	62 ml	4 pc.

Poured **NBB[®]-A** Agar plates should be stacked for cooling in order to prevent the formation of condensation.

NBB[®]-A Agar plates can be stored in a dark place at a temperature of +4 to +8 °C for approx. 1 week. Storing the plates upside down is recommended.

The **NBB[®]-A** Agar plates should be dried in a sterile workbench before use, so that any condensation is removed.

10.2. NBB[®]-A Agar: Detection of beer spoiling bacteria with membrane filtration

10.2.1. Membrane filtration of clear beer samples

Define a suitable volume of a clear beer sample that you want to examine.

Use your membrane filtration system to filter the clear beer samples on a suitable membrane filter. Please note the manufacturer's instructions on the membrane filtration system. For strictly anaerobic beer spoiling bacteria (e.g. *Pectinatus* spp.), the filter should be briefly flushed with CO₂ in order to remove oxygen from the filter.

Label your **NBB[®]-A** Agar plates.


Then transfer the membrane filter to an **NBB[®]-A** Agar plate with sterile tweezers. Place the membrane filter on the surface of the **NBB[®]-A** Agar plate, taking care to avoid air bubbles.

10.2.2. Incubation of NBB[®]-A Agar plates

Incubate the **NBB[®]-A** Agar plates with the membrane filter under anaerobic conditions (e.g. in an anaerobiosis jar) for 5 – 7 days at $27 \pm 2^\circ\text{C}$ in an incubator.

The incubation time will need to be extended for some slow-growing beer spoiling bacteria (e.g.: *Lactobacillus lindneri*).

10.2.3. Evaluation of NBB[®]-A Agar plates with membrane filter

A change in colour from red to yellow shows the growth of typical beer spoiling bacteria on **NBB[®]-A** Agar plates with membrane filters.

Further analysis techniques, such as microscopy and PCR analysis, can be used to identify the bacteria.

10.3. NBB[®]-A Agar: Detection of beer spoiling bacteria using plating or spreading method

10.3.1. Plating or spreading on NBB[®]-A Agar

Label the **NBB[®]-A** Agar plates for this test.

A) Method: plating

Pipette max. 150 µl of your sample onto the **NBB[®]-A** Agar plate in a sterile workbench. Use a Drigalski spatula to plate the sample evenly over the agar.


B) Method: spreading (e.g. cultures)

Take a sterile inoculation loop and immerse it once briefly in a homogeneous (yeast) culture.

Use this inoculation loop to spread the microorganisms onto an **NBB[®]-A** Agar plate.

With the help of the spreading pattern below, you can detect microorganisms individually, i.e. individual colonies. Use the same inoculation loop for all spreading lines. Do not immerse it in the culture again. If using reusable inoculation loops, flame between spreading – or use a new disposable inoculation loop.

Warning:


Always use a new sterile inoculation loop for every new culture or NBB[®]-A Agar plate, or flame the reusable inoculation loop, in order to prevent secondary contamination.

10.3.2. Incubation of NBB[®]-A Agar plates

Incubate the **NBB[®]-A** Agar plates under anaerobic conditions (e.g. in an anaerobiosis jar) for 5–7 days at 27 ± 2 °C in an incubator.

The incubation time will need to be extended for some slow-growing beer spoiling bacteria (e.g.: *Lactobacillus lindneri*).

10.3.3. Evaluation of NBB[®]-A Agar plates

A change in colour from red to yellow shows the growth of typical beer spoiling bacteria on **NBB[®]-A** Agar plates (see image on left: plated samples with negative and positive results).

Further analysis techniques, such as microscopy and PCR analysis, can be used to identify the bacteria.

10.4. NBB[®]-A Agar: Detection of beer spoiling bacteria from air samples

10.4.1. Air sampling with NBB[®]-A Agar

Determine which agar plate diameters can be used with your air sampler. Various air samplers for agar plates are available for purchase.

Label suitable **NBB[®]-A** Agar plates for this test.


Use your air sampler in accordance with the manufacturer's instructions and collect the required volume of ambient air on the plate.

10.4.2. Incubation of NBB[®]-A Agar plates

Incubate the **NBB[®]-A** Agar plates under anaerobic conditions (e.g. in an anaerobiosis jar) for 5 days at 27 ± 2 °C in an incubator.

If you incubate for more than 5 days, unspecific microorganisms, i.e. bacteria that are not obligately beer spoiling, will also be enriched. Therefore, always analyse the **NBB[®]-A** Agar plates after 5 days.

10.4.3. Evaluation of NBB[®]-A Agar plates

A change in colour from red to yellow shows the growth of typical beer spoiling bacteria on **NBB[®]-A** Agar plates.

Further analysis techniques, such as microscopy and PCR analysis, can be used to identify the bacteria.

10.5. NBB[®]-A Agar: Detection of beer spoiling bacteria from CO₂ and compressed air

10.5.1. CO₂ and compressed air sampling with NBB[®]-A Agar

Fill a sterile bottle or a sterile vessel with ca. 50 ml of sterile water.

Decontaminate the outlet of the compressed air or CO₂ system with a 70% alcohol solution.

Then immerse the outlet in the sterile water. Allow compressed air or CO₂ to flow slowly into the water for around 15 seconds, preventing gas from shooting out. Most microorganisms are found in the back of the sampling tap.

10.5.2. Membrane filtration

Use your membrane filtration system to filter the gassed water from 10.5.1. onto a suitable membrane filter. Please note the manufacturer's instructions on the membrane filtration system. For strictly anaerobic beer spoiling bacteria (e.g. *Pectinatus* spp.), the membrane filter should be briefly flushed with CO₂ in order to remove oxygen from the filter.

Label the **NBB[®]-A** Agar plates.

Transfer the membrane filter to an **NBB[®]-A** Agar plate with sterile tweezers. Place the membrane filter on the surface of the **NBB[®]-A** Agar plate, taking care to avoid air bubbles.


10.5.3. Incubation of NBB[®]-A Agar plates

Incubate the **NBB[®]-A** Agar plates with the membrane filter under anaerobic conditions (e.g. in an anaerobiosis jar) for 5–7 days at 27 ± 2 °C in an incubator.

The incubation time will need to be extended for some slow-growing beer spoiling bacteria (e.g.: *Lactobacillus lindneri*).

10.5.4. Evaluation of NBB[®]-A Agar plates with membrane filter

A change in colour from red to yellow shows the growth of typical beer spoiling bacteria on **NBB[®]-A** Agar plates with membrane filters.

Further analysis techniques, such as microscopy and PCR analysis, can be used to identify the bacteria.

10.6. NBB[®]-A Agar: Detection of beer spoiling bacteria using the pour-plate method

10.6.1. Liquefaction of NBB[®]-A Agar

Liquefy an appropriate number of **NBB[®]-A Agar** for your tests in a water bath at approx. 95 °C. Liquefaction of a bottle of **NBB[®]-A Agar** takes approx. 40 – 45 min.

Allow the **NBB[®]-A Agar** to cool to approx. 45 – 50 °C, so that you can touch the bottles without danger of being burned. Make sure that the agar does not solidify again as it cools.

Warning:

Never keep the liquefied NBB[®]-A Agar in a water bath at 45 °C for more than 4 h. This damages the agar structure and changes its solidification properties irrevocably. Therefore, avoid long liquefaction or multiple liquefactions at 95 °C, as the growth components could be damaged otherwise.

Do not use a microwave to liquefy NBB[®]-A Agar. Danger of fire due to the bottle's metal lid and local overheating also destroys some ingredients.

10.6.2. Pouring plates with NBB[®]-A Agar

Take the bottles, now cooled to 45 °C, into a sterile workbench and arrange sufficient empty Petri dishes with vents. Make sure that the liquefied agar is not too hot, as bacteria can be killed by excessive agar temperatures.

Label the Petri dishes.

Use a pipette to add the beer samples to the Petri dishes – do not exceed the maximum sample volumes given below for each Petri dish diameter.

Petri dish diameter	Max. sample volume
9 cm	5 ml
14 cm	10 ml

Homogenise the **NBB[®]-A Agar** by gently swirling the bottle. Fill up the sample volume with **NBB[®]-A Agar**. Use the Petri dish diameter as a guide for the volume of **NBB[®]-A Agar** to be added.

Petri dish diameter	Fill volume with NBB [®] -A Agar	Number of plates from one bottle of NBB [®] -A Agar
5.5cm	10 ml	25 pc.
9 cm	31 ml	8 pc.
14 cm	62 ml	4 pc.

Mix the **NBB[®]-A Agar** pour plates by rotating the Petri dish gently multiple times (in a figure eight). Do not allow the agar to splash onto the lid or run out of the plate. If this does happen, dispose of the plate affected and start again.


Poured **NBB[®]-A Agar** pour plates should be stacked for cooling in order to prevent the formation of condensation.

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A change in colour from red to yellow shows the growth of typical beer spoiling bacteria on **NBB®-A** Agar pour plates.

Further analysis techniques, such as microscopy and PCR analysis, can be used to identify the bacteria.

11. Protocols NBB®-B Broth

DMD® product	Application	Frequency	Protocol
NBB®-B Broth	Yeast sampling	Regular (every production cycle: e.g. after yeast harvesting or before adding yeast)	11.1. & 11.2.
	CO ₂ and compressed air sampling	Monthly/as needed	11.3.

11.1. NBB®-B Broth: Detection of beer spoiling bacteria from yeast samples in tubes

11.1.1. Preparation of NBB®-B Broth tubes

Fill an appropriate test tube to approx. 70% with **NBB®-B** Broth in a sterile workbench using a sterile serological pipette, or use pre-filled **NBB®-B** Broth tubes (art. no. 2.04723.646).

11.1.2. Yeast sampling in NBB[®]-B Broth tubes

A) Low-viscosity yeast suspension

Extract 0.1 – 0.3 ml from a pipettable yeast suspension using a sterile pipette.

Make sure that only the tip of the pipette comes into contact with the yeast suspension in order to prevent secondary contamination.

Then add the yeast suspension to the prepared **NBB[®]-B** Broth tube from Step 11.1.1.

Do not seal the **NBB[®]-B** Broth tube completely to allow CO₂ to diffuse during incubation.

Label the tube.

B) High-viscosity yeast suspension

Extract ca. 0.3 ml of the high-viscosity yeast suspension using a sterile swab.

Transfer the swab to the prepared **NBB[®]-B** Broth tube from Step 11.1.1. Detach the yeast suspension from the swab by stirring and gently shaking.

Then discard the swab. Do not seal the **NBB[®]-B** Broth tube completely to allow CO₂ to diffuse during incubation.

Label the tube.


11.1.3. Incubation of NBB[®]-B Broth tubes

Incubate the **NBB[®]-B** Broth tubes with their screw tops loosely attached under anaerobic conditions (e.g. in an anaerobiosis jar) for 5 days at 27 ± 2 °C in an incubator.

The incubation time will need to be extended for some slow-growing beer spoiling bacteria (e.g.: *Lactobacillus lindneri*).

11.1.4. Evaluation of NBB[®]-B Broth tubes

A change in colour from red to yellow shows the growth of typical beer spoiling bacteria in **NBB[®]-B** Broth. As well as sedimentation and turbidity, gas formation can occur.

Further analysis techniques, such as microscopy and PCR analysis, can be used to identify the bacteria.

11.2. NBB[®]-B Broth: Detection of beer spoiling bacteria from yeast samples in swing top bottles

11.2.1. Addition of yeast into swing top bottles

Use a sterile pipette to place 0.25–0.75 ml of yeast suspension into a transparent swing top bottle (50 ml).

If the yeast cannot be pipetted, you can also transfer it to the bottle with a sterile spatula or by pouring.

11.2.2. Addition of NBB[®]-B Broth in swing top bottles

Pipette approx. 25 ml of **NBB[®]-B** Broth into the swing top bottle that already contains the yeast suspension.

Close the swing top bottle, but do not seal the swing top completely. This prevents explosions due to gas formation.

Label the bottle.

11.2.3. Incubation of NBB[®]-B Broth swing top bottles

Incubate the **NBB[®]-B** Broth swing top bottle for 5 days at $27 \pm 2^\circ\text{C}$ in an incubator. The brief fermentation of the yeast means that enough CO_2 is produced to generate an anaerobic environment. Shortly afterwards, inhibitors reduce the yeast growth.

The incubation time will need to be extended for some slow-growing beer spoiling bacteria (e.g.: *Lactobacillus lindneri*).

11.2.4. Evaluation of NBB[®]-B Broth swing top bottles

A change in colour from red to yellow shows the growth of typical beer spoiling bacteria in **NBB[®]-B** Broth. Turbidity and gas formation can also occur.


Further analysis techniques, such as microscopy and PCR analysis, can be used to identify the bacteria.

11.3. NBB[®]-B Broth: Detection of beer spoiling bacteria from CO₂ and compressed air

11.3.1. Preparation of NBB[®]-B Broth tubes

Fill an appropriate test tube to approx. 70% with **NBB[®]-B** Broth in a sterile workbench using a sterile serological pipette, or use pre-filled **NBB[®]-B** Broth tubes (art. no. 2.04723.646).

11.3.2. CO₂ and compressed air sampling with NBB[®]-B Broth

Sterilise your CO₂ or compressed air outlet (e.g. using a 70% alcohol solution).

Do not allow gas to shoot out, as most bacteria are found in the sampling tap. Allow CO₂ or compressed air to flow into the tube filled with **NBB[®]-B** Broth.

After a defined period (e.g. 15 seconds), shut off the flow of gas.


Close the screw top loosely to allow gas diffusion and label it.

11.3.3. Incubation of NBB®-B Broth tubes

Incubate the **NBB®-B** Broth tubes under anaerobic conditions (e.g. in an anaerobiosis jar) for 5 days at 27 ± 2 °C in an incubator.

The incubation time will need to be extended for some slow-growing beer spoiling bacteria (e.g.: *Lactobacillus lindneri*).

11.3.4. Evaluation of NBB®-B Broth tubes

A change in colour from red to yellow shows the growth of typical beer spoiling bacteria in **NBB®-B** Broth. Turbidity and gas formation can also occur.

Further analysis techniques, such as microscopy and PCR analysis, can be used to identify the bacteria.

12. Protocols NBB[®]-C Concentrate

DMD [®] product	Application	Frequency	Protocol
NBB [®] -C Concentrate	Enrichment of turbid beer samples	Regular (every production cycle)	12.1.
	Enrichment of wort samples (sterile beer addition)	Occasional	12.2.
	Yeast sampling	Occasional	12.3.

12.1. NBB[®]-C Concentrate: Enrichment of beer spoiling bacteria from turbid beer samples

Detection of beer spoiling bacteria with **NBB[®]-C Concentrate** is not based on a change in colour. The samples incubated with **NBB[®]-C Concentrate** are analysed under a microscope (e.g. in a counting chamber). The selectivity of the enrichment can be varied by mixing the beer sample with **NBB[®]-C Concentrate** and adding sterile water, which dilutes the bitter hop compounds. This makes it possible to detect everything from obligate beer spoiling bacteria to indicator microorganisms (see Figure 2 below). For indicator microorganisms air should be left in the bottle head space.

12.1.1. NBB[®]-C Concentrate: Preparation of swing top bottles

Turbid beer samples can be incubated in sterile, transparent swing top bottles for enrichment or for the detection of beer spoiling bacteria.

Warning: In order to produce sterile swing top bottles, fill the bottles with approx. 10 ml of deionised water and autoclave. The water does not have to be removed from the bottle before use.

This small amount of water can be left in the bottle due to the higher hop concentration of unfiltered beer samples.

The swing top on the bottle should be made from stainless steel to prevent corrosion during autoclaving.


Figure 2: Selective detection of bacteria with **NBB®-C** Concentrate. The selectivity of the enrichment can be changed through the addition of sterile water. This allows enrichment of everything from obligate beer spoiling bacteria (e.g. *L. brevis*, *Pectinatus* spp. & *Megasphaera* spp.) to indicator microorganisms by adding air to the bottle headspace (e.g. acidic acid bacteria) (see Back, W., 2005). The volumes stated relate to a 180 ml swing top bottle and may need to be adapted to the actual bottle volume used.

12.1.2. Initiating enrichment with NBB®-C Concentrate in swing top bottles or beer bottles

A) Enrichment in swing top bottles

The beer volume to be investigated is placed into the pre-sterilised swing top bottles in a sterile workbench.

The bottle is then filled to 5% of the total volume with **NBB®-C** Concentrate.

If needed, the selectivity of the enrichment can be modified through the addition of sterile deionised water or sterile tap water, as depicted in Figure 2, p. 32. For example: Adding more than 25% (v/v) water enables the enrichment of potentially beer spoiling bacteria. The bottle should be filled to the top with water. If you want to enrich indicator microorganisms, the headspace has to be filled with air.

Close the sample bottle with the lid. For samples that contain yeast (e.g. yeast containing weiss beer), do not close the bottle tightly in order to allow gas diffusion.

Label the sample.

The lid can be sealed tight after three days of incubation (see 12.1.3, p. 33).

B) Enrichment in beer bottles

Enrichment can also be conducted directly in the beer bottle in order to test beer after filling.

Remove approx. 5% of the bottle volume and fill the bottle almost to the top with approx. 5% **NBB[®]-C** Concentrate.

Close the sample bottle with the lid. For samples that contain yeast (e.g. yeast containing weiss beer), do not close the bottle tightly in order to allow gas diffusion.

Label the sample.

The lid can be sealed tight after three days of incubation (see 12.1.3, p. 33).

12.1.3. Incubation of NBB[®]-C Concentrate

Incubate the bottle with NBB[®]-C Concentrate for 7 – 12 days at 27 ± 2 °C in an incubator or incubation room.

Please note that brief yeast growth can cause CO₂ formation, which puts the bottle under pressure. Yeast growth is quickly suppressed by inhibitors. The CO₂ formation replaces the anaerobic incubation.

The incubation time will need to be extended for some slow-growing beer spoiling bacteria (e.g.: *Lactobacillus lindneri*).

12.1.4. Evaluation of NBB[®]-C Concentrate

NBB[®]-C Concentrate does not contain a colour indicator.

Turbidity can display the growth of typical beer spoiling bacteria in **NBB[®]-C** Concentrate macroscopically. The sample should be constantly checked for the growth of beer spoiling bacteria during incubation. In very turbid samples, a microscope is the only reliable way to determine the growth of beer spoiling bacteria. In the microscopic evaluation (e.g. counting chamber; Neubauer improved), the number of bacteria in 16 microscopic fields (= group square) at a magnification of 600 must be significantly greater than 1. When the chambers have a layer thickness of 0.1 mm, there should be around 4 cells per group square. This corresponds to contamination of approx. 10^6 cells/ml. If there are only isolated cells visible, these are dead or harmless bacteria for the most part, which are almost always found in these kinds of samples.

Ambiguous samples can then be enriched with **NBB[®]-B** Broth to make sure. To do this, mix one drop of **NBB[®]-C** sample sediment with approx. 20 ml **NBB[®]-B** Broth and incubate for one day at 27 ± 2 °C.

Further analysis techniques, such as PCR analysis, can be used to identify the bacteria.

12.2. NBB[®]-C Concentrate: Enrichment of beer spoiling bacteria from wort samples

12.2.1. NBB[®]-C Concentrate: Preparation of swing top bottles

Wort samples with addition of sterile beer can be incubated in sterile, transparent swing top bottles for enrichment or for the detection of beer spoiling bacteria.

Warning:

In order to produce sterile swing top bottles, fill the bottles with approx. 10 ml of deionised water or tap water and autoclave it. The water does not have to be removed from the bottle before use. This small amount of water can be left in the bottle due to the higher hop concentration of wort and unfiltered beer samples. The swing top on the bottle should be made from stainless steel to prevent corrosion during autoclaving.

12.2.2. NBB[®]-C Concentrate: Enrichment of wort samples in swing top bottles

For a 200 ml swing top bottle, which holds around ca. 210 ml when filled to the top, place ca. 20 ml of sterile water in the swing top bottle.

Then fill the swing top bottle with ca. 50 ml wort and add approx. 5% **NBB[®]-C Concentrate** (here ca. 10 ml).

Finally, fill the bottle to the top with pasteurised beer (e.g. samples from the shelf-life cabinet). The beer quantity is around 130 ml.

Loosely close the sample bottle with the lid in order to allow gas diffusion and label it.

12.2.3. Incubation of NBB[®]-C Concentrate

Incubate the bottle with the **NBB[®]-C** Concentrate for 7 – 12 days at 27 ± 2 °C in an incubator or incubation room.

The incubation time will need to be extended for some slow-growing beer spoiling bacteria (e.g.: *Lactobacillus lindneri*).


12.2.4. Evaluation of NBB[®]-C Concentrate

NBB[®]-C Concentrate does not contain a colour indicator. Turbidity can display the growth of typical beer spoiling bacteria in **NBB[®]-C** Concentrate macroscopically. The sample should be constantly checked for the growth of beer spoiling bacteria during incubation. In slightly turbid samples, a microscope is the only reliable way to determine the growth of beer spoiling bacteria. In the microscopic evaluation (e.g. counting chamber; Neubauer improved), the number of bacteria in 16 microscopic fields (= group square) at a magnification of 600 must be significantly greater than 1. When the chambers have a layer thickness of 0.1mm, there should be around 4 cells per group square. This corresponds to contamination of approx. 10^6 cells/ml. If there are only isolated cells visible, these are dead or harmless bacteria for the most part, which are almost always found in these kinds of samples.

Ambiguous samples can then be enriched with **NBB[®]-B** Broth to make sure. To do this, mix one drop of **NBB[®]-C** sample sediment with approx. 20 ml **NBB[®]-B** Broth and incubate for one day at 27 ± 2 °C.

Further analysis techniques, such as PCR analysis, can be used to identify the bacteria.

12.3. NBB[®]-C Concentrate: Enrichment of beer spoiling bacteria from yeast samples (alternative protocol)

12.3.1. NBB[®]-C Concentrate: Preparation of sample enrichment

Autoclave a 200 ml Erlenmeyer flask with a cotton stopper.

In a sterile workbench, add approx. 10 ml of the yeast sample to be tested to the sterile flask.

Then add approx. 90 ml of sterile beer.

Pipette approx. 5 ml **NBB[®]-C** Concentrate into the flask.

Close the flask with the sterile cotton stopper.

Mix the liquids by gently swirling the closed flask.

12.3.2. Incubation of NBB[®]-C Concentrate with yeast samples

Incubate the flask with the yeast sample on a shaker (slow rotation) in an incubator for 5 – 12 days at $27 \pm 2^\circ\text{C}$.

The brief yeast fermentation generates an anaerobic atmosphere in the flask, so anaerobic cultivation is not essential.

12.3.3. Evaluation of NBB[®]-C Concentrate with yeast samples

NBB[®]-C Concentrate does not contain a colour indicator. Turbidity can display the growth of typical beer spoiling bacteria in **NBB[®]-C** Concentrate macroscopically. The sample should be constantly checked for the growth of beer spoiling bacteria during incubation. In very turbid samples, a microscope is the only reliable way to determine the growth of beer spoiling bacteria. In the microscopic evaluation (e.g. counting chamber; Neubauer improved), the number of bacteria in 16 microscopic fields (= group square) at a magnification of 600 must be significantly greater than 1. When the chambers have a layer thickness of 0.1 mm, there should be around 4 cells per group square. This corresponds to contamination of approx. 10^6 cells/ml. If there are only isolated cells visible, these are dead or harmless bacteria for the most part, which are almost always found in these kinds of samples.

Ambiguous samples can then be enriched with **NBB[®]-B** Broth to make sure. To do this, mix one drop of **NBB[®]-C** sample sediment with approx. 20 ml **NBB[®]-B** Broth and incubate for one day at 27 ± 2 °C.

Further analysis techniques, such as PCR analysis, can be used to identify the bacteria.

13. FAQs on application and evaluating the results

13.1. NBB[®]-A Agar

	Cause	Solution
Strictly anaerobic bacteria like <i>Pectinatus</i> and <i>Megasphaera</i> do not grow on the membrane filter.	Especially strict anaerobic beer spoiling bacteria are often exposed to oxygen during the membrane filtration process, which inhibits their growth.	After membrane filtration, either slightly flush the membrane filter with CO ₂ or incubate the membrane filter in NBB[®]-B in an anaerobiosis jar.
After storage in a water bath, NBB[®]-A Agar no longer solidifies completely or at all.	Storing liquid agar at approx. 45°C or at 95°C in a water bath causes chemical damage to the agar structure. This prevents the agar from solidifying.	Never store NBB[®]-A Agar in a water bath at 45°C for more than 4 h. Avoid extended liquefaction times at 95°C or multiple liquefactions.
In the pour-plate method, the colour changes as soon as the samples are added.	You have added too much sample or the sample is excessively acidic to the NBB[®]-A Agar.	Adhere to the maximum sample volumes specified in the protocol (see Protocol 10.6., p. 23). If your sample is very acidic, reduce the sample volume
Bacteria growth with no change in colour.	Microorganisms that do not cause a change in colour on NBB[®]-A Agar should be analysed with a microscope. Indicator microorganisms (e.g. <i>Enterobacter</i> spp.) may grow in case of insufficient anaerobiosis.	For ambiguous samples, you can conduct a bacteria examination under a microscope or a differentiation with PCR methods. You should also check your incubation conditions.

13.2. NBB[®]-B Broth

	Cause	Solution
Turbidity with no change in colour.	If the nutrient medium becomes turbid after 5 days of incubation and there is no colour change, this implies a contamination that may have occurred during sampling.	Please repeat the sampling and comply with the outlined work steps. Please pay attention that microbiological work is carried out in a hygienic manner.
Differentiating between or identifying beer spoiling bacteria.	It is not possible to differentiate or identify beer spoiling bacteria with NBB[®]-B Broth.	Further analysis techniques, such as PCR-based methods, can be applied to differentiate bacteria.
Observing yeast growth.	Shortly after the yeast is added, particularly in large quantities, it may grow briefly in the NBB[®]-B Broth. However, the yeast's growth will be inhibited after a short period. This leads to CO ₂ formation in the short-term and thus provides an ideal anaerobic atmosphere. The yeast quantity may also cause turbidity.	After 5 days, only beer spoiling bacteria are detected in the analysis. Please observe the exact analysis times. If the solution becomes turbid, please also pay attention to a change in colour as this often indicates the presence of beer spoiling bacteria (adding too much yeast prevents the colour change). Unclear results can also be examined under a microscope where beer spoiling bacteria can be recognized in addition to many yeast cell.

13.3. NBB[®]-C Concentrate

	Cause	Solution
The colour change is not visible.	NBB [®] -C Concentrate does not contain a colour indicator.	The beer spoiling bacteria can be determined with a microscope. Only significant growth indicates a result.
Enriching clear beer samples with NBB [®] -C Concentrate.	–	Clear beer samples can also be enriched with NBB [®] -C Concentrate for qualitative observation. Analysis is done by assessing the sediment and turbidity and using a microscope.

14. Appendix

14.1 Brief overview: Detection of beer spoiling bacteria using NBB®-A Agar

	Membrane filtration (clear beer samples/ water)	Plating & spreading method	Air samples	CO ₂ & compressed air	Pour-plate method	
Sample volume	Individual volume (default: 300–500 ml)	Max. 150 µl for plating/ loop volume for spreading	Individual air volume with air sampler or individual time interval with sedimentation process	Membrane filtration of ca. 50 ml of sterile water with CO ₂ /compressed air (15 seconds)	Petri dish diameter	Max. sample volume
					9 cm	5 ml
					14 cm	10 ml
Incubation conditions	Temp.: 27±2°C Time: 5–7 Days	Temp.: 27±2°C Time: 5 Days	Temp.: 27±2°C Time: 5–7 Days			
Analysis	Positive: Colour change, red to yellow, colony formation Negative: No colour change, no colony formation					

14.2 Brief overview: Detection of beer spoiling bacteria using NBB®-B Broth

	Yeast samples (tube)	Yeast samples (swing top bottle 50 ml)	CO ₂ & compressed air
Sample volume	0.1–0.3 ml of yeast or swab sample	0.25–0.75 ml yeast sample	Approx. 10 ml NBB®-B Broth with individual volume of CO ₂ /compressed air (15 seconds)
Incubation conditions	Temp.: 27 ± 2 °C Time: 5 days		
Analysis	Positive: Colour change, red to yellow, clear turbidity (gas formation) Negative: No colour change, no turbidity		

14.3 Brief overview: Detection of beer spoiling bacteria using NBB®-C Concentrate

	Turbid beer samples	Wort samples				Yeast samples (alternative protocol)
Sample volume	Variable (example: at a total volume of 180 ml, approx. 100–170 ml of beer sample, depending on selectivity)	Variable Example for 200 ml swing top bottle. The volume needs to be adapted to the bottle volume.				10 ml yeast suspension (see p. 37)
		Wort	Sterile water	NBB®-C Concentrate	Pasteurised beer	
		ca. 50 ml	ca. 20 ml	ca. 10 ml	ca. 130 ml (filled to the top)	
Incubation conditions	Temp.: 27 ± 2 °C Time: 7–12 days					Temp.: 27 ± 2 °C Time: 5–12 days
Analysis	Positive: Microscopic detection: In 16 microscopic fields at a magnification of 600, significantly more than 1 beer spoiling microorganism in each field of view. Negative: Less than 1 bacteria cell per field of view can be microscopically detected.					

14.4. Glossary (in alphabetical order)

Differentiation	The process of identifying bacteria on the basis of their biochemical or genetic properties down to their genus and species names, e.g. <i>Lactobacillus brevis</i> .
Inert	Is an object that behaves in a stable and non-chemically-reactive way under certain conditions.
Incubation	The process of incubating culture media under defined conditions and parameters (e.g. in terms of time and temperature).
Peleus ball	A rubber ball with valves used for sucking in and discharging liquids from glass or plastic pipettes.
Qualitative detection	Detection based on a yes/no statement, i.e. the sample does or does not contain microorganisms.
Quantitative detection	Detection based on a value-based analysis i.e. the number of microorganisms in a sample can be determined.
Serological pipette	Pipette usually made out of plastic or glass used for transferring liquids. The scale featured on the pipette displays the volume. Plastic pipettes are usually available in individual, sterile packaging.
Sterile	Sterile means free from bacteria and living microorganisms.

15. Supply sources

Döhler provides you with various solutions for your quality control:

Article no.	Product	Packaging unit
2.04709.782	NBB®-A Agar	1x cardboard box (9x 250 ml in bottle)
2.04710.782	NBB®-B Broth	1x cardboard box (9x 250 ml in in bottle)
2.04723.646	NBB®-B Broth (tubes)	1x polystyrene box (20x 10 ml in tubes)
2.04711.782	NBB®-C Concentrate	1x cardboard box (9x 250 ml in bottle)

Sterile wooden swabs are also available as accessories:

Article no.	Product	Packaging unit
2.04725.244	sterile (wooden) swabs (individually packaged)	100 pc./package

All labware can be obtained from a laboratory supplier.

| 16. Information

16.1. Trademarks

NBB® and Doehler Microsafety Design® are both trademarks of Doehler GmbH that are registered and protected worldwide.

15.2. Distribution

You can find more information about our Doehler Microsafety Design® product portfolio on the website www.doehler-dmd.com.

Email contact: dmd@doehler.com

15.3. Recommended further reading

Back, W.: Colour Atlas and Handbook of Beverage Biology, Fachverlag Hans Carl, Nuremberg, 2005.

Back, W. (pub.): [Selected chapters of brewing technology, second updated edition, p. 319 seq., Fachverlag Hans Carl, Nuremberg, 2008.]

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DÖHLER GmbH