
**Revisions**

Rev from	Rev to	ECO #
0905	2010/07	5384-10

**Notes:**

- BD Cat. Number 231044, 231046
- Blank (Sheet) Size : Length: N/A      Width: N/A  
 Number of Pages: N/A      Number of Sheets: N/A  
 Page Size: Length N/A      Width N/A      Final Folded Size: N/A
- Style (see illustrations below): # N/A

- See Specification Control Number VS-88-4070-1 for Material Information
- Ink Colors: Printed two sides  Yes     No  
 No. of Colors: N/A      PMS# N/A
- Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level.

Vendor Spec Controlled by BD Caribe, LTD.

Label Design	Date	COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION	 Becton, Dickinson and Company 7 Loveton Circle Sparks, MD 21152 USA	
Proofer	Date			
Checked By	Date			
Part Number:	Category and Description Package Insert, BBL Taxo Discs for Identification of Neisseria and Pseudomonas		Sheet: 1 of 15 Scale: N/A	<b>A</b>
8840701				

# **BD BBL™ Taxo™ Discs for Identification of *Neisseria* and *Pseudomonas***

English: pages 1 – 3 Italiano: pagine 8 – 10  
Français: pages 3 – 5 Español: páginas 10 – 12  
Deutsch: Seiten 5 – 7

 8840701  
2010/07

Pokyny vám poskytne miestni zástupce spoločnosti BD. / Kontakt den lokale BD repræsentant for at få instruktioner. / Kasutusjuhiste suhtes kontakteeruge oma kohaliku BD esindajaga. / Επικοινωνήστε με τον τοπικό αντιπρόσωπο της BD για οδηγίες. / A használati utasítást kérje a BD helyi képviselőjétől. / Naudojimo instrukcijų teiraukitės vietos BD įgaliotojo atstovo. / Kontakt din lokale BD-representant for mer informasjon. / Aby uzyskać instrukcje użytkowania, skontaktuj się z lokalnym przedstawicielstwem BD. / Contacte o seu representante local da BD para obter instruções. / Instrukcie získate u miestneho zástupcu spoločnosti BD. / Kontakta lokal Becton Dickinson-representant för anvisningar. / Свържете се с местния представител на BD за инструкции. / Contactați reprezentantul dumneavoastră local BD pentru instrucțiuni. / Talimatlar için yerel BD temsilcilerinize danışın. / Obratite se svom lokalnom predstavniku kompanije BD za uputstva. / Для получения инструкций свяжитесь с местным представителем компании BD. / Өзіндік жергілікті БД өкіліне жүгініп нұсқау алыңыз. / Kontaktiraj lokalnog predstavnika BD za upute.

## INTENDED USE

Taxo™ N discs are used in qualitative procedures to distinguish those organisms which produce oxidase. By this reaction, Taxo N discs may be used for the presumptive differentiation of the genus *Neisseria* from other gram-negative cocci and members of the genus *Pseudomonas* from enteric bacilli.

## SUMMARY AND EXPLANATION

Gordon and McLeod described the oxidase test for detecting colonies of *Neisseria*.<sup>1</sup> Kovacs described the oxidase test for *Pseudomonas*.<sup>2</sup>

*Neisseria*, *Pseudomonas* and certain other microorganisms possess an oxidizing enzyme which acts in the presence of air on certain aromatic amines to produce colored compounds.<sup>3,4</sup> Color changes are produced by the indicator in contact with a colony of the organism.

Because the reagent used for the test is very unstable in liquid form, the original test procedures required that the reagent be made fresh, the dye being active only a few hours. The reagent remains quite stable, however, when impregnated into paper discs, dried, and stored in the refrigerator or freezer.

## PRINCIPLES OF THE PROCEDURE<sup>3</sup>

The oxidase test is based on the presence of an intracellular cytochrome oxidase system which activates the oxidation of reduced cytochrome by molecular oxygen; the cytochrome serves as an electron acceptor in the terminal stage of the electron transfer system.

*Neisseria* and nearly all *Pseudomonas* spp. produce an oxidase enzyme which, in the presence of atmospheric oxygen, cytochrome C and the oxidase reagent, *p*-aminodimethylaniline monohydrochloride, oxidizes the reagent to form a colored compound, indophenol.

## REAGENTS

Taxo N discs are 6 mm discs made from high-quality absorbent paper impregnated with 6% *p*-aminodimethylaniline monohydrochloride.

### Warnings and Precautions:

For *in vitro* Diagnostic Use.

Observe aseptic techniques and established precautions against microbiological hazards throughout all procedures. After use, contaminated materials must be sterilized by autoclaving.

**Storage Instructions:** On receipt, store at -20 to +8°C. After use, store vial to protect product integrity at 2 – 8°C. Apply control tests at time of use.

Use the oldest discs first and discard expired discs. Discard discs left out overnight in the laboratory, or else test the discs for performance.

Expiration date applies to product in intact container stored as directed. Do not open until ready to use.

**Product Deterioration:** Discs should be light to medium gray (trace purple acceptable) and uniform. Do not use discs if they show evidence of discoloration or other signs of deterioration. Discs that have absorbed moisture may give false-positive results.

## SPECIMENS

These discs may be used with pure cultures or with specimens containing mixed flora, but the organisms to be tested must be present as separate colonies.

## PROCEDURE

**Material Provided:** Taxo N Discs.

**Materials Required But Not Provided:** Ancillary culture media, quality control organisms and laboratory equipment as required for this procedure.

### Test Procedure:

#### A. Rapid Screening Test

Colonies and broth cultures of *Neisseria* and *Pseudomonas* may both be rapidly and conveniently tested in the following manner. Place a number of **Taxo N** discs in a Petri dish or other convenient container and moisten with distilled or deionized water. Using a platinum wire or sterile wooden applicator stick, apply suspicious colonies to the moistened discs.

Alternatively, place a number of drops of distilled or deionized water in a sterile plate; plates with grids are especially convenient. To each drop add a generous amount of culture from a tube or plate; mix to obtain a heavy suspension. Add a **Taxo N** disc to each drop of suspension. The color change occurs promptly, and can be hastened by incubation at 35°C for 3 – 5 min.

#### B. Plate Tests

Dip **Taxo N** discs in sterile water to thoroughly wet them and place on the surface of cultures which have already grown out during incubation. They may be applied to cultures grown on Chocolate Agar or media for the selective isolation of pathogenic *Neisseria*, **Trypticase™** Soy Agar (*Pseudomonas*), or other suitable media. After application of discs near suspected colorless colonies, return the plates to the incubator for 20 – 30 min; the reaction also takes place, though more slowly, at room and refrigerator temperatures.

### User Quality Control:

1. Examine discs for signs of deterioration as described under "Product Deterioration."
2. At the time of use, check performance with pure cultures of stable control organisms producing known, desired reactions. The following cultures are recommended:

Test Strain	Expected Results
<i>Neisseria gonorrhoeae</i> ATCC™ 43069	Positive; dark purple to black
<i>Pseudomonas aeruginosa</i> ATCC 27853	Positive; dark purple to black
<i>Escherichia coli</i> ATCC 25922	Negative; no color change

Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent NCCLS guidance and CLIA regulations for appropriate Quality Control practices.

## RESULTS

#### A. Rapid Screening Test

Those cultures which contain oxidase-producing organisms will show a reaction within 5 min with the **Taxo N** discs turning dark purple to black. By this means, both *Neisseria* and *Pseudomonas* may be readily differentiated from other bacteria in the culture.

#### B. Plate Tests

Colonies of *Neisseria* and *Pseudomonas* that are within a radius of approximately 10 mm of the disc will turn dark purple and then black, due to the oxidase produced by these organisms. Colonies of oxidase-negative bacteria, such as the gram-positive cocci and coliforms, do not change color, although the medium may be somewhat discolored by the reagent. Continued holding of plates permits further diffusion from the discs and blackening of additional oxidase-positive colonies.

### LIMITATIONS OF THE PROCEDURE<sup>3-6</sup>

**Taxo N** discs are intended only as an aid in the identification of *Neisseria* and *Pseudomonas*.

Not all *Pseudomonas* spp. are oxidase positive. Other gram-negative organisms that may be oxidase positive include species of *Aeromonas*, *Alcaligenes*, *Campylobacter*, *Haemophilus*, *Kingella*, *Moraxella*, *Pasteurella*, *Plesiomonas* and *Vibrio*.

Isolates from selective-differential media for gram-negative bacilli (e.g., MacConkey, EMB, XLD agars) may give false-negative reactions.

A false-positive reaction may result if an iron-containing wire (e.g., nichrome or stainless steel) is used to transfer growth. Use a platinum wire or sterile wooden applicator stick.

A false-negative result may occur with *Neisseria* if a mixed culture contains both *Neisseria* and *Pseudomonas* spp.

Biochemical and serological tests should be performed for complete identification. Appropriate references should be consulted for further information.<sup>4-7</sup>

### PERFORMANCE CHARACTERISTICS

Prior to release, all lots of **BBL™ Taxo™ N** Discs are tested to verify specific product characteristics. Samples are tested by two methods. In the first method, discs are placed in sterile disposable Petri dishes. The discs are then moistened with distilled or deionized water. A loopful of cultures of *Neisseria gonorrhoeae* (ATCC™ 43069 and ATCC 43070), *N. meningitidis* (ATCC 13090), *N. subflava* (ATCC 14799), *Pseudomonas aeruginosa* (ATCC 27853) and *Escherichia coli* (ATCC 25922) are placed on the moistened discs. The discs inoculated with the *Neisseria* and *Pseudomonas*

cultures turn a dark purple to black color indicating a positive reaction. The disc inoculated with *Escherichia coli* is negative and shows no color change.

In the second method, Chocolate Agar plates are inoculated with the above strains of *Neisseria* cultures, and **Trypticase™** Soy Agar with 5% Sheep Blood plates are inoculated with the strains of *Pseudomonas* and *Escherichia* cultures noted above and incubated. After good growth is obtained, moistened discs are placed on the surface of the plates in the area of growth. The plates are then incubated at 35 ±2°C for 20 to 30 min. After incubation, colonies of *Neisseria* and *Pseudomonas* nearest to the disc are positive, exhibiting a dark purple to black color. Colonies of *Escherichia coli* closest to the disc are negative and show no color change.

#### AVAILABILITY

Cat. No.	Description
231044	Taxo™ N Discs, Single vial of 50 discs.
231045	Taxo™ N Discs, Pkg. of 6 vials.

#### REFERENCES

1. Gordon, J., and J.W. McLeod. 1928. The practical application of the direct oxidase reaction in bacteriology. *J. Pathol. Bacteriol.* 31:185.
2. Kovacs, N. 1956. Identification of *Pseudomonas pyocyanea* by the oxidase reaction. *Nature* 178:703.
3. MacFaddin, J.F. 2000. Biochemical tests for identification of medical bacteria, 3rd ed. Lippincott Williams & Wilkins, Baltimore.
4. Shigei, J. 1992. Oxidase test, p.1.19.16. In H.D. Isenberg (ed.), *Clinical microbiology procedures handbook*, vol. 1. American Society for Microbiology, Washington, D.C.
5. Baron, E.J., and S.M. Finegold. 1990. *Bailey & Scott's diagnostic microbiology*, 8th ed. The C.V. Mosby Company, St. Louis.
6. Balows, A., W.J. Hausler, Jr., K.L. Herrmann, H.D. Isenberg, and H.J. Shadomy (ed.). 1991. *Manual of clinical microbiology*, 5th ed. American Society for Microbiology, Washington, D.C.
7. Krieg, N.R., and J.G. Holt (ed.). 1984. *Bergey's manual of systematic bacteriology*, vol. 1. Williams & Wilkins, Baltimore.

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ATCC is a trademark of the American Type Culture Collection.
BD, BD Logo, BBL, Taxo and Trypticase are trademarks of Becton, Dickinson and Company. © 2010 BD.
Benex Limited