-  **(EN) Fraser Broth Base**
-  **(FR) Base Bouillon Fraser**
-  **(DE) Fraser Basis Bouillon**
-  **(IT) Brodo di base Fraser**
-  **(ES) Base de Caldo Fraser**
-  **(NL) Fraser bouillon (basis)**
-  **(SV) Fraser buljong**
-  **(DA) Fraser Bouillon Base**
-  **(NO) Fraser buljong**
-  **(FI) Fraser Elatusainejauhe**
-  **(PT) Fraser Caldo Base**
-  **(EL) Βάση Ζωμού Fraser**
-  **(PL) Bulion do frezowania**
-  **(HU) Fraser Tápleves Alap**
-  **(CS) Fraserův bujón – báze**
-  **(RO) Bulion bază Fraser**
-  **(RU) Бульон Фрейзера**
-  **(TR) Toz Fraser Broth**
-  **(JA) フレーザー基礎培地**
-  **(ZH) Fraser 肉汤基础**
-  **(TH) เฟรเซอร์บรอตเบส**
-  **(KO) Fraser Broth Base**

# Product Instructions

## Fraser Broth Base

### PRODUCT DESCRIPTION

Medium used for the selective enrichment of *Listeria* species including *Listeria monocytogenes*.

### APPLICATIONS

For Laboratory Use Only. Material Safety Data Sheet available.

### COMPOSITION

(typical formula g/L)

Sodium Chloride	20.0 g
Sodium Phosphate, Dibasic CAS # 7558-79-4	9.6 g
Casein Peptone	5.0 g
Meat Peptone	5.0 g
Beef Extract	5.0 g
Yeast Extract	5.0 g
Lithium Chloride	3.0 g
Potassium Phosphate, Monobasic CAS # 7778-77-0	1.35 g
Esculin	1.0 g
Acriflavin	0.025 g
Nalidixic Acid	0.02 g
Final pH 7.2 ± 0.2 at 25°C	

### PREPARATION OF MEDIUM

Suspend 55 g of the powder in 1.0 L of purified water. Mix thoroughly. Autoclave at 121°C for 15 minutes. Cool to room temperature.

**⚠ WARNING!** Irritating to eyes, respiratory system and skin. Can cause central nervous system depression. Contains a chemical or chemicals which can cause birth defects or other reproductive harm.

### ADDITIVE

3M™ Fraser Broth Supplement (**REF:** BP0220010) (5% Ferric Ammonium Citrate) is required for differential enrichment of *Listeria*. If the addition of 3M™ Fraser Broth Supplement is necessary, cool the 3M Demi Fraser Broth Base to room temperature and aseptically add 10 mL of the 3M Fraser Broth Supplement per liter of 3M Fraser Broth Base. Mix thoroughly.

**⚠ WARNING!** Do not add 3M Fraser Broth Supplement (5% Ferric Ammonium Citrate) when used in conjunction with the 3M™ Molecular Detection Assays.

### STORAGE

2-30°C. Hygroscopic. Keep container tightly closed.

### DISPOSAL

Follow established laboratory practices.

**Disclaimer:** 3M disclaims all express and implied warranties including warranties of merchantability or fitness for a particular use. If product is defective, the exclusive remedy is, at 3M's option, replacement, or refund. Except where prohibited by law, 3M will not be liable for any further loss or damage arising from use of this product.

Users are responsible for familiarizing themselves with product instructions and information. Visit our website at [www.3M.com/foodsafety](http://www.3M.com/foodsafety), or contact your local 3M representative or distributor for more information.

## EXPLANATION OF SYMBOLS

 Caution or Warning, consult product instructions.

 Consult product instructions.

**LOT** The lot in a box represents the lot number.

 The hourglass is followed by a year, month and day which represent the expiration date.

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© 2013, 3M. All rights reserved.
3M is a trademark of 3M. Used under license in Canada.
3M Health Care Limited