

MILK CONTROLS

Safety statement

Introduction

The Milk Controls can be used for checking and monitoring a number of performance parameters of antibiotic residue test kits that are in use in the dairy industry for checking cow, farm, truck and process milk for the presence of antibiotic residues.

Composition/information on Milk Controls

The Negative Control contains antibiotic free milk powder derived from commingled cow milk.

The positive control samples contain antibiotic free milk powder derived from commingled cow milk and additional an antibiotic.

The 4 ppb Penicillin positive Control contains 4 +/- 0.5 ppb Penicillin G.

The 5 ppb Penicillin Positive Control contains 5 +/- 0.5 ppb Penicillin G.

The Sulfa Positive Control contains approximately 1000 ppb Sulfadiazine.

The MILK CONTROLS are designed for laboratory use.

For more information please contact:

or our local representative:

DSM Food Specialties B.V.

P.O. Box 1, 2600 MA Delft

The Netherlands

Tel.: +31 15 279 4001

info.dairy-ingredients@dsm.com

www.dsm-foodspecialties.com

Although diligent care has been used to ensure that the information provided herein is accurate, nothing contained herein can be construed to imply any representation or warranty as to the accuracy, currency or completeness of this information. The content of this document is subject to change without further notice. This document is non-controlled and will not be automatically replaced when changed. Please contact us for the latest version of this document or further information.