

Plain Dry Swabs In Tubes CLASSIQSwabs

Sterile plain swabs in labelled tubes.

Code	Applicator	Swab Tip	Color Code
150C	Wood	Cotton	Red
152C	Wood	Cotton/Charcoal coated	Black
154C	Wood	Rayon	White
155C	Plastic	Rayon	White
167C*	Double Plastic	Rayon	White
157C	Plastic	Cotton/Charcoal coated	Black
159C	Plastic	Dacron	Green
160C	Aluminum	Rayon	Orange
162C	Aluminum	Dacron	Yellow
164C	Paper	Cotton	Green
168C	Twisted Wire	Rayon	Blue
184C	Plastic minitip	Rayon	
185C	Plastic minitip	Polyester	
451C	Wood	Vegetal protein coated	
456C	Plastic	Vegetal protein coated	

All swabs are packaged 10x100 swabs/case



* These swabs are packaged 10x50 swabs/case

Swabs In Peel Pouches CLASSIQSwabs

Sterile Swabs individually wrapped in peel pouches

Catalog No.	Applicator	Swab Tip
165KS01	Wood	Cotton
167KS01	Plastic	Rayon
164KS01	Plastic	Dacron
170KS01	Aluminum	Rayon
175KS01	Aluminum	Dacron
184CS01*	Plastic minitip	Rayon
185CS01*	Plastic minitip	Polyester
A905BCS01*	Large plastic	Rayon

* These swabs are packaged 10x50 swabs/case

All swabs are packaged 100 swabs/box, 10x100 swabs/case



Other Copan Products

Code	pack size
706C	Mailing container 26x26 suitable for swab
706CCA	Mailing container 26x26 suitable for swab with absorbing paper