



Scale 1:1

Order number 51343154
former order number -

Specifications
pH range 1 - 11
zero point (pH) 7.00 ± 0.25

temperature range 0 - 80 °C
temperature probe NTC 30 kΩ

shaft material glass
membrane glass LoT
membrane resistance < 250 MΩ

reference system ARGENTHAL™
reference diaphragm open junction
reference electrolyte XEROLYT®PLUS polymer

bridge electrolyte -

electrode head MultiPIN™
cable length -
connector -

cell constant (cm⁻¹) -
cell type -

indicating metal -

Storage FRISCOLYT-B®
(51340053)

Special features
high performance puncture electrode:
- XEROLYT®PLUS polymer electrolyte
- open reference junction
- ARGENTHAL™ reference system
integrated temperature probe NTC 30 kΩ

Choice of cables
All Mettler-Toledo meters have BNC connectors.
Recommended cable: BNC/RCA (Cinch) 1.2 m
(52300009)

Our cable selector guide offers a large choice of MultiPin™ cables for other connectors.