interscience

30, CH. BOIS ARPENTS - 78860 ST NOM LA BRETECHE - FRANCE - T. +33 (0)1 34 62 62 61 - F. +33 (0)1 34 62 43 03 N° Identification TVA: FR52 950 356 220 - E-mail: info@interscience.fr - Website: www.interscience.fr



Scan[®] 4000 Automatic HD colony counter

Ref. 438 000

Scan[®] 4000, automatic HD colour colony counter is adapted to all Petri dish sizes and all media. Its indirect lighting system ensures great user comfort, high accuracy and excellent reproducibility.

In coompliance with CFR 21 part 11, EUCAST, CA-SFM, ISO 7218, ISO 4887-2 standards.

TECHNICAL SPECIFICATIONS

- · All stainless steel body
- Color HD CCD camera
- HD japanese lens
- Digital zoom x64 (enabling to see the letters on the disks)
- Resolution: 1.2 mega-pixels
- Counting: 1000 colonies per second
- Minimum size of colony: 0.05 mm
- Lighting esquipment: long-life white LEDs / Dark Field
- · Lighting system: automatic with 6 combinaisons, top and/or bottom light, white or black background
- · Counting: automatic with manual control
- Results: Scan® file (recountable), PDF report, jpeg, png, bmp, Excel™
- LIMS connection
- USB connection
- Counting on Petri dishes 55 150 mm
- Counting on Spiral® plating
- · Color detection: 7 colors on the same dish
- Counting on chromogenic medium
- Counting on PetriFilm[™]
- Counting on RIDA™ Count
- Counting on Compact Dry™
- · Counting on filtration membranes
- · Automatic separation of confluent colonies
- · Creation of polygonal exclusion areas
- Manual control of adding or removing colonies
- · Inhibition zone reading on 90-100 mm round dishes and 120 mm square dishes
- Inhibition zone measurement: automatic with manual control
- Automatic detection of paper disks (all brands simultaneaously), wells, peni-cylinders (steel, plastic), customizable, add/susbtract manually
- · Maximum number of antibiotics: up to 16 antibiotics on a square 120 mm dish
- Reading time: mesurement of 16 inhibition zones between 2 and 4 seconds maximum
- · Interpretation system: CA-SFM Human Health / EUCAST / CA-SFM Veterinary / customizable list
- Security of data: saving of all data changes in compliance with CFR 21 recommendations / Connection to LIMS / SIL systems
- Results / Traceability: image / sample cumber & antibiotic name / bacteria name / read diameter / result with category SIR by colour code (Sensitive : Intermediate : Resistant) / critical diameter mini and maxi / comments / date / time
- · Software available in english, french, spanish, japanese, chinese, russian
- Power: 100-240V~ / 50-60Hz
- Dimensions: (L x I x H) : 47 x 47 x 64 cm, weight: 24,4 kg

- Guarantee (after registration) : 3 years
- Software update (after registration) : 3 years
- · In compliance with CFR 21 part 11, EUCAST, CA-SFM, ISO 7218, ISO 4887-2 standards
- · Designed and made in France
- Manufactured under: ISO 9001 V2008, CE, WEEE 2002/96/EC, RoHS 2002/95/EC

MINIMUM COMPUTER REQUIREMENTS

- · Processor: (quad core recommended) Intel i5 or i7 or AMD FX 6000 or superior
- System: Windows 7 SP1 or 8 or ulterior
- Frequency: 2.4 Ghz or superior
- RAM: 3 Go
- Ports: CD ROM reader, USB port
- Screen resolution: 1280 x 1024 pixels or more
- · Graphic card: AMD or NVIDIA (chipsets not recommended)

DELIVERED WITH:

- Scan® software CD-ROM
- Connexion cables
- 3 control dishes
- User's manual, IQ/OQ/PQ guidelines
- · Quick user guide

OPTIONS/ACCESSORIES

- Ref. 436 001 DataLink™ traceability system
- Ref. 522 000 Bar-code reader