

CORNING	Life Sciences	Document Number: LSR00018	Rev.: 5
Description (Class): Reference, Customer Technical Data Sheet – Serological Pipets, 2 mL 0.01 mL Increments			Page: 1 of 2

APPLICATIONS:

For tissue culture usage or any non-toxic solution that requires aliquoted or aspirating/dispensing in an aseptic environment.

FEATURES:

- Consistent performance
 - Reliability, accuracy is +/-3% at full volume
- Non-pyrogenic
 - Low endotoxins
- Color coding
 - Ease of identification on packaging and blow-out rings
- Reverse graduations
 - Easy reading of the amount of liquid removed or the amount remaining
- Negative graduations
 - Can measure past the stated capacity
- Crisp, legible, rotogravure printing
 - Easily read sample volumes
- Crystal grade polystyrene
 - Can clearly see fluid sample
- Individually wrapped packaging:
 - Two ways to open: Peel or Pop-through

REGULATORY COMPLIANCE:

Falcon® products are manufactured under the current ISO 9001 and ISO 13485 Standards, and the current FDA Quality System Regulation 21 CFR Section 820.

STERILITY:

Product is gamma irradiated

RESTRICTIONS:

- Refer to "Thermoplastics Properties Chart," available at:
http://catalog2.corning.com/Lifesciences/media/pdf/an_DL_226_Falcon_ThermoPlastics_Broch.pdf

CORNING	Life Sciences	Document Number: LSR00018	Rev.: 5
Description (Class): Reference, Customer Technical Data Sheet – Serological Pipets, 2 mL 0.01 mL Increments			Page: 2 of 2

Falcon® 2 mL, 0.01 mL Increments, Serological Pipet

Material: Pipet - Polystyrene; Plug - Polyester fiber

Related Catalog Numbers:	Cat. No.	Packout Description
Sterile	356507	Individually wrapped, 100/box, 1000/case
Sterile	357507	Individually wrapped, 100/box, 1000/case
Sterile	357508	25/bag, 1000/case

Dimensions: = XXX mm (XXX inches)

TIP:

Length: = N/A
 Interior Diameter: = 1.63 (.064)
 Outside Diameter: = 3.73 (.147)

MOUTHPIECE:

Length: = N/A
 Interior Diameter: = 4.17 (.164)
 Outside Diameter: = 6.27 (.247)
 Plug = N/A

EXTRUSION:

Interior Diameter: = N/A
 Outside Diameter: = N/A
 Length = N/A

PIPET SIZE: 2 mL

INCREMENTS: 0.01


GRADUATIONS: Reverse

OVERFILL CAPACITY: 0.08 mL

COLOR CODE: Green

GENERAL DIMENSIONS:

A = 282.58 (11.125)
 B = 3.73 (.147) Diameter
 C = 6.27 (.247) Diameter
 D = 6.27 (.247) Diameter



Warranty/Disclaimer: Unless otherwise specified, all products are for research use only. Not intended for use in diagnostic or therapeutic procedures.