

JumboMix® 3500 P CC

3500 mL lab blender

Ref 031 230



Standards:



Manufactured under:



Product made for INTERSCIENCE by Interlab, an ISO 9001 certified company

JumboMix® 3500 P CC enables the quick homogenizing of large samples. Avoiding all risks of cross-contamination, it is an easy-to-use and powerful lab homogenizer. To improve blending efficiency and lower preparation time, it is equipped with patented adjustable paddles. Adapted to all kind of applications and with a guarantee of optimal bacterial extraction, it is a great tool for lab analyses.

TECHNICAL SPECIFICATIONS

- All stainless steel 304L
- Plain door
- Multifunction digital screen
- Optimal blending volume: 200 - 3750 mL
- Variable speed (1.5 - 3 - 4.5 - 6 strokes/s)
- Variable blending time (1 s - 1 h or ∞ / count-down)
- **Adjustable blending power (patented)**
- **Side-by-side paddle stop**
- **Click & Clean® removable paddles (patented)**
- **Q-Tight® closing: tight closing for no more leaks from the top of the bag**
- Fully opening door (270°)
- Also available in **CoolMix®** version

- Voltage and frequency: 100-127V~ / 220-240V~ 50-60Hz
- Fuses: T 4A L 240V / T 2A L 240V 5x20 mm
- Max power: 400 W

- Dimensions (w x d x h): 52 x 47 x 47 cm, weight: 53.2 kg
- Box (w x d x h): 62 x 58 x 58 cm, weight: 61 kg

- 3-year warranty
- **Shock absorbers lifetime warranty**

- In compliance with: ISO 7218, ISO 6887-1 and FDA BAM (Bacteriological Analytical Manual)
- Manufactured under CE, RoHS, WEEE - Products made for INTERSCIENCE by Interlab, an ISO 9001 certified company
- Designed and made in France

DELIVERED WITH

- 1 power cable
- 1 user's manual
- 1 quick user guide
- 1 free pack of blender bags

OPTIONS / ACCESSORIES

- **Ref 036 000 DripTray:** Stainless steel security drip tray for **JumboMix®**
- **Ref 037 000 JumboStand:** Ergonomic table for **JumboMix® 3500**

1 / 1

JumboMix® 3500 P CC_FT_0518 Pictures and information are not contractually binding. JumboMix® is trademark. INTERSCIENCE reserves the right to change or improve the specifications of its products without notice. Please visit www.interscience.com for current updates and additional information. RCS 950 356 220 Versailles. INTERSCIENCE SARL F78860

www.interscience.com

interscience

30, ch. Bois Arpents - 78860 St Nom - FRANCE
T: +33 (0)1 34 62 62 61
Email: info@interscience.com

interscience USA & CANADA

32 Cummings Park - Woburn, MA 01801 - USA
P: +1 781 937 0007 - F: +1 781 937 0017
Email: sales.usa@interscience.com

interscience CHINA

上海市徐汇区肇嘉浜路798号坤阳大厦1903室-200030
电话: +86 (0)21-64739390 - +86 189 3097 0733
邮址: sales.china@interscience.com

interscience SOUTH-EAST ASIA

The Centropod, 80 Changi Rd - 05-07 SINGAPORE 419715
T: +65 6909 0825 - M: +65 8118 5178
Email: sales.asia@interscience.com