



Order number 51343153
former order number -

Specifications

pH range 1 - 11
zero point (pH) 7.00 ± 0.25

temperature range 0 - 80 °C
temperature probe -

shaft material glass
membrane glass LoT
membrane resistance < 250 MΩ

reference system ARGENTHAL™
reference diaphragm open junction
reference electrolyte XEROLYT®PLUS polymer

bridge electrolyte -

electrode head S7 (grey)
cable length -
connector -

cell constant (cm⁻¹) -
cell type -

indicating metal -

Storage FRISCOLYT-B®
(51340053)

Special features

high performance puncture electrode:
- XEROLYT®PLUS polymer electrolyte
- open reference junction
- ARGENTHAL™ reference system

Choice of cables

All Mettler-Toledo meters have BNC connectors.
Recommended cable: BNC 1.2m (52300004)

Our cable selector guide offers a large choice of S7 cables for other connectors. Grey S7 connectors are compatible with red S7 connectors. The color is only design-related.