



TECHNICAL DATA SHEET

1) PRODUCT DESCRIPTION

PRODUCT CODE:	902C
GENERAL DESCRIPTION:	SRK 10ML 16X100 PLASTIC
DETAILED DESCRIPTION:	4 Boxes containing 25 individually Paper-Film pouched tube with 10 ml of SRK Solution one plastic applicator

2) PRODUCT COMPONENTS

MEDIUM	COMPOSITION
SRK Solution	Sodium Chloride Potassium Chloride Calcium Chloride Lecithin Sodium thiosulphate Sodium Thyoglycollate Sodium Pyruvate Sodium Bicarbonate Polysorbate 80

APPLICATOR	DIMENSIONAL SPECIFICATIONS	COMPONENTS	RAW MATERIALS
Swab	LENGTH OF THE PLASTIC SHAFT: 144.0 mm \pm 1 mm PLASTIC TIP DIAMETER: 2.5 mm \pm 0.1 STICK DIAMETER HANDLE PART: 2.5 mm \pm 0.2 LENGTH OF THE FIBRE TIP: 16 mm \pm 3 mm TOTAL LENGTH: 146 mm \pm 2mm FIBER TIP DIAMETER: 4.50 mm \pm 1.0	PLASTIC SHAFT	SOLID HIGH IMPACT POLYSTYRENE
		FIBER	RAYON
		GLUE	VINILIC ADHESIVE

TUBE COMPONENTS	DIMENSIONAL SPECIFICATIONS	RAW MATERIALS
CAP	EXTERNAL DIAMETER: 16 mm \pm 1 mm CAP LENGTH: 21 mm \pm 1 mm	POLYETHILEN HIGH DENSITY
TUBE	EXTERNAL DIAMETER: 16 mm \pm 1 mm INTERNAL DIAMETER: 14 mm \pm 1 mm TUBE LENGTH: 101 mm \pm 2 mm TOTAL LENGTH (CAP+TUBE): 103 mm \pm 2 mm	POLYPROPYLEN

This bulletin is for your guidance and is based on information believed to be reliable and should be periodically reviewed.



3) PACKAGING SPECIFICATION

INNER BOX QTY:	25
OUTER BOX QTY:	100
PACKAGING:	PAPER-FILM POUCH
TOP SHEET MATERIAL:	MEDICAL GRADE PAPER PLUS 60GR
BOTTOM SHEET MATERIAL:	PA/PE FILM - 100 μ

4) IMPORTANT INFO

STERILIZATION METHOD	GAMMA IRRADIATION
STORAGE TEMPERATURE	5/25°C – DO NOT FREEZE
SHELF LIFE	24 MONTHS
MEDIUM pH	7.0 \pm 0.5