

	<b>BagMixer® 400 P</b> 400 mL lab blender	Ref. 021 230
	Powerful and reliable. Avoiding all risks of cross-contamination, BagMixer® 400 P is an easy-to-use and powerful lab blender. Adapted to all kinds of applications and with a guarantee of optimal bacterial extraction, it is a great tool for lab analyses.	

## ADVANTAGES

- ▶ Silent and robust
- ▶ Variable blending power
- ▶ Easy cleaning: door opens at 270°

## TECHNICAL SPECIFICATIONS

- ▶ All stainless steel
  - ▶ Plain door
  - ▶ Useful volume: 50 – 400 mL
  - ▶ Fixed blending speed (8 strokes/sec.)
  - ▶ Variable time (30 – 210 secs. or ∞)
  - ▶ Adjustable paddles (patented)
  - ▶ Fully opening door (270°)
  - ▶ Also available in CoolMix® version
- 
- ▶ Power: 110-220 V~ 50-60Hz
  - Machine: size (W x D x H): 26 x 39 x 29 cm, weight: 16.5 kg
  - Box: size (W x D x H): 48 x 33.5 x 33.5 cm, weight: 21 kg
- 
- ▶ 3 year guarantee (after registration)
  - ▶ Shock absorbers life-time guarantee
  - ▶ Life-time service
  - ▶ In compliance with: ISO 7218
  - ▶ Designed and made in France
  - Manufactured under: ISO 9001 V2008, CE, WEEE 2002/96/EC, RoHS 2002/95/EC standards

## DELIVERED WITH

- ▶ Power cord
- ▶ User manual, IQ/OQ/PQ guidelines
- ▶ Quick user guide
- ▶ Free sample pack of bags

## OPTIONS / ACCESSORIES

- ▶ **Ref. 026 000** Drip tray BagMixer®
- ▶ **Ref. 211 040** BagOpen® 400, bag opener stand for holding bags
- ▶ **Ref. 221 040** BagRack® 400, stainless steel stackable storage racks for 10 bags
- ▶ **Ref. 231 040** BagClip® 400, closing clips for lab blender bags 400 mL
- ▶ **Ref. 251 091** BagPipet®, automatic pipet, fixed volume 0.1/0.9/1 ml, ejector for BagTips®

---

## CONSUMABLES

---

- ▶ **Ref. 111 425** BagFilter® 400 P, sterile side filter bag
- ▶ **Ref. 122 025** BagPage® + 400, sterile full-surface filter bag
- ▶ **Ref. 132 025** BagLight® 400 PolySilk, plain blender bag
- ▶ **Ref. 145 040** RollBag® 1300, sterile sampling bag
- ▶ **Ref. 252 019** BagTips® REGULAR, 19 cm sterile tips for BagPipet® or other pipet
- ▶ **Ref. 252 024** BagTips® JUMBO, 24 cm sterile tips for BagPipet® or other pipet