

CORNING	Life Sciences	Document Number: LSR00014	Rev.: 5
Description (Class): Reference, Customer Technical Data Sheet – Serological Pipets, 25 mL Extended 0.5 mL Increments			Page: 1 of 2

APPLICATIONS:

For tissue culture usage or any non-toxic solution that requires aliquoted or aspirating/dispensing in an aseptic environment. Four additional graduations from 32 mL (356525, 357515, 357525) to 36 mL (356535, 357535), which enables researchers the ability to do more aliquoting per pipet.

FEATURES:

- Consistent performance
 - Reliability, accuracy is +/-2% at full volume
- Non-pyrogenic
 - Low endotoxins
- Color coding
 - Ease of identification on packaging and blow-out rings
- Reverse graduations
 - Easy reading of the amount of liquid removed or the amount remaining
- Negative graduations
 - Can measure past the stated capacity
- Crisp, legible, rotogravure printing
 - Easily read sample volumes
- Crystal grade polystyrene
 - Can clearly see fluid sample
- Individually wrapped packaging:
 - Two ways to open: Peel or Pop-through

REGULATORY COMPLIANCE:

Falcon® products are manufactured under the current ISO 9001 and ISO 13485 Standards, and the current FDA Quality System Regulation 21 CFR Section 820.

STERILITY:

Product is gamma irradiated.

RESTRICTIONS:

- Refer to "Thermoplastics Properties Chart," available at:
http://catalog2.corning.com/Lifesciences/media/pdf/an_DL_226_Falcon_ThermoPlastics_Broch.pdf

CORNING	Life Sciences	Document Number: LSR00014	Rev.: 5
Description (Class): Reference, Customer Technical Data Sheet – Serological Pipets, 25 mL Extended 0.5 mL Increments			Page: 2 of 2

Falcon® 25 mL, 0.5 mL Increments, Serological Extended Pipet

Material: Pipet - Polystyrene; Plug - Polyester fiber

Related Catalog Numbers:	Cat. No.	Packout Description
Sterile	356535	Individually wrapped, 50/bag, 200/case
Sterile	357535	Individually wrapped, 50/bag, 200/case

Dimensions: = XXX mm (XXX inches)

TIP:

Length = 39.24 (1.545)
 Interior Diameter = 2.36 (.093)
 Outside Diameter = 4.78 (.188)

MOUTHPIECE:

Length: = 25.91 (1.020)
 Interior Diameter = 5.08 (.200)
 Outside Diameter = 7.62 (.300)
 Plug Length = 14.30 (.563)

EXTRUSION:

Interior Diameter = 13.43 (.529)
 Outside Diameter = 15.16 (.597)
 Length = 274.75 (10.817)

PIPET SIZE: 25 mL

INCREMENTS: 0.50


GRADUATIONS: Reverse

OVERFILL CAPACITY: 7.0 mL

COLOR CODE: Magenta

GENERAL DIMENSIONS:

A = 336.29 (13.240)
 B = 4.78 (.188) Diameter
 C = 7.62 (.300) Diameter
 D = 15.16 (.597) Diameter



Warranty/Disclaimer: Unless otherwise specified, all products are for research use only. Not intended for use in diagnostic or therapeutic procedures.