



FARM
SAMPLING

PROCESS
SAMPLING

PERFORMANCE
ANALYSIS

3 oz Solid Top Vial

One-Piece Flip-Top Vial - Specific for aseptic sampling



Applications

Fields of use

- Specifically designed for field collection of farm milk and related dairy products
- Quality control sample storage containers and food packaging vials
- Veterinary specimen collection, infection control of mastitis
- Environmental sampling containers

Options

- Can be supplied with Bronopol and / or BARCODE label
- Standard Colour is natural white – other colours on request with minimum production run requirements

Benefits

Aseptic single-use vial with flip top construction

- Manufacture in clean room
- Patented « in mould closing » process
- Perfectly tight and aseptic
- Translucent in colour for visual inspection
- Pure Polypropylene suitable for food contact
- 100% Recyclable after use

User-friendly hinged vial concept

- One hand opening and closing for perfect cleanliness and ease of use
- No risk of contamination due to polluted stoppers
- Allows the use of gloves in cold environment

Safety and Quality

- Alimentary : Complies with CE Norm N° 1935/2004
- Traceability of batches of production : a label on each box includes batch number and production date

Characteristics

Internal Volume	75,0 ml
Vial Overall Height	85,6 mm
Bottom Body Diameter	35,1 mm
Top Body Diameter	36,4 mm
Top Lid Diameter	40,3 mm
Body Height	84,5 mm
Inside Top Diameter	34,9 mm
Inside Height	83,5 mm
Filling Level Line	60,0 ml
Vial Weight	10,6 grams

PRODUCTION MOULD # 723

ART. # 03 03 6086 0000 CN

In colour Clear – 400 Pieces per box - 30 Boxes – 12.000 Pieces per Pallet