

# BagMixer<sup>®</sup> 400

The world's best selling lab blender



interscience

# Why choose BagMixer® 400!

## Best-seller

- Sold in more than 80 countries
- 30 years experience

## High-quality

- Optimal analysis with the sterile **BagFilter®** and **BagPage®**
- In compliance with ISO 7218 standard
- Made in France by **interscience**

## Brushless motor

- No wearing parts
- - 50% power consumption
- Extra silent

## Safety

- The handle closes the blender and starts the blending
- **Q-Tight®**: no leaks from the top of the bag
- Security drip tray

## PATENTED

## Click & Clean® removable paddles

- Remove the paddles and easily clean the blending chamber



## Easy cleaning

- All stainless steel
- Door opens at 270°
- Removable **Click & Clean®** paddles

## Guarantee

- Shock absorbers and window door **life-time** guarantee
- 3 year guarantee (after registration)

## PATENTED


## Adjustable blending power

- Up to 28 kg of pressure: the hardest samples don't resist! e.g.: nuts, pills, seeds, dog biscuits etc.



## Optimal extraction

- From as little as 30 seconds
- Easy control with the digital screen (speed, time, power)



More than 10,000 **BagMixer®** in use every day!

Made in France by **interscience**,  
the only manufacturer of blenders and sterile bags for microbiological analysis

Some features are limited to specific models.  
For complete technical info, please check online:  
[www.interscience.fr](http://www.interscience.fr)



# 4 BagMixer® 400 models:

The **BagMixer®** lab blender guarantees optimal extraction by the action of its two paddles. Adapted to any sample (food, agricultural, chemical, medical, environmental) to be analyzed between 50 ml and 400 ml, it guarantees homogenization in 30 to 60 seconds.



## BagMixer® 400 P

> The essential

Ref. 021 230



## BagMixer® 400 W

> Window door

Ref. 022 230



**BEST SELLER**

## BagMixer® 400 VW

> Adjustable speed

Ref. 023 230



## BagMixer® 400 CC

> Click & Clean®

Ref. 024 230

Technical specifications at the back of the leaflet →

# BagMixer® packs:

## Starter Pack



1 BagMixer®

+ Sterile bags: 500 BagFilter®  
+ 1,000 BagLight®

+ Accessories: 50 BagClip®  
+ 1 BagOpen® + 1 BagRack®

## Classic Pack



1 BagMixer®

+ Sterile bags: 500 BagFilter®  
+ 1,000 BagLight® + 500 BagPage®

+ Accessories: 50 BagClip®  
+ 1 BagOpen® + 1 BagRack®  
+ 1 PetriPile® + 1 BagPipet®  
+ 1,000 BagTips® Regular  
+ 6 Clinet® 101

## Premium Pack




1 BagMixer® & 1 BabyGravimat®

+ Sterile bags: 1,000 BagFilter®  
+ 2,000 BagLight® + 1,000 BagPage®

+ Accessories: 50 BagClip®  
+ 1 BagOpen® + 2 BagRack®  
+ 2 PetriPile® + 1 BagPipet®  
+ 1,000 BagTips® Regular  
+ 6 Clinet® 101

## BagSystem® complete range for sample preparation before analysis



**BagFilter®**  
Collect



**Gravimat®**  
Dilute








**BagMixer®**  
Blend



**BagTools®**  
Pipet

# Technical specifications

	BagMixer® 400 P	BagMixer® 400 W	BagMixer® 400 VW	BagMixer® 400 CC®
<b>Reference</b>	021 230	022 230	023 230	024 230
<b>Key features</b>	Removable and autoclavable paddles, <b>Click &amp; Clean®</b> system	-	-	✓
	Blending power indicator	-	-	✓
	Adjustable blending power (up to 28 kg of pressure)	-	-	✓
	Digital screen	-	-	✓
	Leak free <b>Q-tight®</b> system	-	-	✓
	Side-by-side paddle stop	-	-	✓
	Silent brushless motor	-	-	✓
	Window door	-	✓	✓
	Security drip tray	in option	✓	✓
	270° opening door	✓	✓	✓
	All stainless steel	✓	✓	✓
	Compact and ergonomic	✓	✓	✓
<b>Technical specifications</b>	Door weight	1050 g	950 g	
	Blending time	30, 60, 90, 120, 150, 180, 210 sec., ∞		from 1 sec. to 59 min, ∞
	Blending speed	fixed speed: 8 strokes/sec.		variable speed: 4, 6, 8, 10 strokes/sec.
	Dimension (L x l x h)	39 x 26 x 29 cm		42 x 26 x 26.5 cm
	Weight	16.5 kg		
	Power	~100-127V / ~220-240V - 50-60Hz		~100-240V - 50-60Hz
Capacity	50-400 mL			
<b>Guarantees</b>	Made in France	✓	✓	✓
	Life-time guarantee	Unbreakable shock absorbers and window door		
	3 year guarantee (after registration)	✓	✓	✓
<b>Manufactured in compliance with</b>				 

## Discover our complete range for microbiology



Your local distributor:

### interscience INTERNATIONAL

30, Ch. Bois des Arpents - 78860 St Nom - FRANCE  
 T: +33 (0)1 34 62 62 61 - F: +33 (0)1 34 62 43 03  
 info@interscience.fr - www.interscience.fr

### interscience USA & CANADA

199 Weymouth ST - ROCKLAND - MA 02370 - USA  
 T: +1 781 792 2133 - F: +1 781 792 2134  
 info@intersciencelab.com - www.intersciencelab.com