


interscience

30, CH. BOIS ARPENTS – 78860 ST-NOM-LA-BRETECHE – FRANCE – T. +33 (0)1 34 62 62 61 – F. +33 (0)1 34 62 43 03
N° Identification TVA: FR52 950 356 220 – E-mail: info@interscience.fr – Website: www.interscience.fr

	BagMixer® 400 W 400 mL lab blender	Ref. 022 230
<p>Powerful and efficient. Avoiding all risks of cross-contamination, BagMixer® 400 W is an easy-to-use and powerful lab homogenizer. Adapted to all kinds of applications and with a guarantee of optimal bacterial extraction, BagMixer® 400 W, with its window door, ensures permanent control over blending process for perfect results.</p>		

ADVANTAGES

- ▶ Silent and robust
- ▶ Variable blending power
- ▶ Easy cleaning: door opens at 270°

TECHNICAL SPECIFICATIONS

- ▶ All stainless steel
 - ▶ Window door
 - ▶ Useful volume: 50 – 400 mL
 - ▶ Fixed speed (8 strokes/sec.)
 - ▶ Variable time (30 – 210 or ∞)
 - ▶ Variable power, adjustable paddles (patented)
 - ▶ Fully opening door (270°)
 - ▶ Security drip tray
-
- ▶ Power: 110-220 V~ 50-60Hz
 - Machine: size (W x D x H): 26 x 39 x 29 cm, weight: 16.5 kg
 - Box: size (W x D x H): 48 x 33.5 x 33.5 cm, weight: 18.8 kg
-
- ▶ 3 year guarantee (after registration)
 - ▶ Shock absorbers life-time guarantee
 - ▶ Life-time service
 - ▶ In compliance with: ISO 7218
 - ▶ Designed and made in France
 - Manufactured under: ISO 9001 V2008, CE, WEEE 2002/96/EC, RoHS 2002/95/EC standards

DELIVERED WITH

- ▶ Security drip tray
- ▶ Power cord
- ▶ User manual, IQ/OQ/PQ guidelines
- ▶ Quick user guide
- ▶ Free sample pack of bags

OPTIONS / ACCESSORIES

- ▶ **Ref. 211 040** BagOpen® 400, bag opener stand for holding bags
- ▶ **Ref. 221 040** BagRack® 400, stainless steel stackable storage racks for 10 bags
- ▶ **Ref. 231 040** BagClip® 400, closing clips for lab blender bags 400 mL
- ▶ **Ref. 251 091** BagPipet®, automatic pipet, fixed volume 0.1/0.9/1 ml, ejector for BagTips®

CONSUMABLES

-
- ▶ **Ref. 111 425** BagFilter® 400 P, sterile side filter bag
 - ▶ **Ref. 122 025** BagPage® + 400, sterile full-surface filter bag
 - ▶ **Ref. 132 025** BagLight® 400 PolySilk, plain blender bag
 - ▶ **Ref. 145 040** RollBag® 1300, sterile sampling bag
 - ▶ **Ref. 252 019** BagTips® REGULAR, 19 cm sterile tips for BagPipet® or other pipet
 - ▶ **Ref. 252 024** BagTips® JUMBO, 24 cm sterile tips for BagPipet® or other pipet