




PRODUCT CODE: 283462

Ferriin solution 0,025 mol/l (0,025M) for volumetric analysis
redox indicator

$C_{36}H_{24}FeN_6O_4S$

$C_{36}H_{24}FeN_6O_4S$



M.= 692,24

CAS [14634-91-4]

EINECS 238-676-6

TARIC 3822 00 00 00

SYNONYMS: 1,10-Phenanthroline-Iron(II) Sulphate Complex, Tris (1,10-Phenanthroline) Iron(II) Sulphate•

PHYSICAL DATA: Miscible with water D 20/4 1,008 • pH(50 g/l)2,5 •

BIBLIOGRAPHY: Merck Index 12, 7.356

HAZARDOUS:



H: H302 • H412 •

P: P264 • P270 • P273 • P301+P312 • P330 • P501 •

WEIGHT/VOLUME INFORMATION: 1l~1,008 kg 1kg~0,992 l

SPECIFICATIONS:

Suitability: as redox indicator

Passes test