



FARM  
SAMPLING



PROCESS  
SAMPLING



## 10 oz Custody Clear Vial

One-Piece Flip-Top Custody Vial - Specific for aseptic sampling



### Applications

#### Fields of use

- Specific Vial with Tamper Evidence and Lock-Tab for secure sampling
- Hinge can be secured with safety lock tab – exclusive and very easy to use as an anti-counterfeiting device

#### Ideally shaped for

- Field collection of farm milk and related dairy products
- Veterinary specimen collection, infection control of mastitis
- Environmental sampling containers

#### Options

- Can be supplied with Bronopol and / or BARCODE label
- Standard Colour is natural white – other colours on request with minimum production run requirements

### Benefits

#### Aseptic single-use vial with flip top construction

- Manufacture in clean room
- Patented « in mould closing » process
- Perfectly tight and aseptic
- Translucent in colour for visual inspection
- Pure Polypropylene suitable for food contact
- 100% Recyclable after use

#### User-friendly hinged vial concept

- One hand opening and closing for perfect cleanliness and ease of use
- No risk of contamination due to polluted stoppers
- Allows the use of gloves in cold environment

#### Safety and Quality

- Alimentary : Complies with CE Norm N° 1935/2004
- Traceability of batches of production : a label on each box includes batch number and production date

### Characteristics

Internal Volume	343,0 ml
Vial Overall Height	118,5 mm
Bottom Body Diameter	63,0 mm
Top Body Diameter	64,8 mm
Top Lid Diameter	68,5 mm
Body Height	117,4 mm
Inside Top Diameter	62,5 mm
Inside Height	116,0 mm
Filling Level Line	245,0 ml
Vial Weight	37,3 grams

PRODUCTION MOULD # 901

ART. # 03 06 5119 00CU CN

In colour Clear – 90 Pieces per box - 30 Boxes – 2.700 Pieces per Pallet