



Order number 51343110

former order number -

Specifications

pH range 0 - 12

zero point (pH) 7.00 ± 0.25

temperature range 0 - 130 °C

temperature probe -

shaft material glass

membrane glass A41

membrane resistance < 600 MΩ

reference system SteadyForce™

reference diaphragm ceramic

reference electrolyte DPA gel

bridge electrolyte -

electrode head S7 (grey)

cable length -

connector -

cell constant (cm⁻¹) -

cell type -

indicating metal -

Storage 3 mol/L KCl

Special features

high performance electrode:

- SteadyForce™ reference system
- maintenance free
- autoclavable with supplied S7 protective cap

Choice of cables

All Mettler-Toledo meters have BNC connectors.
Recommended cable: BNC 1.2m (52300004)

Our cable selector guide offers a large choice of S7 cables for other connectors. Grey S7 connectors are compatible with red S7 connectors. The color is only design-related.