


pH



<b>Order number</b>	51343050
---------------------	----------

former order number	-
---------------------	---

**Specifications**

pH range	0 - 14
----------	--------

zero point (pH)	7.00 ± 0.25
-----------------	-------------

temperature range	0 - 100 °C
-------------------	------------

temperature probe	-
-------------------	---

shaft material	glass
----------------	-------

membrane glass	HA
----------------	----

membrane resistance	< 600 MΩ
---------------------	----------

reference system	ARGENTHAL™, Ag <sup>+</sup> trap
------------------	----------------------------------

reference diaphragm	ceramic
---------------------	---------

reference electrolyte	3 mol/L KCl
-----------------------	-------------

bridge electrolyte	-
--------------------	---

electrode head	S7 (grey)
----------------	-----------

cable length	-
--------------	---

connector	-
-----------	---

cell constant (cm <sup>-1</sup> )	-
-----------------------------------	---

cell type	-
-----------	---

indicating metal	-
------------------	---

<b>Storage</b>	3 mol/L KCl
----------------	-------------

**Special features**

ARGENTHAL™ reference system with Ag <sup>+</sup> trap
---

**Choice of cables**

All Mettler-Toledo meters have BNC connectors.
--

Recommended cable: BNC 1.2m (52300004)
--

Our cable selector guide offers a large choice of S7 cables for other connectors. Grey S7 connectors are compatible with red S7 connectors. The color is only design-related.