

**FALCON®**

*A Corning Brand*

## Falcon® Pipet Controller

The light, ergonomic, motorized pipet controller is designed for use with glass and plastic serological pipets from 0.5 to 100 mL volume ranges. Conveniently positioned switches change operating modes and speeds quickly to handle different liquid volumes and viscosities. Aspirating and dispensing speed is controlled through the finger triggers. When fully charged, the batteries enable up to 8 hours of continuous use. The large LCD display clearly indicates battery status, pipetting mode and speed.

The Falcon Pipet Controller is supplied with a two-position charging stand, a universal battery charger, and three batteries. Two 0.2 µm and two 0.45 µm hydrophobic polytetrafluoroethylene (PTFE) filters are supplied with the pipet controller. Additional filters are available as standard accessories. The filters, pipet holder and nosepiece are autoclavable.

### Features

Ergonomic design for exceptional comfort

LCD display for continuous visualization of speed, modes and battery status

Easy access to the battery compartment

Easy maintenance




## Falcon® Pipet Controller Ordering Information

Cat. No.	Description
357470	Falcon Pipet Controller, with two 0.2 µm filters and two 0.45 µm filters, 2-position charging stand, universal power supply and set of 3 batteries

### Accessories

357472	Hydrophobic PTFE filter (0.2 µm), set of 5
357473	Hydrophobic PTFE filter (0.45 µm), set of 5
357474	Silicone pipet holder
357486	Battery (set of 3)



Two-position  
Charging Stand



SACCO srl  
Via Manzoni, 29A  
22071 – Cadorago (CO)  
Tel. 0318866611 Fax. 031904596  
Internet: [www.saccosrl.it](http://www.saccosrl.it) Email: [info@saccosrl.it](mailto:info@saccosrl.it)

CORNING | FALCON® cellgro® AXYGEN® PYREX® GOSSELIN®