


Scan® 1200 *High Resolution*

Automatic color colony counter

*Counts up to 7 colors
on the same plate !*




***Instant, complete & traceable results on all media
(Petri dishes, Petrifilm™, RIDA™Count / Sanita-kun™, CompactDry™...)***



Scan[®] 1200: User friendly, simple & high performance

- 1** **Visualisation window**
Detected colonies are highlighted and marked with a cross
- 2** **Parameter window**
Only one cursor to optimize counting.
Customizable settings adaptable to any kind of agar
- 3** **Results window**
All previously recorded data can be checked again



Scan[®] 1200 counts in a second



A full range of accessories



- One cursor to optimize counting
- 6 combinations of lighting and background color; Lighting (upper / lower) / Background (black / white)
- Detects and counts up to 7 colors on the same plate
- Polygonal exclusion areas
- Minimum size of detected colonies: 0,05 mm

- Firewire CCD camera (1280 x 960 pixels)
- Lighting by lifetime lasting LEDs
- English/French/German/Japanese/Chinese/Russian software
- All stainless steel
- Dim : 27 x 27 x 38 cm - 8,7 Kg
- Power : 220 / 110 V - 50 / 60 Hz

Delivered with: Scan[®] software , Firewire PCI card - Option: Firewire PCMCIA card
Minimum requirements: Windows[™] XP or Vista[™] with Firewire connector (or free PCI slot)

ISO 9001 : 2000 / RoHS / WEEE

Please fax this form to interscience or to your local distributor

Scan[®] 1200 colony counter

Quotation

Demonstration

Literature

I would like more information on:

- BagMixer[®] 400 ml blender
- MiniMix[®] 100 ml blender
- JumboMix[®] 3500 ml blender
- BagFilter[®] filter bags
- BagTools[®] accessories
- Gravimat[®] dilutor
- Spiral[®] DS+ plater

Name.....
Company.....
Address.....
.....
Zip/City..... Country.....
Phone..... Fax.....
email.....
web.....

distributor

interscience

Worldwide: 30, Ch Bois Arpents, 78860 St Nom FRANCE
Tel. +33(0)1 34 62 62 61 / Fax. +33(0)1 34 62 43 03
info@interscience.fr / www.interscience.fr

USA/Canada: 199 Weymouth St. ROCKLAND, MA 02370 USA
Tel. +1 781 792 2133 / Fax. +1 781 792 2134
info@intersciencelab.com / www.intersciencelab.com