

# Disinfectant wipes for food-contact surfaces

## DASR -STERICID

---

### Homologation for food contact without rinsing

Authorization product: AMM BTR0307

### Packaging

Canister of 100 blue wipes (12 canisters / carton).

Wipes size: 120 x 200mm.

### Composition

Ethanol (CAS: 64-17-5): 54% m/m.

### Microbiological properties

Ready to use, 20°C:

- Bactericid on *Salmonella typhimurium* :
  - EN1276 in 30 sec. (dirty conditions).
  - EN1276 and EN13697 in 5 min (clean conditions).
- Bactericid on *Listeria monocytogenes*:
  - EN1276 in 30 sec. (clean conditions).
  - EN13697 in 5 min (clean conditions).
- Levuricid on *Candida albicans*:
  - EN1650 in 30 sec (dirty conditions).
  - EN1650 and EN13697 in 5 min (clean conditions).
- Fungicid on *A.brasiliensis* (ex.niger):
  - EN1650, EN13697 on 5 min (Clean conditions).
- Virucid on roativus:
  - EN14476 in 5 min (dirty conditions)

### Direction for use

Pull out the first central wipe.

Ready-to-use and single use wipes. Do not rinse. Use the wipe on the equipment to be disinfected.

Leave on for a minimum of 5 minutes.

Close lid carefully.

### Cautions

Flammable product.

Never use like a skin disinfectant.

Professional use only.

Consult safety data sheet before first use.

Biocid product (type 4 – Directive 98/8).



---

REF : DASR100

Technical data sheet. Version - 04/2016

STERIDIS - ZAC de la Blanche Tâche - 80450 CAMON

Tél : 0033 (0) 3 22 09 41 39 - Fax : 0033 (0) 3 22 50 01 34 - mail : contact@steridis.com