



**Order number** 51344024  
 former order number

**Specifications**  
 pH range -  
 zero point (pH) -

temperature range 0 - 100 °C  
 temperature probe NTC 30 kΩ

shaft material stainless steel  
 membrane glass -  
 membrane resistance -

reference system -  
 reference diaphragm -  
 reference electrolyte -

bridge electrolyte -

electrode head -  
 cable length 1.2 m  
 connector Mini-DIN

cell constant (cm<sup>-1</sup>) 0.105  
 cell type 2 steel poles

indicating metal -

**Storage** dry

**Additional specifications**  
 measuring range 0.001 - 500 μS/cm

Scale 1:1